

APPLICATIONS



Veterinary Equipment



Surgical Equipment



Laboratory Equipment

PRODUCTS



# STERILIZERS

## SELECTION GUIDE

### HARVEY STERILIZERS

	Product Description	Sterilizing Agent	Cycle Types	Dry Cycle	Printer	Page
<i>The New Dry Standard in Sterilization!</i>	PV Dry™	Steam	Liquid 15 min. @121°C Unwrapped 3 min. @ 135°C Wrapped 5 min. @ 135°C Packs 30 min. @ 121°C Special 5 min. @ 134°C	No Yes Yes Yes Yes	Optional	276
<i>Largest Portable Laboratory Model Available!</i>	SterileMax™	Steam	Liquid 15 min. @ 121°C Unwrapped 3 min. @ 135°C Wrapped 10 min. @ 135°C Packs 30 min. @ 121°C	No Yes Yes Yes	Optional	277
<i>There's No Better Way to Extend the Life of Your Instruments!</i>	Chemiclave	Chemical Vapor	Fixed 20 min. @ 132°C	N/A	No	278
<i>Easy-To-Program Tabletop Sterilizers!</i>	MC8/MC10	Steam	Liquid 30 min. @ 121°C Unwrapped 3 min. @ 135°C Wrapped 5 min. @ 135°C Packs 30 min. @ 121°C	No Yes Yes Yes	Optional	279

### WATER PRETREATMENT SYSTEM

	Product Description	Flow Rate LPM	Page
<i>Quick, Pure Water for Your Steam Sterilizer!</i>	Harvey DI+	0.5	280

### BIOLOGICAL INDICATORS AND STERILIZATION BAGS

	Product Description	Page
<i>Fast Indication and Easy to Use!</i>	B/T Sure®/Sterileamp® II	281
<i>Simple Daily Monitor for Your Sterilizer!</i>	Chemitest Indicators/Sterilization Bags	282

CELEBRATING  
**125 YEARS**  
of INNOVATION, QUALITY AND  
SERVICE SINCE 1878

Barnstead International