


Harvey Sterilizers



Chemiclave Chemical Vapor



EC5500



There's No Better Way to Extend the Life of Your Instruments!

- No drying cycle required
- No rust or corrosion of instruments
- ☆ 24 Month Warranty

Product Description

- Automatic system starts sterilization cycle when unit reaches proper sterilization conditions.
- Electronic display times the 20-minute exposure phase and the purge phase.
- Automatic purge for removal of vapor from chamber before door is opened.
- No drying phase required.
- Easy to fill system for sterilant; prevents spillage.
- Built in chemipurge and chemifilter to slash emissions.
- Dry packs and instruments at end of cycle.
- Preheated chamber reduces overall cycle time.
- ASME coded chamber ETL, CSA approved.
- Operates at 132°C and 20-40 psi.
- Sterilization cycle of 20 minutes.
- Easy-to-remove rack for quick cleaning.
- Proprietary solutions ensure consistent quality results.



Control Panel



Vapo-Steril



Built-in Chemifilter

PRODUCT SPECIFICATIONS

	Chemiclave EC5500	Chemiclave EC6000		Chemiclave EC5500	Chemiclave EC6000
Overall Dimensions (W x H x D)	16.75" x 13" x 18.5" (42.5 cm x 33 cm x 47 cm)	19.25" x 16.5" x 20.5" (48.9 cm x 41.9 cm x 52.1 cm)	Heater Wattage	800 Watts	1200 Watts
Chamber Dimensions	8" dia. x 13.25" deep (20.3 cm x 33.7 cm)	10" dia. x 16" deep (25.4 cm x 40.6 cm)	Fixed Cycle Parameters	20 min. @ 132°C	20 min. @ 132°C
Reservoir Capacity	37" fluid oz. (1.1 liters)	37 fluid oz. (1.1 liters)	Operating Pressure	20 - 40 psi (138-276 KPA)	20 - 40 psi (138-276 KPA)
			Operational Weight	59 lb. (26.8 kg)	120 lb. (54.55 kg)
			Shipping Weight	110 lb. (49.9 kg)	161 lb. (73 kg)

ORDERING INFORMATION

Model #	Description	US List Price	Model #	Description	US List Price
266270	EC5500 Harvey Chemiclave Chemical Sterilizer 120V, 60 Hz, 8 amp, 8" chamber	\$4,723.00	266200*	EC5500 Harvey Chemiclave Chemical Sterilizer 230V, 50 Hz, 4 amp, 8" chamber	\$4,723.00
<i>(All units include a chemifilter, chamber tray, tray handle, drain tube, chamber cleaner, indicator bags, tray liners and Chemitest indicators.) Vapo-Steril Solution sold separately.</i>					
<i>*Units have European cord set.</i>					
			266440	EC6000 Harvey Chemiclave Chemical Sterilizer 120V, 60 Hz, 12 amp, 10" chamber	\$5,904.00
			266400*	EC6000 Harvey Chemiclave Chemical Sterilizer 230V, 50 Hz, 6 amp 10" chamber	\$5,904.00

ACCESSORIES

Model #	Description	US List Price	Model #	Description	US List Price
264253	Vapo-Steril® Solution (Case of four 4-packs)	\$145.00	260006701	Door Gaskets for Chemiclave EC6000, EM6000	\$41.80
261823	Dispenser Cap Assembly	\$12.10	280000301	Door Gaskets for Chemiclave EC5500, EM5500 (Replace once per year to assure best seal)	\$21.00
231600201	Tray for Chemiclave EC6000	\$114.70	Tray Liners		
231550201	Tray for Chemiclave EC5500	\$84.60	231600401	Chamber Tray for Chemiclave EC6000 16" (40.6 cm) x 14" (35.6 cm) Box of 250	\$39.00
261569	Chemi-filter (Pkg. of 1)	\$60.40	231550401	Chamber Tray for Chemiclave EC5500 12.75" (32.4 cm) x 11.5" (29.2 cm), Box of 500	\$40.00
264255	Chemi-filter (Pkg. of 4) for Chemiclave EC5500 and EC6000	\$194.40	231603001	See-Thru Tray: 8" (20.3 cm) x 11.125" (28.3 cm), Box of 250	\$33.00