

Harvey Sterilizers



Water Pretreatment System



Harvey DI shown with optional bench stand

Quick Pure Water for Your Steam Sterilizer!

- Removes calcium, magnesium, and other ionic impurities that can harm your steam sterilizer.
- Easy to install
- Low cost operation

Product Description

- Easy to set up—quick fitting connectors connect directly to your faucet—requires no power.
- Compact—fits on the countertop or mounts to the wall for easy installation.
- Fast—produces 0.5 liters of pure water per minute.
- Reliable—Barnstead has been in the pure water market for over 100 years—trust the water expert.
- Low Cost—50 cents per gallon cost based on average tap water (171 ppm).
- No need to store containers of distilled water.
- Connect a Harvey DI system to your faucet and pure water is available at your request.
- The new Harvey DI system incorporates nuclear grade resin enclosed in virgin polypropylene canister to purify tap water.

Operation

- Cartridge resin changes color when it is time to install new cartridges.
- Simple Maintenance—remove cover for easy access when replacing cartridges.
- The Harvey DI water system provides pre-treated water for your steam sterilizer at the push of a button.

ORDERING INFORMATION

Model #	Description	US List Price
D11441	Harvey DI	\$419.00
D50279	Harvey DI cartridge kit	\$64.00
BC1144X2A	Optional bench stand	\$41.60