


Harvey Sterilizers



Biological Indicators



Sterilamp® II

An Easy and Inexpensive Method to Verify Steam Sterilization Processes in Your Lab

- Works with any steam sterilizer
- Results proven in 48 hours
- Save up to 60% of the cost of outside laboratory sterility tests



B/T Sure® Biological Indicator Incubator



B/T Sure Biological Indicator

Product Description

- Verify the sterility of your instruments, vials, media, liquid cycles and other products in less time than it takes to send samples to an outside lab.
- Obvious color change from purple to yellow indicates an ineffective sterilization process.
- Two biological indicators to choose from: 1) for media and liquids 2) for instruments, glassware and other applications.

For Media and Liquids

- Sterilamp II Biological Indicators with spore populations of 10^6 .
- Available in boxes of 25 or 100.

For Instruments, Glassware and Other

- B/T Sure Biological Indicators with spore populations of 10^5 .
- Other spore populations are available upon special request.
- Comply with current AAMI and ISO standards.

**Includes FREE Incubator
With Purchase of 500 Indicator Vials**

ORDERING INFORMATION

HARVEY STERILIZATION PROCESS VALIDATION INDICATORS		
Model #	Description	US List Price
Biological Indicators for Instrument Glassware and Other Product Verification		
AY759X1	B/T Sure Biological Indicator (Box of 5 Indicator vials) Spore population 10^5	\$19.00
AY759X3	B/T Sure Biological Indicator (Box of 100 Indicator vials) Spore population 10^5	\$291.00
AY759X8	B/T Sure Biological Indicator (Box of 100 Indicator vials) Spore population 10^6	\$306.00
Biological Indicators for Media and Liquid Cycle Verification.		
AY759X4	Sterilamp II Biological Indicators (Box of 100) Spore population 10^5	\$303.00
AY759X7	Sterilamp II Biological Indicators (Box of 25) Spore population 10^5	\$95.00
Incubator for Processing Harvey Biological Indicators		
DB104115	B/T Sure Biological Indicator Incubator, 120V	\$286.00
DB104110-33	B/T Sure Biological Indicator Incubator, 240V, CE marked	\$287.00

*Sterilamp® II is a registered trademark of SGM Biotech.