

Thermolyne Stirrers



Cosmo™ Ministirrer



Model
PS61047X1

Low-Cost Personal Size Stirrers!

- Stylish, whimsical design
- Powerful stirring up to 2000 rpm
- Standard, digital and portable versions

Product Description

MiniStirrers

- Combine futuristic design with excellent value.
- Manufactured from a combination of tough Hytrel® bases and clear, chemically-resistant, Nylon® tops.
- Includes a bright power-on LED.
- Powerful enough to stir 1L of liquid at 2000 rpm.
- Digital version has a precise readout of the actual speed in rpm making it ideal for defined standard operating procedures.
- Upon request your own logo or message may be added beneath the clear Nylon top where it is both clearly visible yet protected. (Minimum order of 25 required.)

Portable MiniStir

- Fitted with re-chargeable batteries; runs for up to 16 hours continuously.
- Ideal for use in field situations, fume hoods, glove boxes and incubators.
- All models come complete with instructions, stir bar and line adapter/charger.

ORDERING INFORMATION

Model #	Description	Volts	Overall Dimensions Inches (cm)			Maximum Vessel Capacity (l)	Maximum Vessel Diameter (cm)	Speed Range rpm	Shipping Weight Lb. (kg)	US List Price
			W	H	D					
PS61047X1	Standard MiniStirrer	115	5.6 (14.3)	2.3 (6.6)	5.6 (14.3)	1	8.5	350-2000	2 (1)	\$102.00
PS61048X1	Digital MiniStirrer	115	5.6 (14.3)	2.3 (6.6)	5.6 (14.3)	1	8.5	350-2000	2 (1)	\$133.00
PS61049X1	Portable MiniStirrer	115	5.6 (14.3)	2.3 (6.6)	5.6 (14.3)	1	8.5	350-1400	2 (1)	\$165.00
PS61047X6	Standard MiniStirrer	230	5.6 (14.3)	2.3 (6.6)	5.6 (14.3)	1	8.5	350-1400	3 (1.5)	\$102.00
PS61048X6	Digital MiniStirrer	230	5.6 (14.3)	2.3 (6.6)	5.6 (14.3)	1	8.5	350-2000	3 (1.5)	\$133.00
PS61049X6	Portable MiniStirrer	230	5.6 (14.3)	2.3 (6.6)	5.6 (14.3)	1	8.5	350-1400	3 (1.5)	\$165.00