

## Thermolyne Stirrers



## Student



Models S7800 and S17400

**Includes FREE Stirbar**

### Economical Student Stirrer!

- Motor is designed for hours of continuous operation
- Model S7800 features a solid stainless steel case
- Adjustable stirring speed from 100-1000 rpm

### Product Description

- Perforated steel case on Type S17400 provides cool operation.
- Strong magnetic coupling keeps stir bar synchronized.
- Stirring motor is designed for continuous hours of operation.
- Type S7800 features separate ON/OFF switch permitting starts and stops without changing the speed setting.
- Type S7800 also features a dovetail bracket for mounting on laboratory racks.

### Applications

- General laboratory mixing of aqueous solutions
- Used for stirring while adjusting pH of solutions

### PRODUCT SPECIFICATIONS

Type	Top Plate Dimensions Inches (cm)		Overall Dimensions Inches (cm)		Shipping Weight Lb. (kg)
	Dia.		H	Dia.	
S7800	4 (10.1)		3 (7.6)	4.5 (11.43)	2 (1)
S17400	4 (10.1)		3 (7.6)	4 (10.1)	2 (1)

### ORDERING INFORMATION

Model #	Volts	Electrical (50/60 Hertz)		Stirring Speed Range (rpm)	US List Price
		Amps	Watts		
S7805	120	0.15	18	100-1000	\$218.00
S17415	120	0.15	18	100-1000	\$178.00