

Lab-Line Stirrers



Single Position



Model 1150



Model 1170



Model 1250

Economical Design!

- High speed stirring
- Maximizes bench space
- ☆ 24 Month Warranty



Product Description

- Variable speed heavy duty motors for longer life.
- Compact design maximizes use of available bench space.
- Designed for vessels up to 1 liter capacity.
- Models 1150 and 1250 are sealed to prevent sample contamination of internal components.

Application

- General stirring

PRODUCT SPECIFICATIONS

Model #	Top Plate Dimensions Inches (cm)			Overall Dimensions Inches (cm)		Stirring Speed RPM*
	W	D	W	H	D	
1150	5.4 (14)	3.4 (9)	6.25 (16)	4 (10)	5.1 (13)	60-2250
1170	8 (20)	7 (18)	8.25 (21)	4.4 (12)	8.6 (22)	60-2250
1250	4 (10)	4 (10)	5 (13)	5 (13)	5 (13)	60-2250

*Stirring speed will be affected by liquid viscosity, size of stirring bar used, vessel size and the thickness of the bottom of the vessel.

ORDERING INFORMATION

Model #	Volts	Electrical Hz	Watts	Shipping Weight Lb. (kg)	US List Price
1150 ^{1,2}	120	60	96	6 (3)	\$230.00
1150-1 ^{3,4}	240	50	96	6 (3)	\$230.00
1170	120	50/60	75	11 (5)	\$380.00
1170-1	240	50/60	75	11 (5)	\$380.00
1250	120	50/60	75	6 (3)	\$222.00
1250-1	240	50/60	75	6 (3)	\$222.00

¹ UL, ² CSA, ³ CE, ⁴ ETL