

Lab-Line Stirrers



Multi-Position



Model 1278



Model 1262

Continuous Stirring of Multiple Vessels!

- Simultaneous stirring of 1-6 vessels (up to 1 l flasks)
- Designed for continuous operation
- Efficient operation up to 40°C in an incubator
- ★ 24 Month Warranty



Product Description

- Solid state speed control.
- Fan-cooled, sealed ball bearing motors give synchronous, uniform speeds.
- Sturdy metal housing has Neoprene™ feet for stability and bench top protection.
- Separate power switch on six position models can be turned on or off without changing the stirring dial.
- Stirring speed range: 100-1200 rpm.

Applications

- Ideal for pharmaceutical, chemical, biological and industrial laboratories where long-term mixing of fluids is required
- Biochemical analysis
- Preparation of reagents
- General stirring applications

PRODUCT SPECIFICATIONS

Model #	W	Overall Dimensions Inches (cm)		Shipping Weight Lb. (kg)
		H	D	
1262	11 (28)	6.9 (17)	10 (25)	21 (10)
1262-1	11 (28)	6.9 (17)	10 (25)	21 (10)
1278	18 (46)	7 (18)	14 (36)	36 (16)
1278-1	18 (46)	7 (18)	14 (36)	36 (16)
1288	18.5(47)	7(18)	14(36)	18(8)

ORDERING INFORMATION

Model #	Volts	Electrical		US List Price
		Hz	Watts	
1262	120	60	200	\$528.30
1262-1CE ¹	240	50	200	\$1,100.00
1278	120	60	180	\$1,380.00
1278-1	240	50	180	\$1,380.00
1288 ²	120	60	180	\$1,435.00
1288-1CE ¹	240	50	180	\$1,435.00

¹ CE, ² ETL