



Largest Volume

Largest Volume Stirrer!

- Stirring capacities up to 200 L (53 gallons)
- Microprocessor-based speed control
- ☆ 24 Month Warranty



Model 1297



Model 1295

Product Description

- Speed up to 2400 rpm (maximum speeds vary with fluid viscosity, vessel size and construction and stir bar size).
- LED speed display.
- Digital speed settings.
- Accommodates stir bars from 1.5" to 6" (4 to 15 cm), sold separately.
- Polished stainless steel construction resists corrosion.
- 18" x 18" (46 x 46 cm) top surface has ridged non-slip white rubber mat.
- Model 1297 is remote-controlled by a 10 ft. (3m) cable.

Applications

- Viscous gels
- Latex mixtures
- Polymers
- Starches
- Molasses

PRODUCT SPECIFICATIONS

Model #	Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	
1295	18 (46)	10.25 (26)	18 (46)	45 (20)
1295-1	18 (46)	10.25 (26)	18 (46)	45 (20)
1297	18 (46)	10.25 (26)	18 (46)	50 (23)
	8.4 (21)	Remote 4.4 (11)	5.5 (14)	
1297-1	18 (46)	10.25 (26)	18 (46)	50 (23)
	8.4 (21)	Remote 4.4 (11)	5.5 (14)	

ORDERING INFORMATION

Model #	Volts	Electrical 50/60 Hz		US List Price
		Amps	Watts	
1295	120	1.6	190	\$2,837.00
1295-1	240	0.8	190	\$2,837.00
1297	120	1.6	190	\$3,030.00
1297-1	240	0.8	190	\$3,030.00