



Distribution Pumps



Product Description

- Barnstead's distribution pumps are used with purified water storage tanks when it is necessary to keep the pure water under pressure for delivery to locations above tank level or at greater than gravity flow pressure.
- Choose from four sizes.
- Includes pump bypass with restrictive orifice for tank.
- Includes low-level cutoff switch to shut pump off on storage tank low water conditions.
- Purity is maintained with nonreactive wetted parts.
- Heavy duty, totally enclosed, split-phase induction motor.
- Motor starters not included; all external wiring by customer.
- Bypass piping is precut to fit on tin-lined tanks only.



ORDERING INFORMATION

Description	Model #	US List Price	Model #	US List Price	Model #	US List Price	Model #	US List Price
Assembly Model No. ¹	Tin Bypass H1110	\$2,005.00	H1120	\$2,946.00	H1130	\$3,669.00	H1140	\$4,223.00
	Plastic Bypass H1112	\$2,196.00	H1122	\$2,775.00	H1132	\$3,501.00	H1142	\$4,071.00
Pump Size	AA		BB		CC		DD	
Pump Part No.	01398		01429		01430		01431	
Pump Suction (NPT)	1/2"		3/4"		3/4"		3/4"	
Pump Discharge (NPT)	1/4"		1/2"		1/2"		1/2"	
Distribution Outlet (NPT)	3/4"		3/4"		3/4"		3/4"	
Bypass Return (NPT)	1/4"		1/4"		1/4"		1/4"	
Pump Motor	HP	1/3 HP	1/2 HP		1 HP		1 1/2 HP	
	Speed (rpm)	3450	3450		3450		3450	
	VAC	115/230	115/230		115/230		230/460	
	Frequency (Hz)	50/60	50/60		50/60		50/60	
	Phase	1	1		1		3	

¹Pump starters are not provided with units.