

ERTCO
Thermometers



Non-Mercury Thermometers



Model 24130S

Exact-Temp™ Thermometers for Dry-Block Heaters

- Thermometer is traceable to N.I.S.T.
- Blue spirit-filled
- Thermometer safety-coated with Teflon®
- Housed in 1.5 ml micro tube
- Liquid in tube—mineral oil
- Individually serialized

High-Precision Spirit-Filled Thermometers for Laboratory Instruments

- Red spirit-filled thermometers
- Accurate to within N.I.S.T. tolerances
- Individually serialized
- Glass or Teflon® coated
- Variety of ranges and immersions
- Some with immersion enlargements

NOTICE

"ERTCO Thermometers are available through your local laboratory distributors. Contact Ever Ready Thermometer Company for technical information at 800-453-7826 or 973-812-7474."

ORDERING INFORMATION

Model #	Range	Divisions	Length (mm)	Immersion (mm)	Special Features	Applications	Teflon Coated	US List Price
High-Precision Spirit-Filled Thermometers for Laboratory Instruments								
711S	20° to 100° C	0.5° C	305	76	1, 2	Incubator, Water Bath	No	\$43.90
711-SSC	20° to 100° C	0.5° C	305	76	1, 2	Incubator, Water Bath	Yes	\$49.90
641S	0° to 80° C	0.5° C	250	100	—	Incubator, Water Bath	No	\$32.00
641-SSC	0° to 80° C	0.5° C	250	100	—	Incubator, Water Bath	Yes	\$39.00
639-1S	0° to 70° C	0.5° C	305	150	1	Incubator, Water Bath	No	\$37.40
639-1SSC	0° to 70° C	0.5° C	305	150	1	Incubator, Water Bath	Yes	\$45.90
2457S	25°-37°-56° C	0.5° C	240	76	2	Incubator, Water Bath	No	\$30.80
2457-SSC	25°-37°-56° C	0.5° C	240	76	2	Incubator, Water Bath	Yes	\$39.90
717S	-5° to 20° C	0.5° C	152	Total	2	Blood Bank, Refrigerator	No	\$25.00
717-SSC	-5° to 20° C	0.5° C	152	Total	2	Blood Bank, Refrigerator	Yes	\$32.70
715S	25° to 57° C	0.5° C	175	35	2	Block Heater, Water Bath	No	\$42.50
715-SSC	25° to 57° C	0.5° C	175	35	2	Block Heater, Water Bath	Yes	\$50.00
647S	0° to 110° C	0.5° C	225	35	—	Block Heater, Water Bath	No	\$42.10
647-SSC	0° to 110° C	0.5° C	225	35	—	Block Heater, Water Bath	Yes	\$49.80
Exact-Temp™ Thermometers for Dry Baths and Block Heaters								
DB-030-1	18° to 60° C	0.5° C	125	35				\$55.60
DB-040-1	50° to 110° C	0.5° C	150	35				\$55.60

ERTCO Thermometers



Non-Mercury Thermometers



Model OV-700-01S



Model PS100F

Digital Min-Max Thermometer

- Dual temperature reading - ambient and remote
- °C or °F switchable
- Minimum temperature memory
- Maximum temperature memory
- Min-Max alarm settings
- 9-foot water-resistant sensor cable for remote
- Wall-mounting hole, desk stand and magnetic stand

Exact-Temp™ Certified Chamber Thermometers

- For freezers, refrigerators, incubators and ovens
- Thermometers are graduated in 0.1° C and 0.5° C
- Thermometers traceable to N.I.S.T.
- Furnished with traceable certificate
- Thermometers coated with Teflon®
- Thermometers in plastic bottle with insulating liquid
- Magnets attached to bottle for wall or door mounting



Model R-200-01S

NOTICE

"ERTCO Thermometers are available through your local laboratory distributors. Contact Ever Ready Thermometer Company for technical information at 800-453-7826 or 973-812-7474."

ORDERING INFORMATION

Model #	Range	Divisions	Length (mm)	Immersion	Bottle Liquid	Application	US List Price
0.1 Graduated Thermometers							
F-100-01S	-25° to -5°C	0.1°C	210	Total	Propylene Glycol	Certified Freezer	\$65.90
R-200-01S	-2° to 10°C	0.1°C	180	Total	Ethylene Glycol	Certified Refrigerator	\$65.90
BB-210-01S	-2° to 10°C	0.1°C	180	Total	Ethylene Glycol	Certified Blood Bank	\$84.50
I-300-01S	15° to 30°C	0.1°C	180	Total	Ethylene Glycol	Certified Room/Incubator	\$65.90
I-310-01S	25° to 45°C	0.1°C	190	Total	Ethylene Glycol	Certified Incubator	\$65.90
OV-700-01S	95° to 115°C	0.1°C	200	Total	Vermiculite	Certified Oven	\$65.90
0.5 Graduated Thermometers							
ULF-010S	-90° to 20°C	1.0°C	145	Total	Propylene Glycol	Certified Ultralow Freezer	\$57.40
F-010-1S	-30° to 0°C	0.5°C	115	Total	Propylene Glycol	Certified Freezer	\$45.50
R-020-1S	-5° to 15°C	0.5°C	125	Total	Ethylene Glycol	Certified Refrigerator	\$45.50
BB-100S	-5° to 15°C	0.5°C	145	Total	Ethylene Glycol	Certified Blood Bank	\$70.70
I-030-1S	18° to 50°C	0.5°C	135	Total	Ethylene Glycol	Certified Incubator	\$45.50
OV-070S	20° to 130°C	1.0°C	135	Total	Vermiculite	Certified Oven	\$45.50
Digital Min-Max							
Model #	Range	Remote	Accuracy	Resolution	Display Size	Battery	US List Price
ERT800	Ambient -20° to 70°C and -4° to 158°F	-50° to 70°C and -58° to 158°F	1°C or °F	0.1°C or °F	35mm (1 1/2" x 15mm (0.6"))	AAA	\$29.90