

Thermolyne Vortex Mixers



Vortex Maxi Mix® III



Versatile Vortex Mixer!

- Four interchangeable accessories
- Variable speed control - 100-2200 rpm for gentle or vigorous mixing.
- Continuous, trouble-free shaking of loads up to 5 lb. (2.2 kg)

Product Description

- The Maxi Mix III combines the strength of a heavy cast aluminum base with proven counterbalance system. Spring dampened feet eliminate creeping/walking of unit during large volume mixing.
- Supplied complete with one soft foam top (PT500X6A), a three-wire cord, and plug.

Operation

- Four mixing and shaking accessories perform the function of four separate apparatus, saving valuable bench space in crowded labs and clinics. Choose from foam pad, universal holder, utility tray or flask holder.
- Precise matching of an electronic speed control to a high torque DC type motor provides versatile speed selection ranging from 100 to 2200 rpm.
- Lighted main power switch indicates mixer is operating.



Type M65800 shown with PT500X6A (Foam Pad)

PT500X7A (Utility Tray)

PT500X9A (Universal Holder)

PT500X8A (Flask Holder)

PRODUCT SPECIFICATIONS

Type	Foam Top Mixing Surface Dimensions Inches (cm)		W	Overall Dimensions Inches (cm)		Shipping Weight Lb. (kg)
	D	W		H	D	
M65800	5.5 (13.9)	6.3 (15.8)	6.5 (16.5)	5.1(13.0)	8.5 (21.5)	17.0 (7.7)

Orbital Diameter .156" (4 mm)

ORDERING INFORMATION

Model #	Volts	Electrical Amps	Hz	Operating Speed (rpm)	US List Price	Model #	Volts	Electrical Amps	Hz	Operating Speed (rpm)	US List Price
M65825 ¹	120	0.4	50/60	100-2200	\$924.00	M65820-33 ²	230	0.2	50/60	100-2200	\$924.00

¹ CSA marked
² CE marked

ACCESSORIES

Model #	Description	Mixing Surface Dimensions Inches (cm)		US List Price	Model #	Description	Mixing Surface Dimensions Inches (cm)		US List Price
		W	D				W	D	
PT500X6A	Soft foam rubber top	5.5 (13.9)	6.3 (15.8)	\$63.00	PT500X8A	Flask holder for four 250 ml flasks	—	—	\$209.00
PT500X7A	Utility tray	15.8 (40.1)	7.8 (19.8)	\$126.00	PT500X9A	Universal holder holds a variety of flasks, beakers, bottles and test tube racks. Two soft rubber bars hold vessels in place.	12 (30.5)	5.5 (13.9)	\$309.00