

Lab-Line Water Baths



Circulating Magnestir®



Model 3184

Ideal for Coliform Assays!

- Magnetic stirring circulates water evenly throughout chamber
- Temperature range from slightly above ambient to 65°C (100°C with optional gable cover)
- Temperature uniformity of $\pm 0.2^\circ\text{C}$
- ☆ 24 Month Warranty

Product Description

- Magnetic circulator is driven from below the chamber by a non-sparking induction motor. Constant speed stirring is controlled by a single switch.
- Temperature control of $\pm 0.2^\circ\text{C}$ or better, with $\pm 0.2^\circ\text{C}$ uniformity.
- Easy-to-read LED temperature display.
- Seven gallon (27 l) chamber capacity.
- Corrosion-resistant seamless stainless steel interior.
- High efficiency heating plates on sides of working chamber for maximum heat transfer.
- Three inch (8 cm) fiberglass insulation and double-wall construction increases temperature uniformity and minimizes heat loss.
- Electronic proportional temperature controller for excellent temperature stability.
- Independent hydraulic thermostat provides over-temperature protection.
- Silicone lubrication of motor assures trouble-free operation and quiet, efficient performance.
- Resettable circuit breaker protects unit from electrical overload.
- Stainless steel diffuser plate improves temperature uniformity and provides a level surface for glassware; removal of diffuser plate increases bath depth.

PRODUCT SPECIFICATIONS

Model #	Chamber Dimensions Inches (cm)			Overall Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D	
3184	14 (36)	9.75 (25)	12 (30)	16.5 (42)	17.75 (45)	18.75 (48)	40 (18)
3184-1	14 (36)	9.75 (25)	12 (30)	16.5 (42)	17.75 (45)	18.75 (48)	40 (18)

ORDERING INFORMATION

Model #	Volts	Electrical (60 Hz) Watts	Chamber Capacity Liters	US List Price
3184	120	1700	27	\$1,967.00
3184-1	240	1700	27	\$1,967.00

ACCESSORIES

Model #	Material	US List Price
19900-2	Lexan® Clear	\$101.00

318