



Blow-Hard® O.S. Extra™ Dust Remover

Fast and Easy Cleaning for Hard to Reach Areas

Ozone safe Blow-Hard® Dust Remover provides powerful blasts of triple filtered air for cleaning a wide variety of laboratory equipment which could be damaged by wiping, such as microscopes. A handy extension tube (included) effectively penetrates hard-to-reach places for thorough dust removal. Blow-Hard® O.S. Extra™ Dust Remover is defined as non-flammable (tetrafluoroethane). Warning instructions for safe use are printed on the can. *No Air Shipment.*



CATALOG NO.	SIZE	PRICE/EACH
F17080-0200	10oz can	\$29.00

Protect Valuable Lab Instruments, See Scienceware® Microscope Covers on page 55



Chromerge® Glass Cleaner

The Ultimate Glassware Cleaner

Chromerge® Glass Cleaner leaves surfaces chemically clean, protects glass and aids in creating a perfect meniscus. Each bottle contains 25ml of chromium trioxide in solution which, when mixed with a 4.1kg (9lb) 2.50 liters container of sulfuric acid, makes a highly efficient cleaner. The mixture is then storable and can be re-used until it becomes green, indicating loss of effectiveness. **Not for use with plastic.**



Packaged six 25ml bottles per box, 4 boxes per case. No Air Shipment.

CATALOG NO.	PRICE/BOX
F17089-0000	\$49.45

Cleanware™ Aqua-Clear™ Water Conditioner

Stabilize Your Water Bath

Aqua-Clear™ Water Conditioner stabilizes water in water baths and other water systems. Just 2 ml per liter of water provides long lasting protection against the growth of algae, bacteria and molds which often contaminate standing water. It also protects pumps and thermostats and helps prevent tubing from clogging. In addition, it may be added to water used in developer for color or X-ray films. It will not harm the film and permits longer use of the same solution. Aqua-Clear™ Water Conditioner adds a blue tint to the water which will fade after about 3 months reminding you to change the water and add new conditioner. 100ml polyethylene bottle with cap.



CATALOG NO.	PRICE/EACH
F17093-0000	\$49.00

Extend Life of X-ray Baths with... Cleanware™ Aqua-Clear™ Water Conditioner

NEW Pouch



Aquet® Detergent

Liquid Detergent with Neutral pH

Use Aquet® Detergent for hand washing glass and plastic labware. Used in a 1% dilute solution, this aromatic, neutral pH, phosphate-free and biodegradable polyglycol ether will not etch or mar precision glassware. It drains off glassware readily and naturally dries quickly without wiping. For convenience, Aquet® is available in individual pouches that make one gallon



of cleaner. For frequent users, Aquet® is offered in a 3.8 liter (1gallon) bottle with a 473ml (16oz) dispenser bottle. It is also available in 18.9 liter (5 gallon) and 208 liter (55 gallon) sizes.

CATALOG NO.	SIZE	CASE	BOX	PRICE/EACH
F17094-0020	20 ml pouch makes one gallon	N/A	20 pouches	\$24.95
F17094-0030	1 gallon	12	Each	31.00
F17094-0050	5 gallon Cubitainer™	6	Each	120.00
F17094-0510	55 gallon Drum	N/A	N/A	823.00

Cleanware™ Label-Off™ Label Remover

No More Sticky Residue!

Excellent for use with sample bottles, Label-Off™ Label Remover quickly removes labels, including the gummy residue other cleaners leave behind. The kit contains 59ml (2 oz.) of label adhesive solvent and 20 lint-free 102mm (4") square wipes. *No Air Shipment.*



CATALOG NO.	PRICE/EACH
F17077-0000	\$17.50

Cleanware™ Plastic Labware Cleaner

Concentrated Cleaner Makes Up to 100 Gallons of 1% Solution

Formulated Exclusively for Bel-Art by Alconox, Inc.

Cleanware™ Plastic Labware Cleaner will not cause crazing, clouding or stress cracking on plastics. Labware rinses free of interfering residues. Phosphate free, biodegradable Cleanware™ Plastic Labware Cleaner is concentrated for economical use and works in cold, warm or hot water. Also works well on a wide range of glass, ceramic and metal surfaces.



CATALOG NO.	SIZE	PRICE/EACH
F17084-0000	3.8 liter (1 gallon) bottle	\$43.80

Safely Cleans a Variety of Surfaces