



Dry-Keeper™ Desiccator Cabinet

Small Cabinets Safely Stack Five High

Polystyrene, non-vacuum cabinets are completely transparent from all sides. The door has a foam rubber gasket and two roller latches which provide an air tight seal. When open, the door slides out of the way into a compartment located underneath the cabinet. Inside dimensions are 21 x 18

x 13cm (8 1/4 x 7 x 5") with a door opening of 19 x 11cm (7 1/2 x 4 3/8"). A jar of silica gel desiccant is included.



CATALOG NO.	DIMENSIONS: LENGTH X DEPTH X HEIGHT	PRICE/EACH
H42053-0000	23 x 20 x 17cm (8 1/4 x 7 7/8 x 6 5/8")	\$114.99



Dry-Keeper™ Desiccator Cabinet with Gas Ports

Create and maintain an atmosphere at any concentration. This polystyrene desiccator cabinet is useful for tissue culture work when a controlled carbon dioxide atmosphere is required, or for flooding the chamber with nitrogen to eliminate oxygen. Two gas ports

added to the top and bottom of the H42053-0000 cabinet create a sweep of gas from the inlet to outlet. The gas ports have 4mm stopcocks, a gasketed door and two roller latches which provide an air tight seal.



CATALOG NO.	DIMENSIONS: LENGTH X DEPTH X HEIGHT	PRICE/EACH
H42053-0002	23 x 20 x 17cm (8 1/4 x 7 7/8 x 6 5/8")	\$186.99



Clear

Acrylic Desiccator Cabinets†

Dry, dust free storage for use in refrigerators, freezers or at room temperature. These durable cabinets are made of 6.35mm (1/4") acrylic. Gasketed air tight doors are 9.5mm (3/8") acrylic with full length hinges and stainless

steel twist locking latches. Choose from easy viewing clear units or dark bronze for light sensitive materials. Both afford deep, slide-in trays with perforated bottoms and are available in three sizes.



Bronze



CATALOG NO.	COLOR	SIZE	# OF TRAYS	PRICE/EACH
F42064-0000	Clear	17.78 x 17.78 x 25.4cm (7 x 7 x 10")	2	\$248.00
F42064-0001	Bronze	17.78 x 17.78 x 25.4cm (7 x 7 x 10")	2	\$277.00
F42065-0000	Clear	22.86 x 22.86 x 40.64cm (9 x 9 x 16")	2	\$316.00
F42065-0001	Bronze	22.86 x 22.86 x 40.64cm (9 x 9 x 16")	2	\$366.00
F42066-0000	Clear	30.48 x 30.48 x 30.48cm (12 x 12 x 12")	3	\$349.00
F42066-0001	Bronze	30.48 x 30.48 x 30.48cm (12 x 12 x 12")	3	\$403.00

Find "No Fuss" Desiccant Cartridges on Page 40



Mini Vacuum Desiccator

Ideal for Small Samples

The top section of this vacuum desiccator is clear, molded polycarbonate for observation from any angle. A polypropylene stopcock with a Teflon® TFE plug allows gentle air entrance when releasing vacuum. A peripheral flange around the bottom section fits into two snap locks sealing the

two halves against a Buna-N O-ring. A polypropylene stopcock with 2mm bore Teflon® TFE plug extends from the top center; the bottom section is polypropylene. Vacuum retention rating is 740mm Hg (29"). Interior dimensions are 102mm (4") diameter x 59mm (2 3/8") high. The base is 175mm (6 7/8") diameter.



CATALOG NO.	PRICE/EACH
F41099-0000	\$55.00



Vacuum Chamber and Plate

This Vacuum Chamber is useful for vacuum drying, degassing and teaching demonstrations. The chamber top is molded of clear polycarbonate for full visibility and will withstand a full vacuum at room temperature. A strong ABS plastic plate with an O-ring seals against the chamber's flange. The polypropylene stopcock with Teflon® TFE plug accepts 6.4mm (1/4") I.D. tubing. The chamber I.D.

is 20cm (7 7/8") and maximum inside height is 11cm (4 3/8"). The plate is 44mm (1 3/8") high and 232mm (9 1/8") in diameter.



CATALOG NO.	PRICE/EACH
F42043-0000	\$85.00