



Drum and Carboy Funnel

Provides Optimum Flow

Large funnels with handles are ideal for filling drums, barrels, tanks and carboys providing ample volume and fast flow. The two smaller sizes have a handle at the rim. Molded of natural polypropylene, funnels are autoclavable at 121°C (250°F).



CATALOG NO.	CAPACITY	TOP O.D.	HEIGHT	BOTTOM STEM I.D.	PRICE/EACH
H14712-0200	2.1 Liters (.56 gal)	207mm (8.15")	250mm (9.85")	21mm (.83")	\$10.59
H14712-0250	4.3 Liters (1.15 gal)	260mm (10.23")	309mm (12.5")	23mm (.90")	17.70
H14712-0350	14.1 Liters (3.78 gal)	353mm (13.9")	454mm (17.86")	26mm (1.02")	33.99



Large Industrial Drum Funnel

Safely Fill Drums and Containers

Heavy-duty funnel is molded in one piece, 6.35mm (1/4") thick. The top cylindrical section is 36cm O.D. x 13cm high (15 x 5") with a 2.54cm (1") wide top flange. The 45 degree cone is 11cm (4 1/2") high and a 14cm (5 1/2") long stem tapers from 51mm (2") to 44mm (1 3/4") O.D. Overall height is 39cm (15 3/4"); volume is 25 liters (6.7gal) or 1,550

cubic inches. *Per each, 6 per case*

CATALOG NO.	PRICE/EACH
H14716-0000	\$56.25



For Sampling from Drums and Carboys, See Hand Pumps, Page 143



Heavy Duty Funnel

Durable, Versatile, Available in Assorted Sizes

Sturdy, polypropylene funnels with ridges to prevent vapor lock. Highly recommended for continuous laboratory work. Steam autoclavable at 121°C (250°F).



CATALOG NO.	CAPACITY	STEM I.D.	STEM O.D.	STEM LGTH	TOTAL LGTH	DIA.	CASE	PRICE/BAG
F14694-0000	30ml (1oz)	7.9mm (5/16")	14.3mm (9/16")	38.1mm (1 1/2")	66.7mm (2 1/2")	54mm (2 1/8")	12 bags	\$27.00/6
F14695-0000	90ml (3oz)	6.4mm (1/4")	11.1mm (7/16")	69.9mm (2 3/4")	127.0mm (5")	76.2mm (3")	12 bags	33.00/6
F14696-0000	180ml (6oz)	12.7mm (1/2")	19.1mm (3/4")	44.5mm (1 3/4")	104.8mm (4 1/4")	101.6mm (4")	12 bags	38.00/6
F14697-0000	480ml (16oz)	11.1mm (7/16")	19.1mm (3/4")	79.7mm (3 1/4")	177.8mm (7")	152.4mm (6")	8 bags	23.50/3
F14698-0000	720ml (27oz)	19.1mm (3/4")	25.4mm (1")	69.9mm (2 3/4")	187.3mm (7 3/4")	184.2mm (7 1/4")	8 bags	25.25/3
H14699-0000	3845ml (128oz)	22.2mm (7/8")	34.9mm (1 3/8")	139.7mm (5 1/2")	342.9mm (13 1/2")	266.7mm (10 1/2")	6 each	21.00/Ea

Standard Stem Funnel

Exclusive Design Saves Time

Polypropylene funnels feature the exclusive filter-step internal design for speedier filtrations. External ribs prevent airlock. The parallel stem measures 8mm O.D. Steam autoclavable at 121°C (250°F).



CATALOG NO.	CAPACITY	TOP I.D.	TAKES FILTER PAPER	STEM LENGTH	CASE	PRICE/BAG
F14676-0055	32.8ml (1.1oz)	55mm	90mm	60mm	3 Bags	\$29.90/12
F14676-0065	59.1ml (1.97oz)	65mm	110mm	65mm	3 Bags	22.95/12
F14676-0075	91.8ml (3.1oz)	75mm	125mm	75mm	6 Bags	15.45/6
F14676-0100	209ml (7oz)	100mm	185mm	100mm	6 Bags	19.35/4

Autoclavable Funnels for Demanding Applications



Urbanti® High-Speed Filter Funnel

Internal Ribs Increase Filtration Speed

Unsurpassed for non-vacuum, high-speed filtration, these clear funnels have internal helicoid ribs to increase the filtration speed. Vertical outer ribs prevent airlock. Made from autoclavable polymethylpentene, the funnels are unaffected by corrosive liquids and are heat resistant to 177°C (350°F).



CATALOG NO.	CAPACITY	TOP I.D.	STEM I.D.	STEM O.D.	STEM LENGTH	TAKES FILTER PAPER	CASE	PRICE/BAG
F14640-0000	28.3ml (1oz)	51mm	3mm	7.5mm	150mm	75/90mm	6 Bags	\$14.20/6
F14642-0000	76.ml (2.6oz)	70mm	3mm	7.5mm	150mm	110/125mm	6 Bags	16.50/6
F14644-0000	76.ml (2.6oz)	70mm	7mm	10.0mm	70mm	110/125mm	4 Bags	16.50/6
F14646-0000	224ml (7.6oz)	100mm	7mm	10.0mm	108mm	150/185mm	4 Bags	15.45/4
F14648-0000	586ml (19.8oz)	140mm	10mm	16.0mm	132mm	240/270mm	2 Bags	40.20/3
F14650-0000	1532ml (51.8oz)	196mm	20mm	26.0mm	160mm	320/385mm	2 Bags	35.00/2