

## Heavy Duty Long Stem Funnel



60° Cone

Heavy-duty, polypropylene funnel with smooth walls, 60 degree cone angle and long stem. External ribs prevent air lock. Steam autoclavable at 121°C (250°F). Capacity 70ml (2oz). **6 per bag, 4 bags per case.**



CATALOG NO.	TOP I.D.	STEM I.D.	STEM O.D.	STEM LENGTH	TAKES FILTER PAPER	PRICE/BAG
F14671-0000	70mm	7mm	10mm	70mm	110/125mm	\$13.95



## Powder Funnel

*Specially Designed for Transferring Powders*

Polypropylene funnels with 60 degree angle. Parallel stems minimize bridging of powder. Steam autoclavable at 121°C (250°F).



CATALOG NO.	CAPACITY	TOP I.D.	STEM O.D.	STEM LENGTH	CASE	PRICE/BAG
H14660-0065	49ml (1.7oz)	65mm	15mm	25mm	3 Bags	\$24.15/12
H14660-0080	138ml (4.7oz)	80mm	18mm	25mm	3 Bags	\$33.00/12
H14660-0100	248ml (8.4oz)	100mm	21mm	25mm	4 Bags	18.00/6
H14660-0150	784ml (26oz)	150mm	29mm	25mm	6 Bags	16.20/4

## Wide Stem for Clog-Free Filling



## Powder Funnel Tapered Stem

*Specially Designed for Transferring Powders*

Polypropylene funnels with standard taper stems. Steam autoclavable at 121°C (250°F). **Per each, 12 per case.**



CATALOG NO.	CAPACITY	TOP I.D.	STEM	PRICE/EACH
F14681-0000	169ml (5.7oz)	101mm (4")	with 24/40 $\frac{1}{2}$	\$23.70
F14682-0000	594ml (20oz)	152mm (6")	with 24/40 $\frac{1}{2}$	26.55



## Large Powder Funnel

*Heavyweight High Density Polyethylene. Use for Transferring Powders.*

Economical lightweight funnel is ideal for powders and general use. Cone is 127mm (5") long with 124mm (4") top I.D. Stem is 38mm (1½") long with 38mm (1½") I.D. Capacity 510ml (17.2oz). **Per each, 6 per case.**



CATALOG NO.	PRICE/EACH
H14684-0000	\$12.39

## Utility Type Funnel

*For a Variety of Lab Applications*

Externally ribbed, polypropylene funnels are steam autoclavable at 121°C (250°F). The large capacity, 480ml (16oz) funnel is high-density polyethylene (not steam autoclavable). **12 per bag, except H14784-0000 packed 6 per bag.**



CATALOG NO.	CAPACITY	TOP DIAMETER	PRICE/BAG
H14781-0000	60ml (2.0oz)	63.5mm (2½")	\$13.95/12
H14782-0000	120ml (4.0oz)	82.6mm (3¼")	16.80/12
H14783-0000	240ml (8.11oz)	101.6mm (4")	20.40/12
*H14784-0000	480ml (16oz)	120.65mm (4¾")	18.90/6



## Micro Funnel

*When Working with Small Apertures*

Perfect for micro work, filling burettes, manometers and small bottles. These polypropylene funnels are steam autoclavable at 121°C (250°F). **12 per bag.**



CATALOG NO.	CAPACITY	TOP I.D.	STEM DIAMETER	PRICE/BAG
F14685-0024	2.6ml (.09oz)	24mm	4.5mm	\$17.00
F14685-0028	5.2ml (.17oz)	28mm	6.0mm	14.45
F14685-0035	8.1ml (.27oz)	35mm	8.0mm	13.65
F14685-0043	15.7ml (.53oz)	43mm	8.0mm	30.90

## Offset Stem Funnel

*Easy Access to Hard-to-Reach Filler Caps*

A high density polyethylene funnel with a 51mm (2") diameter, fine-mesh, bronze filter screen. The screen can be removed with a knife point and easily snapped back into place. The funnel resists most acids, alkalies and some organic solvents. It features an anti-slop inward flange and large, strong handle. The stem is offset 204mm (8") center-to-center at a 60 degree angle and tapers to 22mm (¾") O.D. 1 liter (1 quart) capacity, 127 x 152mm (5 x 6") oval top opening with a 102mm (4") deep cone.



CATALOG NO.	PRICE/EACH
H14820-0000	\$41.49



## Funnel Holders

*Supports Several Funnel Sizes*

Non-corrosive, polypropylene funnel holder and clamp includes two bushings

for 12.5mm (½") or 9.5mm (¾") diameter support rod. Outer ring supports 76 to 152mm (3 to 6") top diameter funnels. A removable inner ring holds smaller funnels from 25 to 76mm (1 to 3") top diameter. **Per each, 12 per case.**



CATALOG NO.	DESCRIPTION	PRICE/EACH
F18240-0000	Single - for one funnel	\$9.80
F18242-0000	Double - for two funnels	13.40

**For Hands-Free Operation When Preparing Multiple Samples**