



Mason Jars, Wide Mouth

Durable Versatility for Large Samples or Storage

Highly Shock Resistant polypropylene, heavy-walled Mason Jars are ideal for industrial waste sampling, long term storage and a variety of laboratory uses. They have standard 70mm, G-threaded necks, and screw caps with unique grip lugs and removable leak proof silicone rubber ring seal. Cap and jar are steam autoclavable at 121°C (250°F). Loosen cap when autoclaving. **Bags and case listings below.**



CATALOG NO.	CAPACITY LITERS US	DIAMETER	HEIGHT	CASE	PRICE/BAG
F10913-0000	½ (1 pt)	81mm (3.19")	135mm (5.3")	4 Bags	\$31.00/6
F10914-0000	1 (1 qt)	100mm (3.93")	178mm (7.0")	4 Bags	47.00/6
F10915-0000	2 (½ gal)	125mm (4.92")	213mm (8.38")	4 Bags	35.25/3
F10916-0000	4 (1 gal)	149mm (5.87")	286mm (11.25")	3 Bags	39.25/2
F10917-0000	8 (2 gal)	188mm (7.40")	334mm (13.13")	6 Each	30.00/Ea

Mason Jar Cap

70mm G screw caps fit standard Mason Jars. Grip lugs and removable leakproof seal ring. **12 per bag, 12 dozen per case.**



CATALOG NO.	PRICE/BAG
F12640-0000	\$27.85



1 Gallon Safety Jar with Lid

Seamless, heavy-walled, polyethylene jar is chemically resistant with friction-fit snap-on lid. Capacity of jar is 4 liters (1 gallon). Jar has an inner diameter of 152mm (6") and is 305mm (12") high.



CATALOG NO.	PRICE/EACH
H17205-0000	\$40.80

Write On™ Labeling Tape Is Ideal for Identifying Jar Contents See Page 69



White Polypropylene Jars with Screw Caps

White polypropylene jars and screw caps in a wide assortment of sizes provide excellent storage for liquid or dry stock, or for use with samples. Note: Jars are autoclavable at 121°C (250°F) Screw caps, with PE foam lining are not autoclavable. **12 per bag.**



CATALOG NO.	CAPACITY	DIAMETER	INSIDE DIAMETER	HEIGHT	PRICE/BAG
H17906-0000	14.8ml (½oz)	33mm (1.3")	27.0mm (1.06")	37.0mm	\$7.20
H17908-0000	29.6ml (1oz)	43mm (1.69")	35.5mm (1.40")	37.3mm	9.30
H17910-0000	118.0ml (4oz)	58mm (2.28")	51.0mm (2.00")	66.5mm	13.95
H17911-0000	236.6ml (8oz)	70mm (2.76")	62.0mm (2.44")	86.0mm	21.00
H17912-0000	474.0ml (16oz)	89mm (3.50")	82.0mm (3.23")	101.0mm	30.00



Vial with Cap

This clear 192 ml (6.5 oz) polystyrene vial with white polypropylene screw cap is ideal for small samples. The vial's leak-proof lid has an integral seal ring. It also threads perfectly onto the Scienceware® Samplit™ Scoop. Dimensions: 47 mm I.D. (1.85") x 110 mm (4.36") high. Offered in a box of 25 vials with adhesive labels for sample identification. **25 per box.**



CATALOG NO.	PRICE/BOX
H17590-0000	\$58.65

Sample Vials with Captive Closure

Cylindrical, low-density polypropylene vials with friction-fit, captive closures open and close with a flick of your thumb. Suitable for many applications. **12 per bag, 12 dozen per case.**



CATALOG NO.	CAPACITY	O.D.	HEIGHT	PRICE/BAG
F17570-0000	0.273ml	7.0mm	13.7mm	\$11.00
F17571-0000	1.48ml	11.6mm	24.7mm	11.85
F17572-0000	1.85ml	11.6mm	30.6mm	13.00
F17573-0000	5.00ml	15.00mm	50.0mm	14.20
F17574-0000	7.39ml	16.0mm	56.7mm	15.20

Great for Storing Small Items

Chemical Containers

Useful as sample or weighing bottles, or for storing liquids or solids. Heavy-walled, injection molded, polyethylene jars are easy to clean. Screw cap design features a leakproof inner seal. **6 per bag. Cases listed below.**



CATALOG NO.	CAPACITY	O.D.	HEIGHT	CASE	PRICE/BAG
H17870-0000	5cc	22mm	33mm	6 Bags	\$20.25
H17871-0000	10cc	22mm	54mm	6 Bags	21.00
H17872-0000	30cc	32mm	65mm	6 Bags	27.60
H17874-0000	90cc	52mm	69mm	6 Bags	30.90
H17875-0000	120cc	52mm	84mm	6 Bags	30.90
H17876-0000	180cc	52mm	114mm	4 Bags	32.49

Pipette Jars

For Cleaning and Storage

The leakproof, molded polyethylene Pipette Jars store and soak pipettes and will withstand all common cleaning solutions at room temperature, including sulfuric acid/dichromate. 152mm (6") in diameter. **Per each, 4 per case.**



CATALOG NO.	HEIGHT	PRICE/EACH
F17130-0018	46cm (18")	\$61.80
F17130-0024	61cm (24")	64.90



For Additional Pipette Cleaning Jars, See Page 60

Hydrometer Jar

Clear, unbreakable polycarbonate cylinder, 33mm I.D. x 305mm high (1 ½" x 12") with removable polyethylene base. **Hydrometer Not Included. Per each, 6 per case.**



CATALOG NO.	PRICE/EACH
H17817-0000	\$31.89

