

Double Stack™ Disposable Tube Holders & Needle Holder

Increase Tray Capacity and Improve Organization

Unique die-cut foam blocks accommodate up to 50% more tubes than conventional racks and fit inside the modular compartments of your tray. The well pattern in these blocks positions tubes with spacing and depth so that additional similar sized tubes can be securely held between tubes placed in wells. For easiest cleaning, simply remove old inserts, clean the inside of your tray and pop in fresh, clean DoubleStack™ inserts. **5 per bag.**




Large Tube Holders:

| CATALOG NO. | DESCRIPTION | TUBE SIZE | HOLDS | PRICE/BAG |
|-------------|---------------------------------------|-----------|-------|-----------|
| F18634-0001 | 9 L x 9 W x 5.1cm H (3½ x 3½ x 2) | 16mm | 16 | \$8.25 |
| F18634-0002 | 9 L x 9 W x 3.8cm H (3½ x 3½ x 1½) | 13mm | 25 | 8.25 |
| F18634-0003 | 9 L x 9 W x 3.8cm H (3½ x 3½ x 1½) | 10mm | 25 | 8.25 |

Small Tube Holders:

| CATALOG NO. | DESCRIPTION | TUBE SIZE | HOLDS | PRICE/BAG |
|-------------|--|-----------|-------|-----------|
| F18635-0001 | 9 L x 7 W x 3.8cm H (3½ x 2½ x 1½") | 16mm, | 12 | \$9.00 |
| F18635-0002 | 9 L x 7 W x 3.8cm H (3½ x 2½ x 1½") | 13mm | 20 | 9.00 |
| F18635-0003 | 9 L x 7 W x 3.8cm H (3½ x 2½ x 1½") | 10mm | 20 | 9.00 |
| F18635-0004 | Needle Holder | | 48 | 8.25 |

For Additional Racks, See Pages 103-117



The Collector™ Blood Tray

Everything You Need to Carry from Lab to Patient

Designed for flexibility and versatility, The Collector™ Blood Tray can be rearranged to meet your specific needs. It includes two 40-place test tube supports (for

tubes up to 16mm O.D.), one 18-place slide holder, a large requisition section, and adjustable partitions to hold needles, syringes, cotton and other materials. Sturdy polypropylene construction holds its shape even when fully loaded. Steam autoclavable at 121°C (250°F). Dimensions: 356 x 292 x 51mm H (14 x 11½ x 2"H). The handle is formed by an end mounted, detachable loop that is 178mm (7") high. **Per each, 4 per case.**



| CATALOG NO. | DESCRIPTION | PRICE/EACH |
|----------------|--|------------|
| A. F18631-0715 | Complete Blood Tray Set | \$99.00 |
| B. F18631-0716 | Replacement test tube support. 167 x 109 x 109H mm (6½ x 4¼ x 4¼") 40 places for tubes up to 16mm O.D. | 17.00 |


Custom Organizer Tray

Mix and Match Partitions to Make the Tray Organizer You Need

This rugged injection molded ABS organizer tray is made to hold both long and short items, whichever suits your needs. There are 5 extra partitions which can be used to create custom compartments of various sizes. The tray is 318 x 356mm (12½ x 14") by 51mm (2") deep. The center compartment is 356mm (14") long, with other compartments 178mm (7") long. ABS Plastic. **Per each, 6 per case.**



| CATALOG NO. | PRICE/EACH |
|-------------|------------|
| F18667-0000 | \$21.65 |



A. Lab Fridge™ Tray Racks

Eliminates Clutter, Provides Inventory Refrigeration and Allows Stock Rotation

A highly valuable aid for refrigerated storage of small bottles of blood-typing sera, chemistry controls, enzyme reagents, coagulant reagents, etc. The racks are also useful in drawers for separating pipettes and stirring rods. They provide the most efficient storage when used with the Scienceware® H18663 Lab Fridge™ Tray Cabinet. Made from vacuum-formed white styrene. Overall dimensions, 330 x 330 x 44mm H (13 x 13 x 1¾"H). **Per each, 6 per case.**



| CATALOG NO. | CHANNELS | CHANNEL WIDTH | PRICE/EACH |
|-------------|----------|---------------|------------|
| H18661-0000 | 9 | 25.4mm (1") | \$33.99 |
| H18662-0000 | 8 | 32.0mm (1¼") | 33.99 |

B. Lab Fridge™ Tray Cabinet

Organize Small Bottles in the Lab

This open-front cabinet holds 3 Lab Fridge™ Tray Racks of either style. It will accommodate bottles with a maximum height of 83mm (3¼"). The cabinet is polypropylene with styrene rails to hold trays. Overall dimensions, 36 x 34 x 39cm H (14 x 13½ x 15½"H). (Trays not included.)



| CATALOG NO. | PRICE/EACH |
|-------------|------------|
| H18663-0000 | \$199.50 |

For Space-Saving Storage in the Refrigerator