

**NEW**

### Pipette-EZE™ Pipettor

*Transfer Liquids with Ease*

The Pipette-EZE™ Pipettor combines convenience with ease-of-use in a light weight, comfortable design. Power is continuously supplied by its 120 VAC adapter, eliminating the need for periodic battery charging. The ergonomic design incorporates "trigger-style" pushbutton controls for precise filling and dispensing. The pipettor works with any standard glass or plastic pipettes from 0.5ml to 100ml with a mouth O.D. from 3 to 9.9mm. A

self-contained, replaceable hydrophobic filter prevents internal aspiration and cross-contamination. The pipettor's nosepiece is removable and fully autoclavable. Each unit comes complete with pipettor, hose, vacuum pump and power adapter.

CATALOG NO.	VOLTAGE	PRICE/EACH
F37907-1000	120 VAC	\$265.50

**No Batteries Needed - Runs on Electricity.  
Always Ready to Use - No Charging Required!**



### Max® Electronic Memory Pipettor

US Patent No. 5,509,318

*Built-in Memory Efficiently  
Repeats Your Volumes*

Max® Memory Pipettor is the first portable pipetting aid with a built-in memory. Max® Memory Pipettor has 3 modes of operation:

- 1) Auto stop during aspiration or dispensing
- 2) Step pipetting during dispensing
- 3) Manual pipetting
- 4) Accuracy 5% below, 2% above 2ml

To memorize a volume, simply place in "SET" mode, aspirate or dispense the volume up to 25ml, then switch to "RUN". Max® Pipettor will aspirate and stop at the desired volume automatically or function as a step pipettor, repeatedly dispensing the set volume until the pipette is empty. A thumb wheel dial allows adjustment of the speed, even while pipetting. Max® Memory Pipettor has a unique venting system that protects the unit while pipetting corrosives. An auto sleep function prolongs the battery life by turning the power off when the unit is idle. One charge of the AA NiCad batteries gives a full day's pipetting. Batteries recharge overnight. Max® Memory Pipettor accepts all standard glass and plastic pipettes, from 1ml to 75ml with a mouth O.D. from 3 to 9.9mm. The pipette holder can be autoclaved. A self-contained, replaceable hydrophobic 0.2 micron filter prevents overfilling. Trigger style buttons control filling and dispensing. Max® Pipettor comes complete with a wall battery charger.



CATALOG NO.	DESCRIPTION	PRICE/EACH
F37910-0000	Max® Memory Electronic Pipettor 120 VAC	\$465.00
F37910-0001	Max® Memory Electronic Pipettor 230 VAC	499.50
H37908-0250	Replacement Filters	10.65/2

### Mopet® Portable Pipettor

*Fast and Portable*

The Mopet® Pipettor will pick up or dispense 10ml in three seconds. It takes only ten seconds to fill a 100ml pipette. This powerful fixed speed pipettor is ideal for large volume and viscous liquid pipetting. The Mopet® Pipettor accepts all standard glass and plastic pipettes from 1.0ml to 100ml (mouth O.D. from 3 to 9.9mm). This NiCad battery-powered Mopet® Pipettor aspirates and dispenses liquids with

push-button controls. Built-in replaceable hydrophobic filter prevents liquid from being accidentally drawn into the pipettor and protects against contamination.

The Mopet® Pipettor can be completely recharged in 5½ hours, and a fully charged unit will perform a full day's work.

Available with a plug-in wall charger and with Bel-Art's exclusive charger stand that safely stores the Mopet® Pipettor while charging. Speed not adjustable.



CATALOG NO.	DESCRIPTION	PRICE/EACH
F37908-0000	Mopet® Pipettor with wall charger, 120 VAC	\$297.50
F37908-0001	Mopet® Pipettor with wall charger, 230 VAC	325.00
F37909-0000	Mopet® Pipettor with charger stand, 120 VAC	480.00
F37909-0001	Mopet® Pipettor with charger stand, 230 VAC	480.00
F37908-0100	Charger stand only, 120 VAC	235.00
F37908-0102	Charger stand only, 230 VAC	235.00
H37908-0250	Replacement filters	10.65/2

### Portapet® Air Powered Portable Pipettor

*Completely  
Self-Contained  
Power Source*

The Portapet® Pipettor combines fine control and portability with simplicity and light weight. At only 198gms (7oz), the power source is a low density polyethylene bellows. Just compress the bellows to generate more than 100ml of total pipetting capacity. Portapet® Pipettor is ideal for pipettes from 0.5 to 20ml. By simply removing the bellows, the Portapet®

Pipettor can be attached to a vacuum source for unlimited pipetting power. A self-contained, replaceable hydrophobic filter prevents overfilling. Trigger style buttons control filling and dispensing. Portapet® Pipettor works with all standard glass and plastic pipettes, from 0.5ml to 100ml (mouth O.D. from 3 to 9.9mm).



CATALOG NO.	DESCRIPTION	PRICE/EACH
F37907-0000	Portapet® Pipettor	\$189.00
H37908-0250	Replacement Filters	10.65/2