




No-Wire™ Autoclavable Vial and Bottle Racks

Sizes to Fit Most Small Vials and Bottles

No-Wire™ vial and bottle racks are small enough to hold in one hand and diverse enough to accommodate most available vials and small bottles. Polypropylene rack is available in 5 sizes. Dimensions 128 x 105 x 43mm high (5.08 x 4.15 x 1.70"). **Per each, 8 per case.**



CATALOG NO.	HOLE SIZE	PLACES	ROWS	FITS	PRICE/EACH
F18514-0013	13.38mm	6	6 x 7	13mm vials, bottles	\$12.00
F18514-0016	16.38mm	12	3 x 4	16mm vials, bottles	12.00
F18514-0020	20.37mm	20	4 x 5	20mm vials, bottles	12.00
F18514-0025	25.40mm	30	5 x 6	25mm vials, bottles	12.00
F18514-0030	30.48mm	42	2 x 3	30mm vials, bottles	12.00




In-Line Test Tube Rack

For high-visibility, this 6-place high-density polyethylene rack is end-supported by two columns. The top plate has 25mm diameter holes and the base has corresponding hemispherical wells. Includes six vertical pins for tube drying. Parts snap together for easy assembly. 241 x 76 x 92H mm (9 1/2 x 3 x 3 3/4"). **Per each, 12 per case.**



CATALOG NO.	PRICE/EACH
F18527-0000	\$7.25




Economy Test Tube Rack

Three Levels Offer More Support to Each Tube

Economical, unbreakable and chemical-resistant, this rack has 3 tiers to support tubes. Made from polyethylene, the rack holds 12 tubes up to 19mm in diameter. 194 x 60 x 83mm high (7 3/4 x 2 1/4 x 3 1/4").



CATALOG NO.	PRICE/EACH
F18532-0000	\$14.70



Full-View Test Tube Supports

Hold More Tubes, Waste Less Space

A versatile polypropylene rack with pins for support or draining. Tubes can be stored upright, full or empty, or inverted over pins for draining. Steam autoclavable at 121°C (250°F). **Per each, case listed below**



CATALOG NO.	HOLDS TUBES	PLACES	ROWS	PINS	DIMENSIONS	CASE	PRICE EACH
H18910-0102	10-13mm	80	5 x 16	102	23 x 8.5 x 6.3cm	24	\$15.99
H18910-0189	10-13mm	160	8 x 20	189	29.5 x 13 x 6.4cm	24	32.25
H18911-0000	14-17mm	50	5 x 10	66	18 x 10 x 6.7cm	12	17.70
H18912-0000	17-20mm	40	4 x 10	55	23 x 10.5 x 7cm	12	17.70
H18914-0000	22-25mm	30	3 x 10	44	28 x 10 x 7.5cm	12	17.70
H18915-0000	25-30mm	24	3 x 8	36	29 x 12.7 x 7.5cm	12	20.49

Poxygrid® Hybridization Bottle Rack

Protect and Store Hybridizing Bottles


Heavy-duty rack protects your expensive bottles within individual cells, which keep them from rolling, tipping and striking against one another. The rack holds 150mm and 300mm long bottles and is designed to leave ample space between the bottles to allow for the screw cap diameter and to provide convenient handling. Made of welded steel rods, this rack is epoxy-coated for maximum protection, easy cleaning and is steam autoclavable at 121°C (250°F). Rack is 130mm (5 1/8") high, holes are 42mm square for 40mm diameter bottles.



CATALOG NO.	PLACE	DIMENSIONS	PRICE/EACH
F16964-0006	6	137 x 198mm (5 3/8 x 7 7/8")	\$48.00

Nessler Tube Rack

50ml Nessler tube rack has molded white polypropylene decks supplied in 5 pieces for easy snap together assembly. The flat lower deck creates a reflective surface for the tubes. This rack has 6 places and is 30cm (11 3/4") long, 9.5cm (3 3/4") wide and 14cm (5 1/2") high. Top to bottom deck dimension 10 cm (4"). 27mm hole size. Steam autoclavable at 121°C (250° F). **Per each, 6 per case.**




CATALOG NO.	PRICE/EACH
F18986-0050	\$39.95

Serum Vial Racks

Designed to Withstand Elevated Water Temperatures for Incubation

Just the right size for serum vials (10-13mm diameter), RIA vials and T3/T4 reaction vials (13-16 & 16-20mm diameter). Long-lasting and corrosion-resistant, these polypropylene racks fit most water baths for elevated temperature incubations. Top shelf height is 41mm (1 5/8"). End plates have handles. Steam autoclavable at 121° C (250° F)



CATALOG NO.	LENGTH X WIDTH X HEIGHT	HOLDS TUBES	PLACE	PRICE/EACH
F18513-0072	238 x 111 x 89mm (9 3/8 x 4 3/8 x 9 1/2")	10-13mm	72	\$18.00
F18513-1640	254 x 111 x 89mm (10 x 4 3/8 x 9 1/2")	13-16mm	40	18.00
F18513-2040	254 x 111 x 89mm (10 x 4 3/8 x 9 1/2")	16-20mm	40	18.00

Scintillation Vial Rack

For Start to Finish Sample Preparation

This rugged and stable rack holds 24 scintillation vials (30mm diameter) for convenient storage, handling, and sample preparation. Made of corrosion-resistant, autoclavable polypropylene. Dimensions 298 x 127 x 95mm high (11 3/4 x 5 x 3 3/4"). Shelf height 44mm (1 3/4"). **Per each, 4 per case.**



CATALOG NO.	PRICE/EACH
F18512-0000	\$24.75