



**Safety First!**  
Always remember to wear your safety goggles.



### Safety Goggles

*Fits Comfortably over Prescription Glasses*

Molded of soft, clear, flexible vinyl that conforms to the face and nose, these goggles have 4 vents to reduce fogging. With clear acetate lenses, the goggles can be worn over regular eye-glasses with complete comfort. 13mm

( $\frac{1}{2}$ ") wide adjustable elastic headband. **Per each, 6 per case.**

CATALOG NO.	PRICE/EACH
H24873-0000	\$12.00



### Safety Eyewear Holders\*

*Perfect Between Office and Lab*

Designed to be wall mounted or free standing, these clear acrylic eyewear holders are available with or without a lid. 95.3 x 22.9 x 15.2cm (3 $\frac{3}{4}$  x 9 x 6"). There are pre-drilled holes on the back for wall mounting.



CATALOG NO.	DESCRIPTION	PRICE/EACH
F24877-0000	Holder with lid	\$53.00
F24877-0001	Holder without lid	40.00



### Eyewear Dispenser\*

*Keep Safety Glasses within Reach!*

This wall mounted dispenser makes eyewear readily accessible when needed. Designed to hold either boxed or loose glasses, the dispenser features a large flip-up door for removing glasses and open top for easy filling. Made from clear acrylic, you'll always know when it's time to refill the dispenser. Overall dimensions: 20.3 x 10.2 x 48.3 cm H(8 x 4 x 19"H).



CATALOG NO.	PRICE/EACH
F24876-0000	\$97.75



### Cleanware™ Anti-Fog Lens Cleaner

*Economical, Only One Drop Does the Trick*

Carry Liquid Lens Cleaner in your pocket or purse ready to clean glass and plastic lenses anytime, anywhere. Only one drop on each lens cleans and provides anti-static as well as anti-fogging protection. Supplied in a convenient 59ml (2 oz) cylinder bottle with reclosable dispenser cap, the fluid can be dispensed in a single drop or stream. The polyethylene bottle is 35mm (1 $\frac{1}{8}$ "

diameter x 100mm (4") tall. **2 per bag, 12 bags per case. No Air Shipment.**

CATALOG NO.	PRICE/BAG
F24842-0001	\$11.50

### Rayon Cloth Tissues

*Keep a Pop-Up Dispenser on the Bench*

For use on glass and plastic lenses, these tissues do not contain silicone. 150 sheets in a pop-up dispenser pack. Each sheet is 108 x 127mm (4 $\frac{1}{4}$  x 5"). **Per pack, 10 packs per case.**



CATALOG NO.	PRICE/EACH
F24841-0000	\$9.50

### Disposable Lens Cleaning Station

*All-In-One Station, Tissues and Cleaning Solution*

For economy, convenience and efficiency, this self-contained, disposable station combines cleaning solution and lens tissues. Fabricated of coated paperboard, the 14 x 21 x 21cm high (5 $\frac{1}{2}$  x 8 $\frac{1}{4}$  x 8 $\frac{1}{4}$ ") station has openings for wall mounting. The station contains a generous supply of 600 lens tissues 13 x 20cm (5 x 8") and a 237ml (8 oz) polyethylene spray pump bottle filled with non-silicone fluid that protects, polishes, anti-stats and anti-fogs your lenses. **Per each, 8 per case. No Air Shipment.**



CATALOG NO.	PRICE/EACH
F24835-0000	\$22.00

### Lab Lens® Liquid

Lens cleaning solution in a 473ml (1 pint) polyethylene bottle. One spray pump is included in each case of 12 bottles purchased. **Per each, 12 per case.**



CATALOG NO.	PRICE/EACH
F24833-0000	\$8.80



### Draft Checker™ Air Flow Detection Powder

*Ideal Fume Hood Accessory*

A quick, easy way to test the safety of fume hood exhaust before use. Just a squeeze or two from the plastic bottle dispenses a cloud of ultra-fine, lightweight silica powder, which remains suspended long enough to visually check whether or not fume hood exhaust is functioning properly. It is also useful for detecting drafts, gas leaks and air currents in rooms, ducts and exhaust systems. The amorphous silica powder is inert, non-toxic and non-flammable. 1.5 gm bottle.

CATALOG NO.	PRICE/EACH
F19939-0000	\$13.65