



CATALOG NO.	LABEL	SIZE	PRICE/BAG
F11646-0622	Acetone	500ml (16oz)	\$22.50
F11646-6250	Dichloromethane	500ml (16oz)	27.75
F11646-0631	Deionized Water	500ml (16oz)	22.50
F11646-0620	Distilled Water	500ml (16oz)	22.50
F11646-0637	Ethanol	500ml (16oz)	22.50
F11646-0619	70% Ethanol	500ml (16oz)	22.50
F11646-0612	Ethyl Acetate	500ml (16oz)	22.50
F11646-0624	Isopropanol	500ml (16oz)	22.50
F11646-0638	LYOB (Label Your Own Bottle)	500ml (16oz)	22.50
F11646-0459	Machine Oil	500ml (16oz)	22.50
F11646-0623	Methanol	500ml (16oz)	22.50
F11646-0611	Methyl Ethyl Ketone (MEK)	500ml (16oz)	22.50

Safety Labeled Wash Bottles, Wide Mouth, Color-Coded

Safety Favorites!

Instant Recognition! Informative Labeling! Added Convenience and Safety!

Bel-Art's Safety Labeled wash bottles ensure leakproof dispensing of your most common solvents. Polypropylene closure is 53mm wide with a leakproof design. Easy to read, LDPE solvent resistant printed bottles include: Department of Transportation (DOT) and National Fire Protection Association (NFPA) codes and symbols, Chemical Abstract Service (CAS) numbers, protective equipment symbols and target organ information, making your workplace a safer place. **6 per bag, 2 bags per case.**

LDPE PP

CATALOG NO.	LABEL	SIZE	PRICE/BAG
F11646-0617	Saline Solution	500ml (16oz)	\$22.50
F11646-0614	Soap	500ml (16oz)	22.50
F11646-0627	Sodium Hypochlorite (Bleach)	500ml (16oz)	22.50
F11646-0628	Toluene	500ml (16oz)	22.50
F11646-0621	Water	500ml (16oz)	22.50
F11646-2232	Acetone	1000ml (32oz)	29.00
F11646-6032	Distilled Water	1000ml (32oz)	29.00
F11646-3732	Ethanol	1000ml (32oz)	29.00
F11646-2432	Isopropanol	1000ml (32oz)	29.00
F11646-3832	LYOB (Label Your Own Bottle)	1000ml (32oz)	29.00
F11646-1132	Methyl Ethyl Ketone (MEK)	1000ml (32oz)	29.00

Custom labeling is available LDPE wash bottles are for dispensing use only, not for long term storage.



Safety-Vented™ Labeled Wash Bottles

Patent No. 5,125,543

Eliminate Dangerous Drips

A unique, patented, 2-piece valve prevents pressure build-up of solvents within the wash bottle. The volatiles are vented eliminating solvent drips. Patented valve allows dispensing with bottle upright or inverted. Wash bottles are labeled with six of the most commonly used solvent names. Bottles are low-density polyethylene, closures and venting assembly are polypropylene. Note: Use only with solvents compatible with LDPE. **3 per bag.**

LDPE PP

CATALOG NO.	LABEL	SIZE	PRICE/BAG
F11642-0622	Acetone,	500ml (16oz)	\$17.00
F11642-0637	Ethanol,	500ml (16oz)	17.00
F11642-0612	Ethyl Acetate,	500ml (16oz)	17.00
F11642-0624	Isopropanol,	500ml (16oz)	17.00
F11642-0638	LYOB (Label Your Own Bottle),	500ml (16oz)	17.00
F11642-0623	Methanol,	500ml (16oz)	17.00
F11642-0611	Methyl Ethyl Ketone (MEK),	500ml (16oz)	17.00
F11643-0222	NEW Acetone,	250ml (8oz)	12.50
F11643-0237	NEW Ethanol,	250ml (8oz)	12.50
F11643-0224	NEW Isopropanol,	250ml (8oz)	12.50
F11643-0238	NEW LYOB (Label Your Own Bottle),	250ml (8oz)	12.50
F11643-0223	NEW Methanol,	250ml (8oz)	12.50

LDPE wash bottles are for dispensing use only, not for long term storage.




Safety Labeled Wash Bottles, Variety Six-Pack

Wide Mouth, Color Coded

Consists of one 500ml (16oz) bottle each of: Acetone, Ethanol, Isopropanol, Methanol, Toluene, and Water. **6 per Bag, 2 bags per case.**

CATALOG NO.	PRICE/BAG
F11646-0000	\$22.50



For Bottle Carriers See Pages 11 & 15