



## Scienceware® Economy Sampling Scoops

These new Economy Sampling Scoops are an inexpensive and accurate measure for smaller quantities of dry or liquid measures. Accurate capacity is listed on each scoop and leveling the sample to the flat top edge ensures a repeatable and consistent volume measure. The flat bottom shape and light weight allows for use as a weighing boat. These one-piece scoops are molded of durable white high impact polystyrene. Select 63, 83, 125 or 250ml (2oz, 2.7oz, 4oz or 8oz). **12 per bag.**

CATALOG NO.	CAPACITY	PRICE/BAG
F36933-0021	63ml (2oz)	\$6.00
F36933-0022	83ml (2.7oz)	7.00
F36933-0023	125ml (4oz)	8.00
F36933-0024	250ml (8oz)	11.00



## Assorted Polypropylene Scoops

Six Sizes with Special Bulk Quantity Pricing Available

Molded polypropylene; rigid, strong and smooth. The smaller 60 and 125cc scoops have been improved with a

longer handle that increases the distance between the hand and sample. These scoops have also been redesigned to remain upright when filled or empty allowing them to be used as weighing boats or to be placed on a flat surface. Steam autoclavable at 121°C (250°F). Special prices on yearly contracts for large quantities and colors are available on request.

CATALOG NO.	BOWL SIZE	HANDLE	CAPACITY	CASE	PRICE/BAG
F36747-0000	60 x 29 x 25 (2½ x 1½ x 1")	32 (1¼")	22cc	6 Dozen	\$7.50/12
F36748-0000	76 x 38 x 25 (3 x 1½ x 1")	150 (6")	60cc	6 Dozen	18.80/12
F36749-0000	89 x 38 x 32 (3½ x 1½ x 1¼")	48 (1¾")	75cc	6 Dozen	18.30/12
F36750-0000	92 x 48 x 35 (3¾ x 1¾ x 1¾")	150 (5")	125cc	4 Dozen	18.30/12
F36753-0000	127 x 51 x 51 (5 x 2 x 2")	64 (2½")	200cc	3 Dozen	18.80/12
F36756-0000	203 x 121 x 76 (8 x 4¾ x 3")	102 (4")	1100cc	4 Bags	21.90/6



## Fluo-Kem® Teflon® FEP Coated Aluminum Scoops

FDA Compliant for Food Contact

Coating is anti-stick, corrosion resistant brown Fluo-Kem® Teflon® containing FEP, which complies with FDA regulations for food contact. Handle is not coated so your grip is non-slip.

CATALOG NO.	BOWL LENGTH	CAPACITY liters quarts	DIAMETER	HANDLE LENGTH	PRICE/EACH
H36733-0000	184mm (7¼")	1 1	102mm (4")	83mm (3¼")	\$49.95
H36735-0000	216mm (8½")	2 2	127mm (5")	89mm (3½")	57.30
H36737-0000	260mm (10¼")	4 4	146mm (5¾")	102mm (4")	66.99

## Large Scoop

Large Flat Edge for Easy Scooping

Rectangular shape molded of rigid polypropylene. Excellent for scooping ice cubes. Holds approximately 2 liters (½ gallon) with rounded corners for easy cleaning and a looped handle for hanging. Bowl measures 140 x 229 x 89mm (5½ x 9 x 3½") with a 114mm (4½") handle.



CATALOG NO.  
F36758-0000

PP A  
PRICE/BAG  
\$62.50

## Large Fabricated Scoops

Practical Non-Spill Design in Three Useful Sizes

Large, heavy duty, polyethylene scoops are strong, durable and excellent for scooping dry powders, granules or crystals. They have a closed back, non-spill shape and an ergonomically angled handle.



CATALOG NO.	DESCRIPTION	PRICE/EACH
H36758-0001	Holds approximately 1 liter (1 quart). 125mm diameter (4¾") x 127mm (5") long scoop with 127mm (5") long handle 25mm (1") diameter.	\$24.75
H36759-0000	Holds approximately 2 liters (½ gallon). 152mm diameter (6") x 229mm (9") long scoop with 127mm (5") long handle 32mm (1¼") diameter.	32.25
H36760-0000	Holds approximately 4 liter (1 gallon). 191mm diameter (7½") x 248mm (9¾") long scoop with 127mm (5") long handle 32mm (1¼") diameter.	42.30

## Scienceware® Samplit™ Scoop and Container System

This new collection device incorporates a sturdy sampling scoop that threads into a dedicated leak proof 190ml (6.5oz) container ensuring safe and efficient isolation of material. The rigid scoop is molded of white, high impact polystyrene strong enough for coring semi solids such as cheese, ice cream and even soil or granular materials. The heavily ribbed handle provides an excellent grip for pushing and turning the scoop. The depth and shape of the scoop keeps samples contained while allowing for easy dispensing into sample jars or bags. Vials are crystal clear styrene for easy viewing of contents and are supplied with sample identification labels. Dimensions: 47 mm I.D. (1.85") x 110 mm (4.36") high.



CATALOG NO.	PRICE/BOX
H36915-1000	Samplit™ Scoop and Container Systems \$102.00
H37835-0000	Samplit™ Scoops only 49.20
H17590-0000	Samplit™ Vials with caps only 58.65