

NEW



Scienceware® Mini Sampler Spoons

Disposable Mini Measures

Ideal for retrieving small, specific volume samples of powders or similar materials. Use a separate spoon for each material sampled to avoid cross contamination. Available in 5 sizes. Molded of

white high impact polystyrene. **25 per bag, 10 bags per case.**



CATALOG NO.	DESCRIPTION	PRICE/BAG
F36721-0005	0.05ml (0.0016 oz)	\$8.95
F36721-0010	0.10ml (0.0034 oz)	8.95
F36721-0025	0.25ml (0.0085 oz)	8.95
F36721-0050	0.50ml (0.017 oz)	8.95
F36721-0100	1.00ml (0.034 oz)	8.95

Scienceware® Economy Sampling Spoons



Economical small measures

These economical measures are ideal for retrieving small, specific volume samples of powders or liquids. Their unique flat bottom shape allows them to be used as

weighing boats for additional sampling accuracy. Use a separate spoon for each material sampled to avoid cross contamination. All sizes are molded of white, high impact polystyrene for durability and stiffness when sampling packed powders. **12 per bag.**



CATALOG NO.	CAPACITY	PRICE/BAG
F36740-0021	1.25ml (.04oz)	\$4.80
F36740-0022	2.5ml (.08oz)	4.80
F36740-0023	5.0ml (.17oz)	4.80
F36740-0024	8.0ml (.27oz)	4.80
F36740-0025	20.0ml (.67oz)	4.80



Sampler Spoons

Choose the Bowl Size for Your Needs

These long-handled spoons are convenient for reaching into jars and bottles. Molded of talc-filled polypropylene, they are steam autoclavable at 121°C (250°F). All are 18cm (7") long. **12 per bag, 12 bags per case.**



CATALOG NO.	CAPACITY	WIDTH	PRICE/BAG
F36723-0000	1.23ml (¼ teaspoon)	20mm	\$10.60
F36724-0000	2.46ml (½ teaspoon)	25mm	11.00
F36725-0000	4.93ml (1 teaspoon)	30mm	11.85
F36726-0000	14.79ml (3 teaspoon)	39mm	11.85
F36727-0000	Assortment of 3 of each size		11.60



Lab Spoon

Extra Length for Sampling from Larger Bottles

This large stainless steel lab spoon is 305mm (12") long for reaching into larger bottles. One end is a flat spatula measuring 31mm (1¼") wide; the spoon end is 37mm (1½") wide and holds 5ml. The connecting bar is 14mm (½") wide. **Per each, 12 per case.**



CATALOG NO.	PRICE/EACH
F36711-0012	\$34.55

Stainless Steel Spoons

Dual Bowl Sizes on a Sturdy Handle

Double-ended, stainless steel spoons with a solid, rounded shaft. Available in plain polished stainless steel or with ends

coated with Fluo-Kem® Teflon® FEP which complies with FDA regulations governing direct food contact applications. **2 per pack.**



All Stainless Steel

CATALOG NO.	LENGTH	SMALL SPOON	LARGE SPOON	PRICE/PACK
H36729-0015	15cm	0.5ml	2.0ml	\$28.50
H36729-0018	18cm	0.8ml	2.8ml	31.89
H36729-0021	21cm	1.6ml	5.0ml	36.99
H36729-0025	25cm	0.8ml	3.0ml	40.20
H36729-0030	30cm	2.0ml	8.0ml	42.60

Fluo-Kem® Teflon® FEP coated ends.

CATALOG NO.	LENGTH	SMALL SPOON	LARGE SPOON	PRICE/PACK
F36730-0015	15cm	0.5ml	2.0ml	\$48.70
F36730-0021	21cm	1.6ml	5.0ml	54.85
F36730-0030	30cm	2.0ml	8.0ml	62.50

For Sample Storage Bottles, See Page 7-9



Ellipso-Spoon® Samplers

Unique Shape Means No Missed Sample in Corners

This utensil's specially curved shape permits contact with all corners of bottles and containers; therefore, no loss of material occurs. Crystals, powders, granulated and solid substances, and even thick liquids are easily retrieved and never fall back into the container. The Ellipso-Spoon® Sampler has a spoon end and a spatula end and is made of polished stainless steel. **12 of one size per case.**



CATALOG NO.	DESCRIPTION	DIAMETER OF HANDLE	LENGTH	PRICE/EACH
H36806-0015	15 x 35mm	5mm	15cm	\$16.20
H36806-0018	15 x 35mm	5mm	18cm	17.25
H36807-0021	15 x 35mm	5mm	21cm	18.45
H36807-0025	15 x 35mm	5mm	25cm	21.60
H36808-0030	28 x 65mm	8mm	30cm	32.25
H36809-0040	28 x 65mm	8mm	40cm	39.60
H36809-0050	28 x 65mm	8mm	50cm	47.25



Lab Spoon

Combination Spoon and Spatula with Fluo-Kem® Teflon®

Stainless steel spoon and spatula ends are coated with brown, Fluo-Kem® Teflon® FEP. This coating complies with FDA regulations governing direct food contact and provides an anti-stick, corrosion resistant surface. 4mm (¼") diameter rod is flattened at one end to give an 8 x 51mm (⅝ x 2") spatula. The opposite end is a 13 x 25mm (½ x 1") spoon. Overall length is 229mm (9"). **Per each, 24 per case.**



CATALOG NO.	PRICE/EACH
H36711-0000	\$19.20