

NEW




Scienceware® Mini Tong

These small lightweight tongs are ideal for one time use when sampling small quantities of material. The one piece tongs are molded of polystyrene and have ridges inside the jaws to securely hold wet and dry samples. The overall length of each is 11cm (4 1/4"). **50 per box.**



CATALOG NO.	PRICE/BOX
F37944-0025	\$41.25



Teflon® Coated, Steel Forceps

Meets FDA Regulations Involving Food Contact

Fluo-Kem® Teflon® FEP coating on the length of the tips protects against corrosion at low or high temperatures and provides an anti-stick surface to prevent material build up. Complying with FDA regulations covering direct contact with food, the brown coating

is soft and will not scratch or mar delicate surfaces. All forceps are 98mm (3 7/8") long with medium points. Select the coated length ideal for your use. **2 per pack, 24 packs per case.**



CATALOG NO.	COATED LENGTH	PRICE/PACK
A. H37936-0000	19mm (3/4" tip)	\$25.20
B. H37935-0000	38mm (1 1/2" tip)	25.20
C. H37937-0000	full length, 98mm (3 7/8")	26.49

Hot Hand® Protector for Safely Handling Hot & Cold Objects, See Page 122

Cover Glass Forceps

Wide Tip for Lifting Cover Glasses

Full length coating with brown Teflon® containing FEP which complies with FDA regulations concerning direct food contact. Available in two styles. **2 per bag, 24 bags per case.**



CATALOG NO.	DESCRIPTION	PRICE/BAG
H37942-0000	Bent 102mm (4") long	\$21.69
H37943-0000	Straight 114mm (4 1/2") long	24.45

Stainless Steel Tong

Teflon® Coated Tips for Quick Release

Laboratory type tongs are coated with brown Teflon® containing FEP which complies with FDA regulations governing direct food contact applications. 76mm (3") coated tips. Overall length 23cm (9"). **Per each, 24 per case.**



CATALOG NO.	PRICE/EACH
H37945-0000	\$21.69



A: Policeman

Made of a Teflon® TFE 7.9mm (5/16") diameter rod with flattened ends. 254mm (10") long. Steam autoclavable at 121°C (250°F). **Per each, 12 per case.**



CATALOG NO.	PRICE/EACH
H37738-0000	\$12.00

B: Policeman

Round molded polypropylene rod with one end flattened and cut diagonally. 4.7mm (3/16") diameter, 203mm (8") long. Steam autoclavable at 121°C (250°F). **6 per bag, 12 bags per case.**



CATALOG NO.	PRICE/BAG
F37739-0001	\$9.30

C: Large Policeman *Steel Reinforced Polyethylene*

The handle of this policeman is a rigid 6mm (1/4") steel rod enclosed in polyethylene. The tapered flattened end is 25mm (1 x 1"). Overall length 203mm (8"). Blister packed. **6 per bag, 12 packages per case.**



CATALOG NO.	PRICE/BAG
F37739-0002	\$47.25

E: Stirring Rod *Autoclavable to Avoid Contamination*

Solid, molded polypropylene stirring rod with rounded end, 4.7mm (3/16") diameter x 203mm (8") length is steam autoclavable at 121°C (250°F). **6 per bag, 12 bags per case.**



CATALOG NO.	PRICE/BAG
F37766-0008	\$7.75

Ladles

Safely Sample and Dispense Corrosive Liquids

High-density polyethylene ladles feature a rigid bowl with a handle securely welded at an ergonomically efficient angle. **Per each, 12 per case.**



CATALOG NO.	CAPACITY	HANDLE	PRICE/EACH
F36774-0000	50ml (1 1/2oz)	23cm (9")	\$15.75
F36776-0000	250ml (8 1/2oz)	23cm (9")	19.85