



All **Sterileware®** Sampling Products are shipped with an individual certificate of Sterilization and production number for accurate lot tracing.

NEW

Style 2 & 4oz Scoops!

← Longer Handle and Flat Bottom



Collect, Store and Transport Samples with Any One of These Fine Tools

NEW Improved Design!

Sterileware® Scoop Sampling System

This system eliminates the cost of separate sterile bags or containers, saving time and money in your sample collection process. Both sampling scoops

have been improved with a longer handle that increases the distance between the hand and sample to reduce contamination potential. These scoops have also been redesigned to remain upright when filled or empty allowing them to be used as weighing boats or to be placed on a flat surface. The polystyrene scoop, with a

specially designed detachable lid, is radiation sterilized in its own polyethylene bag. To use, simply remove the scoop from the bag and detach the lid using the convenient tab. After taking the sample, replace the snap-on lid and use the scoop as a container for transport. Sampling data can be marked on the lid with wax pencil, pre-marked label or labeling tape to facilitate tracking and to reduce risk of mix-ups. A lock-tab on the lid can be bent under for a tighter seal. Packed ready for one time use and disposal. **100 individually sterile wrapped per box.**



CATALOG NO.	CAPACITY	BOWL	TOTAL LENGTH	PRICE/BOX
H36922-0000	60ml (2oz)	25 x 76 x 38mm (1 x 3 x 1½")	203mm (8")	\$127.50
H36924-0000	125ml (4oz)	35 x 92 x 48mm (1½ x 3¾ x 1¾")	221mm (8⅞")	156.00

Sterileware® Scoops

3 Convenient Sizes for All Applications

Sampling scoops are the perfect shape and size for taking samples, especially from wide-mouth containers. Designed to hold various capacities, the one-piece scoop is molded of white high impact polystyrene, and each is individually sterilized and sealed in a polyethylene bag. The smaller 2 oz and 4 oz sampling scoops have been improved with a longer handle that increases the distance between the hand and sample. These scoops have also been redesigned to remain upright when filled or empty allowing them to be used as weighing boats or to be placed on a flat surface. Packed ready for one time use and disposal. **100 individually sterile wrapped per box.**

NEW Improved Design! 2 & 4oz



CATALOG NO.	CAPACITY	BOWL	TOTAL LENGTH	PRICE/BOX
H36902-0000	60ml (2oz)	25 x 76 x 38mm (1 x 3 x 1½")	203cm (8")	\$72.99
H36904-0000	125ml (4oz)	35 x 92 x 48mm (1½ x 3¾ x 1¾")	221cm (8⅞")	93.99
H36906-0000	250ml (8oz)	51 x 127 x 51mm	216cm (8½")	