



Komet™ Spinbar® Magnetic Stirring Bar

Improved Efficiency

The Teflon® Komet™ Spinbar® Magnetic Stirring Bar is made from a high-energetic magnetic material. That, combined with an octagonal shape, produces exceptional power, although its length is only 50mm (1.969") with a 21mm (.827") diameter. The Samarium-Cobalt magnet is encapsulated in Teflon® PTFE. Its eight edges

generate strong turbulence so liquids are thoroughly mixed and its stability is not impaired in vessels with curved bottoms. In addition, the risk of demagnetization caused by external magnetic fields is completely eliminated. Komet™ Spinbar® Magnetic Stirring Bar transmits torque loads 2 to 3 times larger than those of conventional stirring bars of this length, significantly improving efficiency. **Per each, 10 per case.**

CATALOG NO.	SIZE	PRICE/EACH
F37174-0000	50mm x 21mm (2 x 13/16")	\$62.40



Egg-Shaped Spinbar® Magnetic Stirring Bars

Teflon® PTFE

This egg-shaped Spinbar® Magnetic Stirring Bar provides its own pivot point without a separate ring or abrupt change in contour. Designed especially for use in round-bottom vessels, it spins unrestricted by the vessel contour. Most sizes have egg-shaped magnets. The three smallest sizes have cylindrical magnets with egg-shaped Teflon® PTFE coating. **Per each, 50 per case.**

Per each, 50 per case.

CATALOG NO.	FITS ROUND BOTTOM FLASKS	SIZE	PRICE/EACH
F37130-0038	-----	9.5 x 4.7mm (3/8 x 3/16")	\$6.25
F37130-0058	-----	15.9 x 6.35mm (5/8 x 1/4")	7.00
F37130-0034	10,25,50,100ml	19.1 x 9.5mm (3/4 x 3/8")	7.20
F37130-0001	100 and 200ml	25.4 x 12.7mm (1 x 1/2")	7.15
F37130-0114	300 and 500ml	31.8 x 15.9mm (1 1/4 x 5/8")	9.00
F37130-0112	300 and 500ml	38.1 x 15.9mm (1 1/2 x 5/8")	9.80
F37130-0158	500 and 1000ml	41.3 x 19mm (1 1/4 x 3/4")	12.25
F37130-0002	2000 and 3000ml	50.8 x 19.1mm (2 x 3/8")	14.65
F37130-0212	3000 and 5000ml	63.5 x 19mm (2 1/2 x 3/8")	18.55
F37130-0003	12000 and 22000ml	76 x 19.1mm (3 x 3/8")	21.35



Flea Micro Spinbar® Magnetic Stirring Bars

Teflon® PTFE, White

Per each, 50 per case.

CATALOG NO.	SIZE	PRICE/EACH
F37119-0002	2 x 2mm (.079" x .079")	\$3.90
F37119-0003	3 x 3mm (.118" x .118")	3.85
F37119-0005	5 x 2mm (.197" x .079")	3.50
F37119-0635	6.35 x 3mm (.250" x .118")	3.85
F37119-0007	7 x 2mm (.276" x .079")	3.50
F37119-0008	8 x 1.5mm (.315" x .059")	3.50
F37119-1083	8 x 3mm (.315" x .118")*	4.35
F37119-0083	8 x 3mm (.315" x .118")	3.85
F37119-0010	10 x 3mm (.394" x .118")	3.85
F37119-0127	12.7 x 3mm (.500" x .118")	3.85
F37119-0015	15 x 1.5mm (.591" x .059")	3.45



Flea Micro Spinbar® Magnetic Stirring Bars

Teflon® PTFE, Color-Coded

Per each, 50 per case.

CATALOG NO.	COLOR	SIZE	PRICE/EACH
F37121-0007	red	2 x 2mm (.079" x .079")	\$4.30
F37121-0008	yellow	2 x 2mm (.079" x .079")	4.30
F37121-0009	blue	2 x 2mm (.079" x .079")	4.30
F37121-0016	red	3 x 3mm (.118" x .118")	4.25
F37121-0017	yellow	3 x 3mm (.118" x .118")	4.25
F37121-0018	blue	3 x 3mm (.118" x .118")	4.25
F37121-0010	red	5 x 2mm (.197" x .079")	3.85
F37121-0011	yellow	5 x 2mm (.197" x .079")	3.85
F37121-0012	blue	5 x 2mm (.197" x .079")	3.85
F37121-0019	red	6.35 x 3mm (.250" x .118")	4.25
F37121-0020	yellow	6.35 x 3mm (.250" x .118")	4.25
F37121-0021	blue	6.35 x 3mm (.250" x .118")	4.25
F37121-0013	red	7 x 2mm (.276" x .079")	3.85
F37121-0014	yellow	7 x 2mm (.276" x .079")	3.85
F37121-0015	blue	7 x 2mm (.276" x .079")	3.85
F37121-0001	red	8 x 1.5mm (.315" x .059")	3.85
F37121-0002	yellow	8 x 1.5mm (.315" x .059")	3.85
F37121-0003	blue	8 x 1.5mm (.315" x .059")	3.85
F37121-0022	red	10 x 3mm (.394" x .118")	4.25
F37121-0023	yellow	10 x 3mm (.394" x .118")	4.25
F37121-0024	blue	10 x 3mm (.394" x .118")	4.25
F37121-0025	red	12.7 x 3mm (.500" x .118")	4.25
F37121-0026	yellow	12.7 x 3mm (.500" x .118")	4.25
F37121-0027	blue	12.7 x 3mm (.500" x .118")	4.25
F37121-0004	red	15 x 1.5mm (.591" x .059")	3.85
F37121-0005	yellow	15 x 1.5mm (.591" x .059")	3.85
F37121-0006	blue	15 x 1.5mm (.591" x .059")	3.85



Round Spinbar® Magnetic Stirring Bar with Tapered Ends

Teflon® PTFE

Tapered design acts as a pivot point, eliminating the need for a separate pivot ring. End diameters are 1mm less than center diameter. **Per each, 50 per case.**

CATALOG NO.	SIZE	PRICE/EACH
F37113-0208	20 x 8mm (.787" x .315")	\$2.85
F37113-0309	30 x 9mm (1.181" x .354")	3.40
F37113-0409	40 x 9mm (1.575" x .354")	3.70
F37113-0509	50 x 9mm (1.969" x .354")	4.30