



**Octagon Spinbar®  
Magnetic Stirring  
Bars**

Teflon® PTFE

With Molded-On Pivot Ring. Per each, 50 per case.

CATALOG NO.	SIZE	PRICE/EACH
F37110-1218*	12.7 x 3.2mm (½ x ¼")	\$2.85
F37110-1256	12.7 x 8mm (½ x ⅝")	3.20
F37110-1238	12.7 x 9.5mm (½ x ⅜")	3.70
F37110-0058	15.9 x 8mm (⅝ x ⅝")	3.20
F37110-5838	15.9 x 9.5mm (⅝ x ⅜")	4.15
F37110-0078	22.2 x 8mm (⅞ x ⅝")	3.45
F37110-7838	22.2 x 9.5mm (⅞ x ⅜")	4.30
F37110-0001	25.4 x 8mm (1 x ⅝")	3.50
F37110-0138	25.4 x 9.5mm (1 x ⅜")	4.30
F37110-0118	28.6 x 8mm (1 ⅛ x ⅝")	3.70
F37110-1388	34.9 x 9.5mm (1 ⅜ x ⅜")	5.10
F37110-0112	38.1 x 8mm (1 ½ x ⅝")	3.95
F37110-1128	38.1 x 9.5mm (1 ½ x ⅜")	5.10
F37110-1122	38.1 x 12.7mm (1 ½ x ½")	12.95
F37110-0158	41.3 x 8mm (1 ⅝ x ⅝")	4.55
F37110-0178	47.6 x 9.5mm (1 ⅞ x ⅜")	6.45
F37110-0002	50.8 x 8mm (2 x ⅝")	4.85
F37110-0238	50.8 x 9.5mm (2 x ⅜")	6.85
F37110-0212	63.5 x 8mm (2 ½ x ⅝")	7.75
F37110-2128	63.5 x 9.5mm (2 ½ x ⅜")	8.00
F37110-0003	76.2 x 12.7mm (3 x ½")	14.20

\*No pivot ring on this size.



**Round Spinbar®  
Magnetic Stirring  
Bars**

Teflon® PTFE

With Removable Pivot Ring. Per each, 50 per case.

CATALOG NO.	SIZE	PRICE/EACH
F37111-0012*	12.7 x 8mm (½ x ⅝")	\$4.10
F37111-1238	12.7 x 9.5mm (½ x ⅜")	5.25
F37111-0034*	19.1 x 8mm (¾ x ⅝")	4.15
F37111-0001	25.4 x 8mm (1 x ⅝")	4.30
F37111-0138	25.4 x 9.5mm (1 x ⅜")	5.60
F37111-0114	31.8 x 8mm (1 ¼ x ⅝")	5.45
F37111-1148	31.8 x 9.5mm (1 ¼ x ⅜")	6.00
F37111-0112	38.1 x 8mm (1 ½ x ⅝")	5.65
F37111-1128	38.1 x 9.5mm (1 ½ x ⅜")	7.15
F37111-0158	41.3 x 9.5mm (1 ⅝ x ⅜")	8.00
F37111-0134	44.5 x 8mm (1 ⅞ x ⅝")	7.00
F37111-1348	44.5 x 9.5mm (1 ⅞ x ⅜")	7.65
F37111-0002	50.8 x 8mm (2 x ⅝")	7.35
F37111-0238	50.8 x 9.5mm (2 x ⅜")	8.40
F37111-0218	54.0 x 9.5mm (2 ⅛ x ⅜")	9.25
F37111-2125	63.5 x 16mm (2 ½ x ⅝")	32.75
F37111-0003	76.2 x 12.7mm (3 x ½")	16.50
F37111-0004**	101.6 x 16mm (4 x ⅝")	41.95
F37111-0005**	127.0 x 16mm (5 x ⅝")	48.00
F37111-0006**	152.4 x 19mm (6 x ⅜")	118.80

\*No pivot ring on this size

\*\*10 bars per case

**Metric Polygon Spinbar® Magnetic Stirring  
Bars**

With Center Molded Pivot Ring. Teflon® PTFE

The molded pivot ring minimizes the contact area between the bar and the vessel, thus reducing friction and chattering. Marring of plastic containers is also diminished. Multifaceted surface adds turbulence relative to similar size smooth cylindrical stirring bars. **Per each, 50 per case.**

CATALOG NO.	SIZE	PRICE/EACH
F37122-0012	12 x 4.5 mm (½ x ⅝")	\$3.40
F37122-0020	20 x 6 mm (¾ x ¼")	4.10
F37122-0025	25 x 6 mm (1 x ¼")	4.20
F37122-0035	35 x 6 mm (1 ⅜ x ¼")	4.35
F37122-0040	40 x 8 mm (1 ⅞ x ⅝")	4.95
F37122-0045	45 x 8 mm (1 ¾ x ⅝")	5.40
F37122-0050	50 x 8 mm (2 x ⅝")	5.40
F37122-0060	60 x 10 mm (2 ⅜ x ⅝")	5.50
F37122-0070	70 x 10 mm (2 ¾ x ⅝")	9.25

**Giant Polygon Spinbar®  
Magnetic Stirring Bars**

With Center Molded Pivot Ring Teflon® PTFE

Designed for stirring large volumes in large vessels such as drums and tanks. The molded pivot ring minimizes the contact area between the bar and the vessel, thus reducing friction and chattering. Marring of plastic containers is also diminished. Multifaceted surface adds turbulence relative to similar size smooth cylindrical stirring bars. Ideal for use with Scienceware® Cool Stir® Large Volume Stirrer, H37028-0000, which is capable of stirring volumes up to 190 L (50 gallons). **Per each, 10 per case.**

CATALOG NO.	SIZE	PRICE/EACH
F37122-0005	127 x 16 mm (5 x ⅝")	\$92.50
F37122-0006	152 x 19 mm (6 x ⅜")	102.80



**Octagon Shaped  
Spinbar® Magnetic  
Stirring Bars**

Color-Coded Teflon® PTFE

All PTFE coded Spinbar® Magnetic Stirring Bars are available in red, yellow or blue. **Per each, 50 per case.**

CATALOG NO.	COLOR	SIZE	PRICE/EACH
F37109-0034*	red	12.7 x 3.2mm (½ x ⅝")	\$3.90
F37109-0035*	yellow	12.7 x 3.2mm (½ x ⅝")	3.90
F37109-0036*	blue	12.7 x 3.2mm (½ x ⅝")	3.90
F37109-0001	red	12.7 x 8mm (½ x ⅝")	3.90
F37109-0002	yellow	12.7 x 8mm (½ x ⅝")	3.90
F37109-0003	blue	12.7 x 8mm (½ x ⅝")	3.90
F37109-0004	red	15.9 x 8mm (⅝ x ⅝")	4.15
F37109-0005	yellow	15.9 x 8mm (⅝ x ⅝")	4.15
F37109-0006	blue	15.9 x 8mm (⅝ x ⅝")	4.15

CATALOG NO.	COLOR	SIZE	PRICE/EACH
F37109-0007	red	22.2 x 8mm (⅞ x ⅝")	\$4.40
F37109-0008	yellow	22.2 x 8mm (⅞ x ⅝")	4.40
F37109-0009	blue	22.2 x 8mm (⅞ x ⅝")	4.40
F37109-0010	red	25.4 x 8mm (1 x ⅝")	4.65
F37109-0011	yellow	25.4 x 8mm (1 x ⅝")	4.65
F37109-0012	blue	25.4 x 8mm (1 x ⅝")	4.65
F37109-0019	red	38.1 x 8mm (1 ½ x ⅝")	5.70
F37109-0020	yellow	38.1 x 8mm (1 ½ x ⅝")	5.70
F37109-0021	blue	38.1 x 8mm (1 ½ x ⅝")	5.70
F37109-0028	red	50.8 x 8mm (2 x ⅝")	6.10
F37109-0029	yellow	50.8 x 8mm (2 x ⅝")	6.10
F37109-0030	blue	50.8 x 8mm (2 x ⅝")	6.10
F37109-0031	red	76.2 x 12.7mm (3 x ½")	19.00
F37109-0032	yellow	76.2 x 12.7mm (3 x ½")	19.00
F37109-0033	blue	76.2 x 12.7mm (3 x ½")	19.00

\*No pivot ring on this size.