

# Safe-Lab® Stoppers


Patent numbers D280,080 and D289,378

## How Safe-Lab® Stoppers Work

Safe-Lab® Stoppers feature an exclusive, patented safety device for easy and simple separation of fused glassware. This innovative extracting mechanism minimizes expensive glassware breakage and potential bodily injuries when separating glassware with frozen ground joints.


No More...

- tapping glassware against a bench top
- running glassware under warm water to loosen
- risking injury with dangerous removal techniques



Safe-Lab® stoppers are easily removed from standard tapered glassware by turning the red extraction nut that applies a safe, gentle and even force against the mouth of the vessel, helping to eliminate accidents and breakage.

“Easy on, easy off” Safe-Lab® stoppers save time and prevent accidents with standard tapered glassware.



### Safe-Lab® Teflon® PTFE Stoppers

Use with Separatory Funnels And Volumetric Flasks

Solid Teflon® PTFE stopper with threaded polypropylene extracting nut for easy removal of stopper from flask. Interchangeable with penny head stoppers. *Per pack listed below.*



CATALOG NO.	STOPPER SIZE	PACK	PRICE/PACK
F20191-1465	8	3	\$34.00
F20191-1470	13	3	36.35
F20191-1480	16	3	41.20
F20191-1490	19	3	60.50
F20191-1495	22	3	63.50
F20191-1505	27	3	64.50
F20191-1520	32	2	65.00
F20191-1545	38	Each	53.00



### Safe-Lab® Solid Teflon® PTFE Stoppers

Use with Boiling and Distilling Flasks

Solid Teflon® PTFE stopper with extracting nut for flasks with ground joints. *Per pack listed below.*



CATALOG NO.	JOINT SIZE	PACK	PRICE/PACK
F20196-1485	14/20	3	\$60.00
F20196-1500	19/22	3	63.50
F20196-1502	24/40	2	62.85

### Safe-Lab® Hollow Teflon® PTFE Stoppers

Use these hollow Teflon® PTFE stoppers for boiling and distilling flasks with ground joints. Stopper comes with polypropylene extracting nut for easy removal from flask. *Per pack listed below.*



CATALOG NO.	STOPPER SIZE	PACK	PRICE/PACK
F20197-1515	24/40	2	\$65.95
F20197-1530	29/42	2	75.00
F20197-1540	34/45	Each	51.00
F20197-1550	45/50	Each	99.95
F20197-1555	55/50	Each	118.50

### Safe-Lab® Glass Stoppers

Interchangeable with Penny Head Stoppers

Lightweight hollow borosilicate glass stopper. Flat hex head for easy gripping. Threaded polymethylpentene or polypropylene extracting nut above ground surface allows easy and safe removal of stopper from flask.



CATALOG NO.	JOINT SIZE	THREAD SIZE	NUT	PRICE/EACH
H20199-1300	14/20	14	PMP	\$10.50
H20199-1320	24/40	24	PMP	15.90
H20199-1330	29/42	29	PP	18.00