

**NEW**

### Stackable Petri Dish Incubation Tray

Increase overall capacity of your incubator with Stackable Petri Dish Incubation Trays. Designed to hold five 100mm diameter or fifteen 60mm diameter Petri dishes, each polypropylene tray is provided with two large label areas for critical sample identification by marker or label. **3 per bag, 4 bags per case.**

PP

CATALOG NO.	DIMENSIONS	PRICE/BAG
F18983-1000	250.8 x 236.5 x 35 cm H (97/8 x 95/16 x 13/8" H)	\$33.50



### Petri Dish Incubation Tray

Organize Your Incubator

The Petri dish incubator tray holds 100mm (4") Petri dishes securely on an incubator rack. It is an excellent aid for keeping cell cultures organized and has 28 holes 19mm (3/4") in diameter in the bottom to allow air passage. Made from 1.6mm (1/16") thick, white epoxy-coated steel. 349 x 137 x 22mm H (13 3/4 x 5 1/2 x 7/8").

STEEL A

CATALOG NO.	PRICE/EACH
F18983-0000	\$39.75

### Compartment Trays

A practical and functional variety of vacuum-formed, high impact, heavy walled, white polystyrene trays that will eliminate lost and broken lab equipment. All trays measure 36 x 43cm (14 x 17") and fit into standard lab bench drawers. The trays are great for organizing:

- Long Items
- Thermometers
- Vials
- Stirring Bars and Retrievers
- Gadgets
- Stoppers
- Round Containers
- Medication Cups
- and more

PS

#### COMPARTMENT TRAY

Superb for holding pipettes and other slim, fragile items up to 406mm (16") long. Two compartments are 102mm (4") wide, with a third compartment 140mm (5 1/2") wide all are 64mm (2 1/2") deep. **Per each, 6 per case.**

CATALOG NO.	PRICE/EACH
H18648-0000	\$33.99

#### BEAKER, FLASK OR JAR TRAY

Round containers up to 98mm (3 7/8") in diameter nest into molded wells making this tray perfect for holding flasks and beakers up to 600ml. Features 20 wells of various diameters. **Per each, 6 per case.**

CATALOG NO.	PRICE/EACH
H18650-0000	\$33.99

#### GADGET ORGANIZER

A 12 compartment tray that is invaluable for storing small items. Each compartment is 98 x 98 x 73mm deep (3 7/8 x 3 7/8 x 2 7/8"). **Per each, 6 per case.**

CATALOG NO.	PRICE/EACH
H18652-0000	\$36.39

### Stak-A-Tray™ System

Design Your Own Stackable Storage

Modular rack and tray system maximizes storage for Petri dishes and cell culture flasks in incubators, cold rooms or on the bench. Stack five racks in one square foot for up to 20 levels of storage. Perforated trays allow air flow and slide easily in the rack. The system accommodates either one large or two small trays on each of four tiers. All parts are made of stainless steel and are autoclavable at 121°C (250°F). Rack overall dimension: 37.5 x 10.8 x 35.8cm H (14 1/4 x 4 1/4 x 14 1/16"H). Large tray: 35.6 x 35.6 cm (14 x 14"), small tray: 17.8 x 35.6 cm (7 x 14"). Racks and trays sold separately.

STAINLESS STEEL A

CATALOG NO.	DESCRIPTION	PRICE/EACH
H18610-0420	Rack with two center supports (for larger culture flasks) 5cm (1.96") clearance between trays. Rack and 2 supports only.	\$231.75
H18610-0440	Rack with four center supports (for Petri dishes and small culture flasks) 2.4cm (.96") clearance between trays. Rack and four supports only.	231.75
H18610-0470	Small tray, 17.8 x 35.6cm (7 x 14")	43.29
H18610-1440	Large tray, 35.6 x 35.6cm (14 x 14")	60.30

### Stack-a-Tray™ System, Storage Capacity

ITEM SIZE	LARGE TRAY HOLDS	SMALL TRAY HOLDS
35 x 100mm dish	64	32
5 x 12mm dish	36	18
60 x 12mm plastic dish	36	17
60 x 15mm glass dish	25	12
100 x 10mm dish	13	6
100 x 15mm dish	13	6
100 x 20mm dish	13	6
150 x 20mm dish	4	2
150 x 25mm dish	4	2
125 x 15mm square	4	2



#### FIVE COMPARTMENT TRAY

Store stirring bars and retrievers in separate areas of the same tray. The center compartment is 102 x 406mm (4 x 16") for long items. Four side compartments are all 191mm (7 1/2") long and 64mm (2 1/2") deep. Two of the compartments are 102mm (4") wide and two are 140mm (5 1/2") wide. **Per each, 6 per case.**

CATALOG NO.	PRICE/EACH
H18654-0000	\$36.39



#### OPEN TRAY

This open, flat tray holds a variety of large objects or light bulk materials. 64mm (2 1/2") deep. **Per each, 6 per case.**

CATALOG NO.	PRICE/EACH
H18655-0000	\$36.39



#### FOUR COMPARTMENT TRAY

Four compartments, each 197 x 152 x 57mm (7 3/4 x 6 x 2 1/4") deep. **Per each, 6 per case.**

CATALOG NO.	PRICE/EACH
H18658-0000	\$33.99