



The Collector™ Blood Tray

Everything You Need to Carry from Lab to Patient

Designed for flexibility and versatility, The Collector™ Blood Tray can be rearranged to meet your specific needs. It includes two

40-place test tube supports (for tubes up to 16mm O.D.), one 18-place slide holder, a large requisition section, and adjustable partitions to hold needles, syringes, cotton and other materials. Sturdy polypropylene construction holds its shape even when fully loaded. Steam autoclavable at 121°C (250°F). Dimensions: 356 x 292 x 51mm H (14 x 11½ x 2"H). The handle is formed by an end mounted, detachable loop that is 178mm H (7"H). **Per each, 4 per case.**



CATALOG NO.	DESCRIPTION	PRICE/EACH
A. F18631-0715	Complete Blood Tray Set	\$99.00
B. F18631-0716	Replacement test tube support. 167 x 109 x 109H mm (6½ x 4¼ x 4¼") 40 places for tubes up to 16mm O.D.	17.00



Benchtop Working Tray†

The polypropylene Benchtop Working Tray is designed to hold the most popular sizes of safety mats to help contain spills. It can be used in combination with biohazard,

radioactive or general purpose mats. The lip of the tray assures containment of spilled liquids that are not absorbed by the mats. Inside dimensions of the tray are 46 x 56 x 1.3cm (18⅛ x 22⅛ x ½).



CATALOG NO.	PRICE/EACH
F24675-0060	\$62.50



Trays†

Offered in Three Popular Styles

Our lab trays are made from two materials; choose the one that suits your needs. Polypropylene is autoclavable and offers resistance to temperatures up to 135°C (274°F). Economical high-density polyethylene is available with a spigot.



CATALOG NO.	DESCRIPTION	LENGTH x WIDTH x HEIGHT	PRICE/EACH
F16300-1015	PP	254 x 381 x 76mm (10 x 15 x 3")	\$65.50
F16300-1216	PP	305 x 406 x 203mm (12 x 16 x 8")	95.75
F16300-1520	PP	381 x 508 x 76mm (15 x 20 x 3")	85.25
F16300-1818	PP	457 x 457 x 102mm (18 x 18 x 4")	97.75
F16300-2020	PP	508 x 508 x 152mm (20 x 20 x 6")	127.00
F16300-2226	PP	559 x 660 x 102mm (22 x 26 x 4")	130.00

CATALOG NO.	DESCRIPTION	LENGTH x WIDTH x HEIGHT	PRICE/EACH
F16301-1015	HDPE	254 x 381 x 76mm (10 x 15 x 3")	\$50.00

CATALOG NO.	DESCRIPTION	LENGTH x WIDTH x HEIGHT	PRICE/EACH
F16302-1818	HDPE w/spigot	457 x 457 x 102mm (18 x 18 x 4")	89.00



Heavy Duty Rectangular Tanks with Top Flanges

Several Styles to Suit a Variety of Applications

Seamless, rotational-molded polyethylene tanks are lightweight and easy to handle. Available with or without a faucet, these tanks are resistant to most acids and all alkalis. All sizes are made with wide top flanges, include covers, and can be supplied with spigots, fittings, dams, partitions, etc. Contact Bel-Art® for quotations on special sizes and fabrication in polypropylene and PVC.



Without Faucet

CATALOG NO.	CAPACITY	DIMENSIONS (L X W X H)	WALL	PRICE/EACH
H34080-0000	19 liters (5gal)	28 x 28 x 25cm (11 x 11 x 10")	6.4mm (¼")	\$63.00
H34082-0000	39 liters (10gal)	46 x 33 x 25cm (18 x 13 x 10")	6.4mm (¼")	84.00
H34085-0000	57 liters (15gal)	39 x 30 x 48cm (15.25 x 12 x 19")	9.5mm (¾")	104.49
H34090-0000	78 liters (20gal)	53 x 41 x 36cm (21 x 16 x 14")	9.5mm (¾")	114.90

With Faucet

CATALOG NO.	CAPACITY	DIMENSIONS (L X W X H)	WALL	PRICE/EACH
H34092-0000	19 liters (5gal)	28 x 28 x 25cm (11 x 11 x 10")	6.4mm (¼")	\$87.90
H34093-0000	28 liters (7.5gal)	42 x 28 x 25cm (16.5 x 11 x 10")	6.4mm (¼")	99.00
H34094-0000	39 liters (10gal)	46 x 33 x 25cm (18 x 13 x 10")	6.4mm (¼")	114.90
H34095-0000	57 liters (15gal)	39 x 30 x 48cm (15.25 x 12 x 19")	9.5mm (¾")	129.00
H34096-0000	78 liters (20gal)	53 x 41 x 36cm (21 x 16 x 14")	9.5mm (¾")	160.50



Polly-Crocks with Lids

Several Styles to Suit a Variety of Applications

Heavy-duty, rotational molded polyethylene tanks will safely hold chemicals and hot liquids up to 52°C (125°F) and will meet FDA requirements for food handling. The open-top tanks are stress-free and self-supporting with a seamless wall about 4.8 mm (3/16") thick. Resistant to rust, corrosion and dents, each container's capacity is marked on both the tank and the lid. Lids have a central grip knob for easy removal and replacement. Also available with the Scienceware® F30852 faucet, (for replacement faucet see page xx). Polly-Crocks are also offered with a sight glass as a special order.



CATALOG NO.	CAPACITY	DIAM. TOP X BOTTOM	HEIGHT	WEIGHT	PRICE
LITERS GALLONS				KG LB	EACH
H35103-0000	11.4 3	28 x 24cm (11 x 9½")	28cm (11")	2.0 4.5	\$58.50
H35106-0000	22.7 6	31 x 28cm (12¼ x 11")	41cm (16")	3.0 6.6	72.00
H35110-0000	37.9 10	38 x 34cm (15 x 13½")	46cm (18")	4.0 8.8	90.00
H35115-0000	56.3 15	43 x 39cm (16¾ x 15¼")	52cm (20½")	5.1 11.3	108.00

With Faucet

CATALOG NO.	CAPACITY	DIAM. TOP X BOTTOM	HEIGHT	WEIGHT	PRICE
LITERS GALLONS				KG LB	EACH
H35223-0000	11.4 3	28 x 24cm (11 x 9½")	28cm (11")	2.0 4.5	\$87.00
H35226-0000	22.7 6	31 x 28cm (12¼ x 11")	41cm (16")	3.0 6.6	105.00
H35230-0000	37.9 10	38 x 34cm (15 x 13½")	46cm (18")	4.0 8.8	118.50
H35235-0000	56.3 15	43 x 39cm (16¾ x 15¼")	52cm (20½")	5.1 11.3	135.00