



## Complete Fitting Assortment

74 Pieces to Make the Right Connection

Stop wasting time hunting for the right tubing connection. This generous assortment of 74 fittings has everything you need for routine laboratory set-ups. A clear, hinged-top, polystyrene storage box makes it easy to find what you need, when you need it.

Kit includes two of each size of the following:

### Polypropylene Connectors:

- 10 Straight Tubing connectors
- 10 "T" Connectors
- 12 Straight-Stepped Tubing Connectors
- 12 "Y" Connectors



### Polyethylene Connectors:

- 6 Quick Disconnects
- 24 Cup-Type Stoppers



CATALOG NO.	DESCRIPTION	PRICE/BOX
F19625-0000	Fitting Assortment	\$85.00



## Stepped Tubing Connectors

Polypropylene. Autoclavable.

12 per bag



CATALOG NO.	TUBING I.D.	PRICE/BAG
H19558-0000	4.7 x 6mm (5/16 to 1/4")	\$16.20
H19559-0000	4.7 x 9.5mm (5/16 to 3/8")	16.20
H19560-0000	4.7 x 12.7mm (5/16 to 1/2")	16.20
H19561-0000	6 x 9.5mm (1/4 to 3/8")	16.20
H19562-0000	6 x 12.7mm (1/4 to 1/2")	16.20
H19563-0000	9.5 x 12.7mm (3/8 to 1/2")	16.20



## Straight Connectors

Polypropylene. Autoclavable.

12 per bag.



CATALOG NO.	TUBING I.D.	PRICE/BAG
H19503-0000	3.2mm (1/8")	\$7.80
H19507-0000	4.7mm (3/16")	16.20
H19508-0000	6mm (1/4")	16.20
H19509-0000	9.5mm (3/8")	16.20
H19510-0000	12.7mm (1/2")	16.20

## Connector Assortment

Always Have the Right Size Available...  
Keep an Assortment on Hand!

An assortment of polypropylene straight and straight-stepped tubing connectors. 22 per bag (2 of each size, 19558-19563 and 19503-19510, as above).



CATALOG NO.	PRICE/BAG
H19570-0000	\$26.70

## Quick Disconnects

Interchangeable Parts to Make Reducing Combinations



Full range of sizes to fit tubing from 3 to 12mm (1/8 to 1/2") I.D. Consists of two serrated polyethylene connectors which assemble tightly together with a male-female center taper. Very useful in connecting and disconnecting vacuum lines and other tubing assemblies. Use on filtering flasks to avoid accidents. All sizes are interchangeable to make reducing combinations. 12 per bag, 6 bags per case.



CATALOG NO.	FITS TUBING I.D.	PRICE/BAG
H19727-0000	3 to 6mm (1/8 to 1/4")	\$19.80
H19728-0000	6 to 9mm (1/4 to 3/8")	19.80
H19729-0000	9 to 12mm (3/8 to 1/2")	19.80

## Quick Disconnects

For Same or Different Sized Tubing



Polyethylene connectors are ideal for vacuum applications, water lines, manometers, etc. The F19731-0552 connects different size tubing. 12 per bag, 6 bags per case.



CATALOG NO.	MALE FITTING ID	FEMALE FITTING ID	OVERALL LENGTH	PRICE/BAG
F19730-0001	4.7 to 6.4mm (5/16 to 1/4")	4.7 to 6.4mm (5/16 to 1/4")	60mm (2 3/8")	\$30.65
F19730-0002	6.4 to 9.5mm (1/4 to 3/8")	6.4 to 9.5mm (1/4 to 3/8")	60mm (2 3/8")	17.25
F19731-0552	4.7 to 6.4mm (5/16 to 1/4")	6.4 to 7.9mm (1/4 to 5/16")	60mm (2 3/8")	29.65

## Check Valves

Use in Vacuum Lines to Filter Flasks



These polyethylene check valves operate positively and instantly in any position. Arrows indicate direction of free flow. Neoprene disk, no metal. 12 per bag, 6 bags per case.



CATALOG NO.	BORE I.D.	OVERALL LENGTH	TAPER	PRICE/BAG
F19713-0001	2.4mm (3/32")	70mm (2 3/4")	7.9 to 6.4mm (5/16 to 1/4")	\$37.00
F19713-0002	6.4mm (1/4")	81mm (3 1/8")	12.7 to 9.5mm (1/2 to 3/8")	51.00

## "5 in 1" Connectors

Adapts to 15 Combinations as Connector or Reducer



The "five-in-one" fits 5 sizes of tubing and will connect the same size tubing on both ends or will act as a reducer for any combination within the range. The inside diameter increases with notch size when unneeded length is cut off for unrestricted flow. The molded polypropylene connectors are steam autoclavable at 121°C (250°F). Ring diameters are 14, 12.5, 11, 9.5 and 7.9mm (1/2, 1/2, 1/2, 3/8 and 5/16"). 12 per bag, 6 bags per case.



CATALOG NO.	DESCRIPTION	DIMENSIONS	PRICE/BAG
F19655-0000	Straight connector 5-in-1	90mm (3 3/8") long	\$23.25
F19657-0000	"T" connector 5-in-1	118mm (4 5/8") long x 67mm (2 5/8") branch	46.25