




Insert gloved hand into bag.



Pick up waste materials with hand inside bag.



Peel bag over the hand, turning it inside out while maintaining grip on waste material.



Waste is now inside bag and ready for autoclaving.



Bench-Top Biohazard Bag

Easy Open, Convenient Size

Autoclavable bags feature one long flap for ease in opening. When used with the

Poxygrid® H13193-1000 desk-size bag holder, they provide a convenient, space saving disposal unit for use in fume hoods, or on benchtops and desks. The bright orange-red bags are made of High Molecular Weight, High-Density (HMHD) polyethylene and are printed with the biohazard symbol and a sterilization indicator patch. When picking up individual pieces of waste, the bag is turned inside out, over the hand, then turned right side out with waste safely inside the bag. Measures 22 x 28cm (8½ x 11") and is .038mm (.0015") thick. Steam autoclavable at 121°C (250°F). **Per box, 100 or 1000.**



CATALOG NO.	DESCRIPTION	PRICE/BOX
F13166-0000	Benchtop Biohazard bag	\$11.85 /100
F13166-0001	Benchtop Biohazard bag	108.00 /1000

Poxygrid® Desk Model Bag Holder

Handy Disposal Holder for Used Tips and Tubes

Economical and convenient, this rectangular frame holds a polyethylene bag 22 x 28cm (8½ x 11") which is just the right size for collecting pipette tips, microcentrifuge tubes, etc. during daily use. The steel wire holder is coated in durable biohazard orange epoxy and measures 8.5 x 12.7 x 21.6cm H (3¼ x 5 x 8½"). Use with Scienceware® Bench Top Biohazard Bags H13182-0812 and F13166-0000.



CATALOG NO.	PRICE/EACH
H13193-1000	\$15.99

Easy to use clips on page 4



A. Biohazard Sharp Object Safety Pouch

Hard-to-Pierce Pouch Offers Greater Protection

Avoid the danger of being cut by contaminated sharp objects. Use the Biohazard Sharp Object Safety Pouch for bench-top collection and disposal. Made of 10mil thick solid paperboard which opens to an oval shape by folding bottom flaps. After filling, top flaps close and lock in place for safe disposal. Bold biohazard symbol for easy visibility. Measures 14 x 33cm (5½ x 13") flat. **200 per pack, 5 packages per case.**



CATALOG NO.	PRICE/PACK
H13234-0000	\$124.50



B. Poxgrid® Safety Pouch Stand

The bright biohazard orange epoxy coated steel wire stand provides a stable support to keep the Biohazard Sharp Objects Safety Pouch upright and open on the bench for easy access. 10 x 13 x 20cm H (4 x 5½ x 8").



CATALOG NO.	PRICE/EACH
F13236-0000	\$13.95

Clavies® Autoclavable Bags

Safe, Convenient Biohazard Disposal

For maximum protection and convenience when disposing of biohazardous wastes, Clavies® transparent, autoclave bags are made of .05mm (.002") thick polypropylene. Fill the bag with used Petri dishes, membrane filters, culture plates, pipettes, etc., then autoclave and discard. The three largest sizes hold 65, 200 and 240 standard Petri dishes respectively. Steam autoclavable at 121°C (250°F). **100 per box, 10 boxes per case.**



CATALOG NO.	SIZE	PRICE/BOX
H13182-0810	20 x 25cm (8 x 10")	\$26.70
H13182-0812	20 x 30cm (8 x 12")	31.89
H13182-0815	20 x 38cm (8 x 15")	39.00
H13182-1010	25 x 25cm (10 x 10")	28.80
H13182-1012	25 x 30cm (10 x 12")	34.50
H13182-1015	25 x 38cm (10 x 15")	39.99
H13185-1224	30 x 61cm (12 x 24")	44.25
H13185-2430	61 x 76cm (24 x 30")	105.00
H13185-2436	61 x 91cm (24 x 36")	123.60

Odo-Clave® Deodorant Pads Neutralize Organic Waste Odors in the Autoclave, See Page 4

Polyethylene Utility Bags

All Purpose Bags, Available in Two Weights

Reusable, transparent polyethylene bags fill a variety of needs. These are ideal for shipping, storage and are especially useful as sample bags. Not autoclavable. FDA approved material. **The H13176 lightweight bags are .038 mm (.0015") thick, packed 100 per box, 15 boxes per case. The H13178 heavyweight bags are .1016 mm (.004") thick, packed 100 per box, 10 boxes per case.**



Order by complete catalog number and size code.



CATALOG NO.		SIZE CODE	SIZE	PRICE/BOX	
.038MM	.1016MM			.038MM	.1016MM
—	H13178	-0046	102 x 152mm (4 x 6")	\$—	\$11.25
—	H13178	-0068	152 X 216mm (6 x 8.5")	—	14.10
—	H13178	-0610	152 x 254mm (6 x 10")	—	16.95
H13176	H13178	-0612	152 x 305mm (6 x 12")	10.50	18.99
—	H13178	-0810	203 x 267mm (8 x 10.5")	—	16.95
H13176	H13178	-0812	203 x 305mm (8 x 12")	11.70	18.00
H13176	—	-0814	203 x 356mm (8 x 14")	13.29	—
H13176	H13178	-0816	203 x 406mm (8 x 16")	14.70	23.25
H13176	—	-0818	203 x 457mm (8 x 18")	15.00	—
H13176	H13178	-1012	254 x 305mm (10 x 12")	12.99	23.70
H13176	—	-1014	254 x 356mm (10 x 14")	13.20	—
H13176	H13178	-1016	254 x 406mm (10 x 16")	16.35	27.75
—	H13178	-1018	254 x 457mm (10 x 18")	—	31.95
—	H13178	-1214	305 x 356mm (12 x 14")	—	26.25
—	H13178	-1216	305 x 406mm (12 x 16")	—	29.49
H13176	H13178	-1218	305 x 457mm (12 x 18")	16.95	32.49
H13176	H13178	-1220	305 x 508mm (12 x 20")	19.95	35.25

Great for Sample Collection and Storage!
Other Bag Sizes Available. Contact Us for Quotation.