



Graduated Beakers

"No-Drip" Spout for Easy Pouring


Griffin low form polypropylene beakers are translucent, chemical resistant and will stack without sticking together. Raised numeric graduations show approximate volume ($\pm 5\%$). All sizes are steam autoclavable at 121°C (250°F). **Bag and case quantities listed below.**



CATALOG NO.	CAPACITY	CASE	PER BAG	PRICE/BAG
F26210-0000	50ml	4	12	\$21.25
F26211-0000	100ml	4	12	20.00
F26212-0000	150ml	4	12	31.00
F26213-0000	250ml	6	6	18.00
F26214-0000	400ml	6	6	18.25
F26215-0000	600ml	6	4	19.25
F26216-0000	1000ml	4	3	14.25
F26217-0000	2000ml	6	1	11.75
F26219-0000	4000ml	6	1	22.75

Beaker Specifications

CATALOG NO.	GRADUATION	HEIGHT	TOP OD.	BOTTOM OD.
F26210-0000	Every 10ml	48.26mm (1.9")	50.80mm (2")	40.64mm (1.6")
F26211-0000	Every 25ml	66.04mm (2.6")	64.26mm (2.53")	50.80mm (2")
F26212-0000	Every 25ml	76.20mm (3")	73.15mm (2.88")	58.42mm (2.3")
F26213-0000	Every 50ml	88.90mm (3.5")	83.31mm (3.28")	66.80mm (2.63")
F26214-0000	Every 100ml	111.13mm (4.37")	96.52mm (3.8")	74.68mm (2.94")
F26215-0000	Every 100ml	127mm (5")	111mm (4.37")	83.82mm (3.3")
F26216-0000	Every 200ml	152mm (6")	130.18mm (5.125")	101.60mm (4")
F26217-0000	Every 500ml	184.15mm (7.25")	165.10mm (6.5")	133.35mm (5.25")
F26219-0000	Every 1000ml	234.95mm (9.25")	212.60mm (8.37")	171.45mm (6.75")



This "spec" chart answers our most commonly asked beaker questions.

Scienceware® Graduated Pitchers Are Also Ideal for Pouring Large Volumes, See Page 101



Large Volume Beakers

5 and 10 Liter Capacity Ideal for Pilot Plant Work

These non-graduated beakers feature flared rims as well as spouts for quick, easy pouring when large volumes are handled. Made of polypropylene, they offer good chemical resistance and are steam autoclavable at 121°C (250°F).



CATALOG NO.	CAPACITY	BOTTOM OD.	TOP OD.	HEIGHT	PRICE/EACH
F26219-0005	5 liter	16.5cm (6½")	21cm (8¼")	26cm (10¼")	\$49.50
F26219-0010	10 liter	22.2cm (8¾")	26.6cm (10½")	31.1cm (12¼")	91.50

Tubing Holder

Clips on Beakers and Pitchers

Hold flexible tubing steady in beakers or pitchers with this polypropylene clip. It fits securely onto the edge and has a 51mm (2") long arm with a hole to accept tubing up to 9.5mm (¾") O.D. Autoclavable. Holders fit onto beakers or pitchers up to ¼" thick. **6 per bag, 6 bags per case.**



CATALOG NO.	PRICE/BAG
F18235-0000	\$12.00

Dialysis Bag Holder

Suspend dialysis bags in large beakers with these strong clips. Just slide the bag into the slot and clip onto the edge of the beaker. The molded polypropylene surface may be marked with identifying sequence numbers. Holders fit onto beakers or pitchers up to ¼" thick. Autoclavable. **6 per bag, 6 bags per case.**



CATALOG NO.	PRICE/BAG
F18237-0000	\$12.00

Thermometer Holder

Protect Expensive Thermometers

This convenient holder clips onto the edge of a beaker or other vessel and holds a thermometer vertically near the side of the container for safe, easy reading. The thermometer is held in place by small "fingers" that will expand to accommodate thermometers ranging from 3.2mm (¼") to 11mm (7/16") in diameter. Made from polypropylene, it withstands temperatures up to 121°C (250°F). Autoclave. Holders fit onto beakers or pitchers up to ¼" thick. **3 per bag, 12 bags per case.**



CATALOG NO.	PRICE/BAG
F18241-0000	\$9.55

For Other Scienceware® Clips, See Our Clamps and Holders Section on Pages 21-23



Weighted Beaker/Flask Holders

Eliminate Drifting and Tipping in Water Baths

Weighted holder grips the beaker or flask firmly in four stainless steel prongs flared out for easy use. The Vikem™ vinyl-coated lead ring secured beneath the holder keeps the beaker immersed in water, provides a low center of gravity and allows water to circulate underneath the beaker. Can be used in boiling liquid up to 135°C(275°F) for 1 hour.



CATALOG NO.	FITS BEAKER	WEIGHT	DIAMETER	PRICE/EACH
F18309-0100	100ml	165gm (.36lb)	60mm (2.4")	\$50.00
F18309-0250	250ml	237gm (.5lb)	80mm (3.3")	52.50
F18309-0600	600ml	454gm (1lb)	95mm (3.7")	56.00