

NEW



Precisionware® Polypropylene Wide Mouth Bottles - Autoclavable

Made of FDA Compliant Material

Choose Precisionware® Polypropylene Wide Mouth Bottles for pharmaceutical or similar packaging of both liquid and dry materials. Bottles are Boston round in style, water-clear and molded from FDA compliant material. Matched polypropylene caps have a valve-seal design for a perfect, leak proof fit. Deep square threads withstand vibration and help prevent closure back-off during shipping. Ribbed caps are designed for torque-closure and easy opening. Inner bottle surfaces are smooth and rounded for easy cleaning and thorough emptying. Underside of bottle lists volume, material and includes an indexing notch. Recycling codes listed on 500 and 1000 ml bottles. A protruding ring under the cap holds shrink wrapping securely in place. Autoclavable 121°C (250°F) *Do not autoclave with caps in place. 12 per bag except 1000ml, 6 per bag.*



CATALOG NO.	CAPACITY	CLOSURE SIZE	CASE	PRICE/BAG
F10632-0003	30ml (1 oz)	28mm	6 Bags	\$12.00
F10632-0004	60ml (2 oz)	28mm	6 Bags	12.50
F10632-0005	125ml (4 oz)	38mm	6 Bags	17.50
F10632-0006	250ml (8oz)	43mm	6 Bags	26.00
F10632-0007	500ml (16oz)	53mm	4 Bags	37.00
F10632-0008	1000ml (32oz)	63mm	4 Bags	32.00



Precisionware® Low-Density Polyethylene Wide Mouth Bottles

Wide Mouth Bottles for Easy Filling and Pouring

Translucent, leak proof design, flexible and easy to fill wide mouth bottles. Made of FDA compliant material, ideal for general storage of liquids and powders. Low-density polyethylene with matched polypropylene closures. Bottles 250ml (8oz) and larger, feature recycling code molded into base. *12 per bag except 1000ml, 6 per bag.*



CATALOG NO.	CAPACITY	CLOSURE SIZE	CASE	PRICE/BAG
F10626-0005	125ml (4oz)	38mm	6 Bags	\$16.25
F10626-0006	250ml (8oz)	45mm	6 Bags	21.65
F10626-0007	500ml (16oz)	53mm	4 Bags	30.65
F10626-0018	1000ml (32oz)	53mm	4 Bags	28.00



Precisionware® High-Density Polyethylene Wide Mouth Bottles

Wide Mouth Bottles for Easy Filling and Pouring

Rigid leak proof design with easy to fill, wide mouth. These high-density polyethylene bottles with matched polypropylene closures and made of FDA compliant material, are ideal for packaging and storing of liquids and powders. Bottles 250ml (8oz) and larger, feature recycling code molded into base. *12 per bag except 1000ml, 6 per bag.*



CATALOG NO.	CAPACITY	CLOSURE SIZE	CASE	PRICE/BAG
F10625-0005	125ml (4oz)	38mm	6 Bags	\$14.95
F10625-0006	250ml (8oz)	45mm	6 Bags	21.25
F10625-0007	500ml (16oz)	53mm	4 Bags	30.65
F10625-0018	1000ml (32oz)	53mm	4 Bags	27.95



Wide Mouth Round Bottles - Heavy Duty Closure

These popular bottles offer proven dependability for storing, shipping, and handling materials in both liquid and powder form. Molded of low-density polyethylene, they feature a heavy threaded, self sealing screw cap closure. **Note: the 5,000 and 10,000ml sizes are made of high-density polyethylene and have an octagonal shape. These bottles are not leakproof and are best used for dry materials.*



CATALOG NO.	CAPACITY	DIAMETER	HEIGHT	NECK I.D.	PRICE/BAG
F10906-0050	50ml	39mm	78mm	24mm	\$26.50/12
F10906-0100	100ml	48mm	100mm	24mm	30.00/12
F10906-0200	200ml	59mm	115mm	31mm	38.50/12
F10906-0250	250ml	64mm	127mm	31mm	39.50/12
F10906-0500	500ml	77mm	160mm	39mm	30.00/6
F10906-1000	1000ml	95mm	210mm	50mm	48.00/6
F10906-1500	1500ml	108mm	230mm	50mm	36.00/3
F10906-2000	2000ml	120mm	255mm	50mm	37.00/3
F10906-5000*	5000ml	192mm	235mm	110mm	42.00/EACH
F10906-0010*	10,000ml	232mm	310mm	110mm	45.00/EACH



Buttruss Cap Bottles

Low-density polyethylene bottles with water tight buttruss screw caps. PE caps have a 38/430 finish. *12 per bag.*



CATALOG NO.	CAPACITY	DIAMETER	HEIGHT	CAP	CASE	PRICE/BAG
F10611-0016	500ml (16oz)	79mm (3.10")	153mm (6.02")	38mm	3 Bags	\$33.50
F10611-0032	1000ml (32oz)	94mm (3.71")	220mm (8.66")	38mm	2 Bags	48.95



Small Bottle

This small polyethylene bottle with screw cap is ideal for containment of critical quantities. *12 per bag, 6 bags per case.*



CATALOG NO.	CAPACITY	DIAMETER	HEIGHT	CAP	PRICE/BAG
F10605-0005	5.0ml (0.167oz)	22mm (.86")	41mm (1.60")	15mm	\$20.60



Square Bottles

Square shaped, leak proof design and chemical resistant. Tamper-proof tear-lock closure. These bottles with heavy-threaded necks make good use of limited shelf space and are ideal for storing or shipping samples. They are molded of high-density polyethylene with deep-shelled, blue polypropylene closures. *6 per bag.*



CATALOG NO.	NECK CAPACITY	WIDTH	HEIGHT	OPENING	PRICE/BAG
F10904-0500	500ml	75mm	158mm	43mm	\$34.00
F10904-1000	1000ml	85mm	207mm	43mm	43.00

Square Bottles Help to Save Valuable Work Space