



PRODUCT SPECIFICATION SHEET
BUTYL ACETATE
REAGENT ACS GRADE
Catalog No: 304000ACS

Available in 55 gallon Drums only

Test	Specification	Typical Lot Analysis
Assay	NLT 99.5%	99.84%
Titration Acid	0.0016meq/g	<0.0016meq/g
Specific Gravity @ 20/20	0.882 - 0.884	0.883
Residue After Evaporation	NMT 0.001%	<0.001%
Color	NMT 10	5
Substance Darkened by H ₂ SO ₄	To pass test	Pass
Water	NMT 0.1%	0.0278%
n-Butyl Alcohol	NMT 0.2%	0.11 %
Acetic Acid	0.01%, max.	Pass
n-Butyl Formate	NMT 0.1%	Pass
n-Butyl Propionate	NMT 0.1%	Pass

Form Butyl Acetate-ACS, # 101, Rev. 2.0, 08/01

Approved by: Paul DiMarco, Vice President Quality Management

FOR INDUSTRIAL OR LAB USE ONLY. Not intended as an active substance in Food or Drug. Not to be considered a Medical Device. Not intended for use as a Disinfectant as defined by the EPA. The expiration date of this product is three years from the date of manufacture. (Rev. # Disclaimer only: Rev 2.0 11/19/01 PJM)

This document has been electronically generated and is therefore not signed. A signed copy can be obtained from Pharmco's Quality Control office.

58 Vale Road, Brookfield, CT 06804, USA
Telephone: (203) 740-3471 Fax (203) 740-3481
E-Mail: paul@pharmco-prod.com
Web Site: www.pharmco-prod.com



Approved by: Paul DiMarco, Vice President Quality Management

FOR INDUSTRIAL OR LAB USE ONLY. Not intended as an active substance in Food or Drug. Not to be considered a Medical Device. Not intended for use as a Disinfectant as defined by the EPA. The expiration date of this product is three years from the date of manufacture. (Rev. # Disclaimer only: Rev 2.0 11/19/01 PJM)

This document has been electronically generated and is therefore not signed. A signed copy can be obtained from Pharmco's Quality Control office.

58 Vale Road, Brookfield, CT 06804, USA
Telephone: (203) 740-3471 Fax (203) 740-3481
E-Mail: paul@pharmco-prod.com
Web Site: www.pharmco-prod.com