



PRODUCT SPECIFICATION SHEET

TETRAHYDROFURAN
REAGENT GRADE ACS
Catalog Number: 346000ACS
Available in the following sizes:

4x4 liter Amber Glass, 5 gallon Metal w/Rieke Spout, 55 gallon New Metal Drum

PRODUCT SPECIFICATIONS	LIMITS	TYPICAL LOT ANALYSIS
Assay (GC), min	99.0%	99.99%
Inhibitor (BHT)	250-400 ppm	305ppm
Color (APHA), max	20	5
Peroxide (as H ₂ O ₂)	0.015%	0.002%
Residue After Evaporation	0.03%	<10 ppm
Water, max	0.05%	0.0037%
Specific Gravity	0.886-0.889	0.888
Max Total Organic Impurities	0.08%	<0.007%


Form: Tetrahydrofuran,ACS,#201,Q517, Rev. 2.0

FOR INDUSTRIAL OR LAB USE ONLY. Not intended as an active substance in Food or Drug. Not to be considered a Medical Device. Not intended for use as a Disinfectant as defined by the EPA. The expiration date of this product is three years from the date of manufacture. (Rev. # Disclaimer only: Rev 2.0 11/19/01 PJM)

Approved by: Paul DiMarco, Vice President Quality Management

This document has been electronically generated and is therefore not signed. A signed copy can be obtained from Pharmco's Quality Control office.

58 Vale Road, Brookfield, CT 06804, USA
Telephone: (203) 740-3471 Fax (203) 740-3481
E-Mail: paul@pharmco-prod.com
Web Site: www.pharmco-prod.com



Approved by: Paul DiMarco, Vice President Quality Management

This document has been electronically generated and is therefore not signed. A signed copy can be obtained from Pharmco's Quality Control office.

58 Vale Road, Brookfield, CT 06804, USA
Telephone: (203) 740-3471 Fax (203) 740-3481
E-Mail: paul@pharmco-prod.com
Web Site: www.pharmco-prod.com