



Product	Grade Certifications	C6H5CH3	CAS 108-88-3	FW 92.14
---------	----------------------	---------	--------------	----------

Toluene - Lab Grade - Hospital Clinical Grade

<u>Catalog Number</u>	<u>Package Sizes</u>	<u>Product Specifications</u>	<u>Limits</u>
347000LAB		Assay (GC), min	99.5%
		Water, max	0.10%
	4 x 1 gallon Poly	Sulfur Compounds, max	0.15%
	5 gallon Metal	Appearance	Clear
	5 gallon POLY	Residue After Evaporation, max	10 ppm
	55 gallon Metal	Boiling Range, within 2C, inc.	Passes Test
		110.6C +/- 0.1C	
		Color	Clear
7.26 lbs./gallon	700 lbs. net	Not Intended For Use As A Medical Device	

WHEN ORDERING THIS PRODUCT THE ACS GRADE AND LABEL IS SUBSTITUTED

Approved by Paul DiMarco, Vice President Quality Management

58 Vale Road, Brookfield, CT 06804, USA
Telephone: (203) 740-3471 ext. 134, 800-243-5360 ext. 134, Fax (203) 740-3481
E-Mail: paul@pharmco-prod.com
Website: www.pharmco-prod.com