

2005 List Price



UNICO 1000 Series Spectrophotometer

Item #	Description	List Price
S-1000	UNICO Model 1000 Spectrophotometer 20 nm Bandpass complete with 10 mm Test Tube Cuvettes (Box of 12), 10 mm Square Cuvette Adapter, Analog Output, Dust Cover, User Manual. Wavelength Range of 400-1000 nm. Voltage preset at 110V	\$695.00
S-1000E	Same as Model 1000 but Voltage preset at 220V	\$695.00
Accessories		
S-1100-101	Experiment Manual includes safety in the lab, 10 experiments, & instructor's guide	\$40.00
S-1100-102	Test Tube Sample Holder for 10 mm diameter test tubes (Standard with instrument)	\$35.00
S-1100-103	Square Cuvette Adapter for 10 mm cuvettes (Standard with instrument)	\$10.00
S-1100-113	Test Tube Sample Holder for 3/4 inch test tubes	\$35.00
S-1100-114	Test Tube Sample Holder for 1 inch test tubes	\$35.00
S-1100-115	COD Vial Sample Holder	\$35.00
S-1100-220	Power cord, European plug	\$6.00
S-2100-115	Holmium Oxide Filter, 10 mm x 10 mm x 45 mm (Requires S-1100-103)	\$60.00
S-2100-116	Didymium Filter, 10 mm x 10 mm x 45 mm (Requires S-1100-103)	\$60.00
Glassware		
S-90-301	Test Tube Cuvettes, 10mm diameter, Box of 12	\$35.00
S-90-302P-100	Disposable Cuvettes, Polystyrene, 10 mm pathlength, Box of 100 pcs	\$30.00
S-90-302P-500	Disposable Cuvettes, Polystyrene, 10 mm pathlength, Case of 500 pcs	\$80.00
S-90-304G	Square Cuvettes, Optical Glass, Set of 2	\$40.00
S-90-305P	Disposable semi-micro Cuvettes, Polystyrene, 10 mm pathlength, 500 pcs	\$100.00
Replacement Items		
S-1100-505	Tungsten-Halogen Lamps, Package of 1 (6V/10W, G4 pin)	\$19.00
S-1000-511	User Manual	\$8.00
S-1100-512	Dust Cover	\$14.00
S-1100-513	Fuse, 3 A, Package of 1	\$5.00
S-1100-521	Allen Wrench	\$4.00

Actual Size: 406 mm X 305 mm X 203 mm

Net Weight: 14 lbs.

Shipping Size: 496 mm X 410 mm X 286 mm

Shipping Weight: 16 lbs.



Model 1000

2005 List Price

UNICO®**UNICO 1100 Series Spectrophotometer**

Item #	Previous Item #	Description	List Price
S-1100	1100	UNICO Model 1100 Spectrophotometer 20nm Bandpass Complete with 10 mm Test Tube Cuvettes (Box of 12), 10mm Square Cuvette Adapter, Analog Output, Dust Cover, User Manual Wavelength Range of 335-1000 nm. Voltage preset at 110V	\$795.00
S-1100E	1100-E	Same as Model 1100 but preset at 220V	\$795.00
S-1100RS	1100RS	UNICO Model 1100RS Spectrophotometer 10 nm Bandpass Complete with 10 mm Test Tube Cuvettes (12pcs), 10 mm Square Cuvette Adapter, RS-232C Serial Port, Dust Cover, User Manual Wavelength range of 335-1000 nm. Voltage Preset at 110V	\$995.00
S-1100RSE	1100RS-E	Same as Model 1100RS but preset at 220V	\$995.00
Software			
S-1100-401	1100-401	UNICO Application Software for PC's, Window's 95, 98, 2000, Windows ME or Windows XP required. Programs include Standard Curve, Abs.%T/Conc., and Abs.vs.Time. Includes S-2100-226 RS232 serial cable. For 1100RS ONLY	\$150.00
Accessories			
S-1100-101	1100-101	Experiment Manual includes safety in the lab, 10 experiments, & instructor's guide	\$40.00
S-1100-102	1100-102	Test Tube Sample Holder for 10 mm diameter test tubes (Standard with instrument)	\$35.00
S-1100-103	1100-103	Square Cuvette Adapter for 10 mm cuvettes (Standard with instrument)	\$10.00
S-1100-113	1100-113	Test Tube Sample Holder for 3/4 inch test tubes	\$35.00
S-1100-114	1100-114	Test Tube Sample Holder for 1 inch test tubes	\$35.00
S-1100-115	1100-115	COD Vial Sample Holder	\$35.00
S-2100-115	1100-109	Holmium Oxide Filter, 10mm x 10mm x 45mm (Requires S-1100-103)	\$60.00
S-2100-116	1100-110	Didymium Filter, 10mm x 10mm x 45mm (Requires S-1100-103)	\$60.00
Output Devices			
S-1100-206	1100-206	RS-232 Printer (Requires S-1100-207 Cable)	\$390.00
S-1100-207	1100-207	Printer Cable for use with RS-232 printer (Male 25-pin to Female 9-pin, Null modem)	\$33.00
S-1100-220		Power cord, European plug	\$6.00
S-2100-226	1100-226	Computer RS-232 Serial Cable for use with software (9-pin to 9-pin, Female/Female, Null Modem)	\$25.00
Glassware			
S-90-301	1100-301	Test Tube Cuvettes, 10 mm diameter, Box of 12	\$35.00
S-90-302P-100	1100-302	Disposable Cuvettes, Polystyrene, 10 mm pathlength, Box of 100	\$30.00
S-90-302P-500	1100-303	Disposable Cuvettes, Polystyrene, 10 mm pathlength, Case of 500	\$80.00
S-90-304G	1100-304	Square Cuvettes, Optical Glass, Set of 2	\$40.00
S-90-305P	1100-305	Disposable semi-micro Cuvettes, Polystyrene, 10 mm pathlength, 500pcs.	\$100.00
Replacement Items			
S-1100-505	1100-505	Tungsten-Halogen Lamps, Package of 1 (6V/10W, g4 pin)	\$19.00
S-1100-508	1100-508	Printer Paper, Package of 3	\$15.00
S-1100-511	1100-511	User Manual	\$8.00
S-1100-512	1100-512	Dust Cover	\$14.00
S-1100-513	1100-513	Fuse, 3 A, Package of 1	\$5.00
S-1100-520	1100-520	Service Manual	\$80.00
S-1100-521	1100-521	Allen Wrench	\$4.00

**Model 1100**

Actual Size: 408 mm X 308 mm X 185 mm

Net Weight: 13 lbs.

Shipping Size: 496 mm X 410 mm X 286 mm

Shipping Weight: 16 lbs.

**Model 1100RS**

2005 List Price**UNICO 1200 Series Spectrophotometer**

Item #	Previous Item #	Description	List Price
S-1200	1200	UNICO Model 1200 Spectrophotometer 5nm Bandpass Complete with 10 mm single Square Cuvette Holder, V-type Holder for 5~10 mm round tube, a pair of Glass Cuvettes, 12 pcs of 10 mm Test Tubes RS-232C Port, Dust Cover, User Manual Wavelength range: 325~1000 nm. Voltage preset at 110V	\$1,195.00
S-1200-E	1200-E	Same as Model 1200 but preset at 220V	\$1,195.00
S-1201	1201	UNICO Model 1200 Spectrophotometer 5nm Bandpass Complete with 4-position Cell Holder for 10 mm Square Cuvettes, a pair of Glass Cuvettes, RS-232C Port, Dust Cover, User manual Wavelength 325~1000 nm. Voltage preset at 110V	\$1,195.00
S1201-E	1201-E	Same as Model 1201 but preset at 220V	\$1,195.00
Software			
S-1200-401	1200-401	UNICO Application Software for PC's, Window's 95 or 98, 2000, Windows ME Windows XP required. Programs include Standard Curve, Abs.%T/Conc., and Abs.vs.Time. Includes S-2200-226 RS-232 serial cable.	\$150.00
Accessories			
S-1200-101	1200-101	Test Tube V-type Cell Holder with base for 8~25 mm diameter test tubes	\$80.00
S-1200-102	1200-102	Single Square Cuvette Holder with base for 10 mm Cuvettes	\$80.00
S-1200-103	1200-103	4-position Cell Holder for 10 mm square cuvettes (Standard in model 1201)	\$110.00
S-1200-104	1200-104	Longpath 4-Cell holder for up to 50 mm rectangular cuvettes	\$120.00
S-1200-105	1200-105	Water-jacketed cell holder with base and panel for temperature controlled application	\$130.00
S-1200-106	1200-106	Universal Base for V-type holder, single cell holder and water-jacketed holder	\$40.00
S-1200-111	1200-111	Black Block to set 0%T	\$10.00
S-2100-115	1100-109	Holmium Oxide Filter, 10mmx10mmx45mm	\$60.00
S-2100-116	1100-110	Didymium Filter, 10mmx10mmx45mm	\$60.00
Output Devices			
S-1100-206	1200-206	RS-232 Printer (Requires S-1100-207 Cable)	\$390.00
S-1100-207	1200-207	Printer Cable for use with RS-232 printer (Male 25-pin to Female 9-pin, Null modem)	\$33.00
S-1200-220		Power cord, European plug	\$6.00
S-2100-226	1200-226	Computer RS-232C Cable for use with software (9-pin to 9-pin, Female/Female, Null Modem)	\$25.00
Glassware			
S-90-301	1100-301	Test Tube Cuvettes, 10 mm diameter, Box of 12	\$35.00
S-90-302P-100	1100-302	Disposable Cuvettes, Polystyrene, 10 mm pathlength, Box of 100	\$30.00
S-90-302P-500	1100-303	Disposable Cuvettes, Polystyrene, 10 mm pathlength, Case of 500	\$80.00
S-90-304G	2100-304G	Square Cuvettes, Optical Glass, Set of 2	\$40.00
S-90-305P	1100-305	Disposable Semi-Micro Cuvettes, Polystyrene, 10 mm pathlength, 500pcs.	\$100.00
Replacement Items			
S-1100-101	1100-101	Experiment Manual includes safety in the lab, 10 experiments, & instructor's guide	\$40.00
S-1200-505	1200-505	Tungsten-Halogen Lamps, Package of 1 (6V/10W, G4 pin)	\$19.00
S-1100-508	1100-508	Printer Paper, Package of 3	\$15.00
S-1200-511	1200-511	User Manual	\$8.00
S-1200-512	1200-512	Dust Cover	\$14.00
S-1200-513	1200-513	Fuse, 3 A, Package of 1	\$5.00
S-1200-520	1200-520	Service Manual	\$80.00

**Model 1200**

Actual Size: 408 mm X 308 mm X 180 mm

Net Weight: 15 lbs.

Shipping Size: 496 mm X 410 mm X 286 mm

Shipping Weight: 20 lbs.

2005 List Price



UNICO 2100 Series Spectrophotometer

Page 1-2

Item #	Previous Item #	Description	List Price
S-2100	2100	UNICO Model 2100 Spectrophotometer 5 nm Bandpass, Complete with 4-position Cell Holder, 4 Optical Glass Cuvettes, RS-232C Port, Dust Cover, User Manual Wavelength Range: 325~1000 nm. Automatic Wavelength Change. Voltage preset at 110V	\$1,350.00
S-2100-S		Same as Model S-2100 with scanning feature and bundled software, scan at 1 nm interval	\$1,650.00
S-2100-E	2100-E	Same as Model 2100 but preset at 220V	\$1,350.00
S-2100-SE		Same as Model S-2100-S but preset at 220V	\$1,650.00
S-2100UV	2100UV	UNICO Model 2100UV Spectrophotometer 5 nm Bandpass, Complete with 4-position Cell Holder, 4 Optical Glass Cuvettes, 2 Quartz Cuvettes, RS-232C Port, Dust Cover, User Manual Wavelength Range: 200~1000 nm, Automatic Wavelength Change, Voltage preset at 110V	\$3,250.00
S-2100UV-S		Same as Model S-2100UV with scanning feature and bundled software, scan at 1 nm interval	\$3,550.00
S-2100UV-E	2100UV-E	Same as Model 2100UV but preset at 220V	\$3,250.00
S-2100UV-SE		Same as Model S-2100UV-S but preset at 220V	\$3,550.00
Software			
S-2100-401	2100-401	UNICO Application Software for PC's. Window's 95, 98,2000, ME or XP required Programs include Standard Curve, Abs/%T/Conc., and Abs. Vs. Time. Includes CD software, manual and item # S-2100-226 serial cable. (for Model 2100 and 2100UV)	\$150.00
Accessory Packages			
S-2100-101P	2100-101P	Test Tube Holder Kit for 8~20 mm diameter test tubes. Includes universal base, one V-type tube holder and top hat for test tubes	\$90.00
S-2100-102P	2100-102P	Rectangular Long Path Cell Holder Kit for cells up to 100 mm pathlength. Includes universal base and one holder	\$75.00
S-2100-103P	2100-103P	Single 10 mm Square Cuvette Holder Kit. Includes universal base and square cuvette holder	\$75.00
S-2100-104P	2100-104P	Cylindrical Cell Holder Kit for cells up to 100 mm pathlength (22mm dia.). Includes universal base and one holder	\$75.00
S-2100-105P	2100-105P	Water-Jacketed Cuvette Holder Kit for 10 mm single square cuvette Includes universal base, water-jacketed holder with flow-thru panel (Water Bath required)	\$150.00
S-2100-106P	2100-106P	Micro-Cell Holder Kit for use with 10 mm microcuvettes. Includes adjustable x-y base and single cuvette holder for 10 mm square cuvette	\$150.00
S-2100-107P	2100-107P	Peltier Package for continuous temperature control: 15°C~40°C. Includes item # S-2100-107A peltier thermal electric 10 mm cell holder and item # S-2100-107B peltier controller unit	\$1,225.00
S-2100-108P	2100-108P	Flow-thru Package for ambient flow through. Includes item # S-2100-109 flow-thru controller and item # S-2100-108B flow-thru panel (Requires flowcell and tubing)	\$1,100.00
S-2100-109P	2100-109P	Peltier/Sipper Package for flow-thru and continuous temperature control: 15°C~40°C Includes item # S-2100-109 peltier/flow-thru controller and item # S-2100-107A peltier cell holder with panel (Requires flowcell and tubing)	\$1,325.00
Components and Accessories			
S-2100-101	2100-101	Test Tube Holder without base (Requires Universal Base)	\$35.00
S-2100-102	2100-102	Rectangular Long Path Cell Holder without base (Requires Universal Base)	\$35.00
S-2100-103	2100-103	Single 10 mm Cuvette Holder without base (Requires Universal Base)	\$35.00
S-2100-104	2100-104	Cylindrical Cell Holder without base (Requires Universal Base)	\$35.00
S-2100-105	2100-105	Water-Jacketed Single Cuvette Holder with flow-thru panel but without base (Requires Universal Base, Water Bath)	\$100.00
S-2100-107A	2100-107A	Peltier thermal electric 10 mm Cell Holder with panel	\$275.00
S-2100-107B	2100-107B	Peltier controller unit for controlling temperature: 15°C~40°C	\$950.00
S-2100-108A	2100-110	Flow-thru Panel for use with ambient flowcell system (Requires flowcell, tubing and item # S-2100-109 flow-thru/peltier controller)	\$50.00
S-2100-109	2100-109	Peltier Thermal Temperature / Flow-thru Controller Only	\$1,050.00
S-2100-111	2100-111	4-position Manual Cell Holder (Standard with the instrument)	\$150.00
S-2100-112	2100-112	Universal Base. Accommodates up to 2 same or different type of holders	\$60.00

continued on next page

2005 List Price



UNICO 2100 Series Spectrophotometer

Page 2-2

Item #	Previous Item #	Components and Accessories (continued)	List Price
S-2100-113	2100-113	Top Hat for test tubes with height of 75~150 mm	\$25.00
S-2100-114	2100-114	Spacers for 5 mm cuvettes	\$45.00
S-2100-115	2100-115	Holmium Oxide Filter	\$60.00
S-2100-116	2100-116	Didymium Filter	\$60.00
S-2100-120	2100-120	Haake Model DC10-B3 water bath for temperature: 25°C~100°C for use with S-2100-105	please call
S-2100-121	2100-121	Haake Model DC10-K10 water bath for temperature: -10°C~100°C for use with S-2100-105	please call
Output Devices			
S-1100-206	2100-206	RS-232 Printer (Requires item # S-1100-207 printer cable)	\$390.00
S-1100-207	2100-207	Printer Cable for use with RS-232 printer (Male 25-pin to Female 9-pin, Null Modem)	\$33.00
S-2100-220		Power cord, European plug	\$6.00
S-2100-226	2100-226	Computer RS-232 Cable for use with application software (9-pin to 9-pin Female/Female, Null Modem)	\$25.00
Glassware			
S-90-301	2100-301	Test Tube Cuvettes, 10 mm Diameter, Box of 12	\$35.00
S-90-302P-100	2100-302P	Disposable Cuvettes, Polystyrene, 10 mm Pathlength, Box of 100	\$30.00
S-90-302P-500	2100-303P	Disposable Cuvettes, Polystyrene, 10 mm Pathlength, Case of 500	\$80.00
S-90-305P	2100-305P	Disposable Cuvettes, Polystyrene, 10 mm Pathlength, Case of 500	\$100.00
S-90-304G	2100-304G	Square Cuvettes, Optical Glass, Set of Two	\$40.00
S-90-309Q	2100-309Q	Square Cuvettes, Quartz, 10 mm, Set of Two	\$100.00
S-90-308-50G	2100-308-50G	Cylindrical Cell, 50 mm pathlength	\$95.00
S-90-326G	2100-326G	Rectangular Cell, Optical Glass, 50 mm pathlength, with stopper	\$89.00
S-90-344FG	2100-344FG	Flowcell, Optical Glass, 4x12 mm window, 0.48 ml	\$290.00
S-90-346FQ	2100-346FQ	Flowcell, Quartz, 4x12mm window, 0.48 ml	\$400.00
S-90-350G	2100-350G	Square Cuvettes, Optical Glass, 4 mm width 1.0 ml	\$42.00
S-90-351Q	2100-351Q	Square Cuvettes, Quartz, 4 mm width, 1.0 ml	\$100.00
S-90-352Q	2100-352Q	Square Cuvettes, Quartz, 4 mm width, with stopper, 0.5 ml	\$130.00
S-90-353G	2100-353G	Square Cuvettes, Optical glass, 2 mm width, 0.5 ml	\$55.00
S-90-354Q	2100-354Q	Square Cuvettes, Quartz, 2 mm width, 0.5 ml	\$135.00
S-90-355Q	2100-355Q	Square Cuvettes, Quartz, 2 mm width, with stopper, 0.5 ml	\$155.00
S-90-358Q	2100-358Q	Square Cuvettes, Quartz, 2x5 mm window, 0.100 ml	\$195.00
S-90-359Q	2100-359Q	Square Cuvettes, Quartz, 2x2.5 mm window, 0.050 ml	\$195.00
For more cuvettes and glassware, please refer to "Glassware and Cuvettes" Price Sheet			
Miscellaneous and Replacement Items			
S-1100-508	2100-508	Printer Paper, Package of 3, for use with RS-232 printer	\$15.00
S-2100-510	2100-510	User Manual	\$8.00
S-2100-511	2100-511	Dust Cover	\$14.00
S-2100-515	2100-515	Tungsten-Halogen Lamp for Model 2100 and 2100UV, Package of 1 (6V/10W, G4 pin)	\$19.00
S-2100-525	2100-525	Deuterium Lamp for Model 2100UV	\$290.00
S-2100-535	2100-535	Fuse, 3 A, Package of 1	\$5.00



Model 2100UV

Actual Size: 18.5 in x 15.7 in x 5.5 in (470 mm x 400 mm x 140 mm)
 Net Weight: 26 lbs/32 lbs. (Model 2100/2100UV)
 Shipping Size: 22.4 in x 19.7 in x 11.4 in (570 mm x 500 mm x 290 mm)
 Shipping Weight: 29/35 lbs. (Model 2100/2100UV)

2005 List Price



Glassware and Cuvettes for UNICO Spectrophotometers

Item #	Previous Item #	Materials	Shape	Pathlength	Interior Width	Height	Capacity	Wavelength Range	Qty	List Price
Regular or Macro Cuvettes (10 mm pathlength)										
S-90-304-1G	2100-304-1G	Glass	Square	10 mm	10 mm	45 mm	3.5 ml	335 to 900 nm	1	\$20.00
S-90-304G	2100-304G	Glass	Square	10 mm	10 mm	45 mm	3.5 ml	335 to 900 nm	2	\$40.00
S-90-309-1Q	2100-309-1Q	Quartz	Square	10 mm	10 mm	45 mm	3.5 ml	200 to 1200 nm	1	\$50.00
S-90-309Q	2100-309Q	Quartz	Square	10 mm	10 mm	45 mm	3.5 ml	200 to 1200 nm	2	\$100.00
S-90-301	2100-301	Glass	Tube	10 mm	10 mm	100 mm	4.0 ml	335 to 1100 nm	12	\$35.00
S-90-302P-100	2100-302P	Polystyr.	Square	10 mm	10 mm	45 mm	3.5 ml	335 to 1100 nm	100	\$30.00
S-90-302P-500	2100-303P	Polystyr.	Square	10 mm	10 mm	45 mm	3.5 ml	335 to 1100 nm	500	\$80.00
S-90-305P	2100-305P	Polystyr.	Square	10 mm	4 mm	45 mm	1.5 ml	335 to 1100 nm	500	\$100.00
Non-10 mm pathlength Cuvettes										
S-90-308-50G	2100-308-50G	Glass	Cylindr.	50 mm	n/a	22 mm	14.1 ml	335 to 2500 nm	1	\$95.00
S-90-308-100G	2100-308-100G	Glass	Cylindr.	100 mm	n/a	22 mm	28.2 ml	335 to 2500 nm	1	\$130.00
S-90-310Q	2100-310Q	Quartz	Rectan.	5 mm	w/ stopper	48 mm	1.7 ml	170 nm-2500 nm	1	\$125.00
S-90-311Q	2100-311Q	Quartz	Rectan.	50 mm	w/ stopper	48 mm	17.5ml	170 nm-2500 nm	1	\$165.00
S-90-312Q	2100-312Q	Quartz	Rectan.	30 mm	w/ stopper	48 mm	10.5ml	170 nm-2500 nm	1	\$146.00
S-90-313Q	2100-313Q	Quartz	Rectan.	20 mm	w/ stopper	48 mm	7.0 ml	170 nm-2500 nm	1	\$125.00
S-90-314Q	2100-314Q	Quartz	Rectan.	100 mm	w/ stopper	48 mm	3.0 ml	170 nm-2500 nm	1	\$315.00
S-90-320G	2100-320G	Glass	Rectan.	20 mm	open top	45 mm	7.0 ml	335 to 2500 nm	1	\$40.00
S-90-321G	2100-321G	Glass	Rectan.	50 mm	open top	45 mm	17.5 ml	335 to 2500 nm	1	\$55.00
S-90-322G	2100-322G	Glass	Rectan.	100 mm	open top	45 mm	35.0 ml	335 to 2500 nm	1	\$100.00
S-90-325G	2100-325G	Glass	Rectan.	20 mm	w/ stopper	48 mm	7.0 ml	335 to 2500 nm	1	\$65.00
S-90-326G	2100-326G	Glass	Rectan.	50 mm	w/ stopper	48 mm	17.5 ml	335 to 2500 nm	1	\$89.00
S-90-327G	2100-327G	Glass	Rectan.	100 mm	w/ stopper	48 mm	35.0 ml	335 to 2500 nm	1	\$175.00
S-90-330G	2100-330G	Glass	Rectan.	1 mm	10 mm	45 mm	0.35 ml	335 - 1000 nm	1	\$90.00
S-90-331Q	2100-331Q	Quartz	Rectan.	1 mm	10 mm	45 mm	0.35 ml	170 - 2500 nm	1	\$165.00
S-90-332G	2100-332G	Glass	Rectan.	2 mm	10 mm	45 mm	0.70 ml	170 - 2500 nm	1	\$90.00
Flowcell, Macro and Micro										
S-90-340FG	2100-340FG	Glass	Square	Frosted	10x7x40 mm	48 mm	3.0 ml	320 to 2500 nm	1	\$280.00
S-90-341FG	2100-341FG	Glass	Square	Frosted	10x4x40 mm	48 mm	1.8 ml	320 to 2500 nm	1	\$295.00
S-90-342FQ	2100-342FQ	Quartz	Square	Frosted	10x7x40 mm	48 mm	3.0 ml	170 to 2500 nm	1	\$430.00
S-90-343Q	2100-343Q	Quartz	Square	Frosted	10x4x40 mm	48 mm	1.8 ml	170 to 2500 nm	1	\$475.00
S-90-344FG	2100-344FG	Glass	Square	Black	10x4x12 mm	48 mm	0.48 ml	320 to 2600 nm	1	\$290.00
S-90-345FG	2100-345FG	Glass	Square	Black	10x3 mm (dia.)	48 mm	0.07 ml	320 to 2500 nm	1	\$310.00
S-90-346FQ	2100-346FQ	Quartz	Square	Black	10x4x12 mm	48 mm	0.48 ml	170 to 2500 nm	1	\$400.00
S-90-347FQ	2100-347FQ	Quartz	Square	Black	10x3 mm (dia.)	48 mm	0.07 ml	170 to 2500 nm	1	\$415.00
Semi-Micro, Sub-micro Cuvettes										
S-90-350G	2100-350G	Glass	Lid	Black	10x4	45 mm	1.0 ml	320 to 2,600 nm	1	\$42.00
S-90-351Q	2100-351Q	Quartz	Lid	Black	10x4	45 mm	1.0 ml	170 to 2,600 nm	1	\$100.00
S-90-352Q	2100-352Q	Quartz	Stopper	Black	10x4	45 mm	1.0 ml	170 to 2,600 nm	1	\$130.00
S-90-353G	2100-353G	Glass	Lid	Black	10x2	45 mm	0.5 ml	320 to 2,600 nm	1	\$55.00
S-90-354Q	2100-354Q	Quartz	Lid	Black	10x2	45 mm	0.5 ml	170 to 2,600 nm	1	\$135.00
S-90-355Q	2100-355Q	Quartz	Stopper	Black	10x2	45 mm	0.5 ml	170 to 2,600 nm	1	\$155.00
S-90-356G	2100-356G	Glass	Lid	Black	10x1	45 mm	0.25 ml	320 to 2,600 nm	1	\$75.00
S-90-357Q	2100-357Q	Quartz	Lid	Black	10x1	45 mm	0.25 ml	170 to 2,600 nm	1	\$175.00
S-90-358Q	2100-358Q	Quartz	Stopper	Black	10x2 x 5 mm	45 mm	0.100 ml	170 to 2,600 nm	1	\$195.00
S-90-359Q	2100-359Q	Quartz	Stopper	Black	10x2 x2.5 mm	45 mm	0.050 ml	170 to 2,600 nm	1	\$195.00