

# SELECTION GUIDE

## Refrigerators/Freezers



**HAZARDOUS ENVIRONMENT?  
SAFETY FIRST!**



**SAFE COLD STORAGE!**



**SIZES FOR YOUR EVERY NEED!**



**OUR LARGEST COLD STORAGE  
UNITS!**

REFRIGERATORS	Product	Total Volume cu. ft.	Cabinet Style	Page
<b>HAZARDOUS ENVIRONMENT? SAFETY FIRST!</b>	Explosion Proof	5.6 20	Undercounter Upright	266
<b>SAFE COLD STORAGE!</b>	Flammable Material Storage	5.6 20	Undercounter Upright	268
<b>SIZES FOR YOUR EVERY NEED!</b>	General Purpose	5.6 20	Undercounter Upright	270
<b>READY TO USE!</b>	Chromatography	23, 47, 69	Upright	273
<b>OUR LARGEST COLD STORAGE UNITS!</b>	Large Capacity	23, 49, 72	Upright	275
<b>FREEZERS</b>				
<b>HAZARDOUS ENVIRONMENT? SAFETY FIRST!</b>	Explosion Proof	5.6 6, 18 20	Undercounter Chest Upright	266
<b>SAFE COLD STORAGE!</b>	Flammable Material Storage	5.6 20	Undercounter Upright	268
<b>SIZES FOR YOUR EVERY NEED!</b>	General Purpose	1.9 5.6 6, 18 20	Undercounter Chest Upright Benchtop	270
<b>REACH IN FOR EASY ACCESS!</b>	Large Capacity	23, 49, 72	Upright	274
<b>REFRIGERATOR/FREEZER COMBINATIONS</b>				
<b>HAZARDOUS ENVIRONMENT? SAFETY FIRST!</b>	Explosion Proof	5.6 10, 11	Upright Undercounter	266
<b>SAFE COLD STORAGE!</b>	Flammable Material Storage	5.6 10, 11	Upright Undercounter	268
<b>SIZES FOR YOUR EVERY NEED!</b>	General Purpose	5.6 10, 14.4	Undercounter Upright	270

# Refrigerators / Freezers

## EXPLOSION - PROOF

**Lab-Line**

Model 3560



### HAZARDOUS ENVIRONMENT? SAFETY FIRST!

- — Engineered to prevent triggering an explosion inside or outside the unit
- — Suitable for use in Class I, Group C and D hazardous environments
- — Meets OSHA and National Fire Protection Association safety requirements



### PRODUCT DESCRIPTION

- Available in a wide range of sizes and styles, including uprights, combination refrigerator/freezers, undercounter and horizontal chest units.
- Thermostat and compressor components are isolated within a vapor-proof enclosure.
- Electrical connections are hardwired to prevent sparking and accidental disconnection.
- Magnetic vinyl gaskets for positive door seal.
- CFC-free insulation and refrigerant.

#### Model 3560

- 20 cubic ft. corrosion-resistant model designed to store caustic materials and acids.
- Five adjustable stainless steel shelves.
- Enameled steel interior and exterior with stainless steel door liner for greater corrosion resistance.
- Built-in key lock.
- UL listed, CSA certified.

#### Model 3551

- 11 cubic ft. combination refrigerator/freezer with separate doors.
- Three adjustable shelves plus four door shelves for maximum storage space.
- Hydraulic thermostat; adjustable natural air flow vent.
- Manual defrost.
- UL listed, CSA certified.

#### Model 3552

- 20 cubic ft. upright freezer.
- Built-in evaporators in shelves for even cooling.
- Four permanently positioned wire shelves.
- Door equipped with key lock.
- Manual defrost.
- UL listed, CSA certified.

#### Model 3557

- 5.6 cubic ft. undercounter refrigerator.
- Hydraulic thermostat.
- Three adjustable wire shelves plus three door shelves.
- Off cycle self defrost.
- UL listed.

#### Model 3557-4

- 5.6 cubic ft. undercounter refrigerator/freezer.
- Hydraulic thermostat.
- Two adjustable wire shelves, three door shelves.
- ABS plastic interiors with enameled steel exterior.
- Manual defrost.

#### Model 3557-2

- 5.6 cubic ft. undercounter freezer.
- Hydraulic thermostat.
- Three fixed shelves, three door shelves.
- Aluminum evaporator coil.
- Off cycle defrost.

#### Model 3555

- 18 cubic ft. horizontal chest freezer.
- Storage basket.
- Built-in key lock in lid.
- Manual defrost.

#### Model 3555-6

- 7 cubic ft. horizontal chest freezer.
- Storage basket.
- Cover with key lock.
- UL listed, CSA certified.

#### Model 3566

- 20 cubic ft. one-door refrigerator.
- Copper evaporator coil for efficient thermal transfer.
- Five adjustable shelves.
- Door with built-in key lock.
- UL listed, CSA certified.

#### Model 3559

- 10 cubic ft. single door refrigerator/freezer with interior freezer.
- Three adjustable shelves plus four door shelves.
- Hydraulic thermostat.
- Manual defrost.
- UL listed, CSA certified.

Model 3551



# Refrigerators / Freezers

## EXPLOSION-PROOF

**Lab-Line**

### PRODUCT SPECIFICATIONS

Model #	Chamber Dimensions Inches (cm) Refrigerator Section			Chamber Dimensions Inches (cm) Freezer Section			Total Volume Cubic ft.	Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D		W	H	D	
3551, 3551-1	20.75 (53)	41.5 (105)	14 (36)	19 (48)	12.75 (32)	15 (38)	11	23.25 (59)	63 (160)	24.5 (62)	150(68)
3552, 3552-1		N/A		27 (69)	57 (145)	18 (46)	20	32 (81)	71.5 (182)	33 (84)	240 (109)
3555, 3555-1		N/A		50.25 (128)	28 (71)	22.25 (57)	18	60.25 (153)	35 (89)	29.5 (75)	200 (91)
3555-6, 3555-7		N/A		26.25 (67)	28 (71)	16.5 (42)	6	31.5 (80)	35 (89)	26.25 (67)	125 (57)
3557	19 (48)	27.5 (70)	13 (33)		N/A		5.6	23.5 (59)	33.5 (85)	23.75 (60)	120 (55)
3557-2		N/A		18 (46)	27.25 (70)	15 (38)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3557-4	18 (46)	21 (53)	15.25 (39)	15.25 (39)	5.25 (13)	11 (28)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3559, 3559-1	19.75 (50)	32 (81)	17.5 (44)	17.25 (44)	8 (20)	12 (30)	10	23 (58)	58 (147)	23.5 (60)	150 (68)
3560, 3560-1	27 (69)	55 (140)	20 (51)		N/A		20	32 (81)	71.5 (182)	28.5 (72)	280 (127)
3566, 3566-1	27 (69)	55 (140)	20 (51)		N/A		20	32 (81)	71.5 (182)	28.5 (72)	280 (127)

### ORDERING INFORMATION

Model #	Electrical			Temperature Range °C (°F)		Freezer Section	US List Price
	Volts	Hertz	Watts	Refrigerator Section	Refrigerator Section		
3551	120	60	360	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)		\$3274.00
3551-1	240	50	360	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)		\$3932.00
3552	120	60	516	N/A	-23° to -5° (-9° to 23°)		\$3769.00
3552-1	240	50	516	N/A	-23° to -5° (-9° to 23°)		\$4447.00
3555	120	60	516	N/A	-28° to -18° (-18° to 0°)		\$4162.00
3555-1	240	50	516	N/A	-28° to -18° (-18° to 0°)		\$4753.00
3555-6	120	60	276	N/A	-21° to -4° (-5° to 24°)		\$3233.00
3555-7	240	50	276	N/A	-21° to -4° (-5° to 24°)		\$3907.00
3557	120	60	160	2° to 7° (36° to 45°)	N/A		\$2295.00
3557-2	120	60	200	N/A	-26° to -15° (-15° to 5°)		\$2295.00
3557-4	120	60	160	2° to 7° (36° to 45°)	-26° to -6° (-15° to 21°)		\$2295.00
3559	120	60	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)		\$3101.00
3559-1	240	50	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)		\$3749.00
3560	120	60	480	-5° to 15° (23° to 59°)	N/A		\$5212.00
3560-1	240	50	480	-5° to 15° (23° to 59°)	N/A		\$5885.00
3566	120	60	480	-5° to 15° (23° to 59°)	N/A		\$2963.00
3566-1	240	50	480	-5° to 15° (23° to 59°)	N/A		\$3198.00

Model 3555



# Refrigerators / Freezers

## FLAMMABLE MATERIAL STORAGE

**Lab-Line**

Model 3552-10



### SAFE COLD STORAGE!

- — Designed for safe cold storage of volatile materials
- — Engineered with no internal electrical components which could trigger an explosion
- — Meets OSHA and National Fire Protection Association safety requirements



### PRODUCT DESCRIPTION

- Provide a wide variety of units to keep your lab safe from potentially destructive explosions.
- Available in many sizes and styles, including uprights, combination refrigerator/freezers, and undercounter models.
- Thermostat located on exterior of cabinet to eliminate interaction with stored materials.
- Magnetic vinyl gaskets for positive door seal.
- CFC-free insulation and coolant.

#### Model 3551-10

- 11 cubic ft. refrigerator/freezer with separate doors.
- Three adjustable shelves and four door shelves.
- Hydraulic thermostat; adjustable natural air flow vent.
- ABS plastic interior with enameled steel exterior.
- Manual defrost.
- UL listed.

#### Model 3552-10

- 20 cubic ft. upright freezer.
- Built-in evaporators in shelves for even cooling.
- Four permanently positioned wire shelves.
- Enameled steel exterior with ABS plastic interior.
- Door equipped with key lock.
- Manual defrost.
- UL listed.

#### Model 3556-4

- 5.6 cubic ft. undercounter freezer.
- Hydraulic thermostat.
- Three fixed wire shelves, three door shelves.
- Aluminum evaporator coil.
- Off cycle defrost.

#### Model 3556-6

- 5.6 cubic ft. undercounter refrigerator/freezer.
- Hydraulic thermostat.
- Two adjustable wire shelves, three door shelves.
- Manual defrost.

#### Model 3556

- 5.6 cubic ft. undercounter refrigerator.
- Hydraulic thermostat.
- Three adjustable wire shelves plus three door shelves.
- Aluminum evaporator coil.
- Off cycle self defrost.
- Key locking door.
- UL listed.

#### Model 3559-10

- 10 cubic ft. single door refrigerator with interior freezer.
- Three adjustable shelves plus four door shelves.
- Hydraulic thermostat.
- ABS plastic interior with enameled steel exterior.
- Manual defrost.
- UL listed.

Model 3556-4



#### Model 3566-10

- 20 cubic ft. one-door refrigerator.
- Five adjustable shelves.
- Enameled steel exterior.
- Manual defrost.
- UL listed.

#### Model 3560-10

- 20 cubic ft. corrosion-resistant model designed to store caustic materials and strong acids.
- Stainless steel evaporators.
- Five adjustable stainless steel shelves.
- Enameled steel interior and exterior with stainless steel door liner for greater corrosion resistance.
- Door comes with built-in key lock.
- Manual defrost.
- UL listed.

# Refrigerators/Freezers

FLAMMABLE MATERIAL STORAGE

Lab-Line

## PRODUCT SPECIFICATIONS

Model #	Chamber Dimensions Inches (cm) Refrigerator Section			Chamber Dimensions Inches (cm) Freezer Section			Total Volume Cubic ft.	Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D		W	H	D	
3551-10, 3551-11	20.75 (53)	41.5 (105)	14 (36)	19 (48)	12.75 (32)	15 (38)	11	23.25 (59)	63 (160)	24.5 (62)	150 (68)
3552-10, 3552-11		N/A		27 (69)	57 (145)	18 (46)	20	32 (81)	71.5 (182)	33 (84)	240 (109)
3556	19 (48)	27.5 (70)	13 (33)		N/A		5.6	23.5 (59)	33.5 (85)	23.75 (60)	120 (55)
3556-4		N/A		18 (46)	27.5 (70)	15 (38)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3556-6	18 (46)	21 (53)	15.25 (39)	15.25 (39)	5.25 (13)	11 (28)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3559-10, 3559-11	19.75 (50)	32 (81)	17.5 (45)	17.25 (44)	8 (20)	12 (30)	10	23 (58)	58 (147)	23.5 (60)	150 (68)
3560-10, 3560-11	27 (69)	55 (140)	20 (51)		N/A		20	32 (81)	71.5 (182)	28.5 (72)	280 (127)
3566-10, 3566-11	27 (69)	55 (140)	20 (51)		N/A		20	32 (81)	71.5 (182)	28.5 (72)	280 (127)

## ORDERING INFORMATION

Model #	Electrical			Temperature Range °C (°F)		Freezer Section	US List Price
	Volts	Hertz	Watts	Refrigerator Section	Refrigerator Section		
3551-10	120	60	360	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)		\$2815.00
3551-11	240	50	360	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)		\$3463.00
3552-10	120	60	516	N/A	-23° to -5° (-9° to 23°)		\$3203.00
3552-11	240	50	516	N/A	-23° to -5° (-9° to 23°)		\$3896.00
3556	120	60	160	2° to 7° (36° to 45°)		N/A	\$2397.00
3556-4	120	60	200	N/A	-26° to -15° (-15° to 5°)		\$2397.00
3556-6	120	60	160	2° to 7° (36° to 45°)	-26° to -6° (-15° to 21°)		\$2397.00
3559-10	120	60	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)		\$2621.00
3559-11	240	50	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)		\$3269.00
3560-10	120	60	480	-5° to 15° (23° to 59°)		N/A	\$5059.00
3560-11	240	50	480	-5° to 15° (23° to 59°)		N/A	\$5738.00
3566-10	120	60	480	-5° to 15° (23° to 59°)		N/A	\$2866.00
3566-11	240	50	480	-5° to 15° (23° to 59°)		N/A	\$3101.00

Model 3559-10



# Refrigerators / Freezers

## GENERAL PURPOSE

**Lab-Line**

Model 3767



### SIZES FOR YOUR EVERY NEED!

- — Available as all-refrigerator, combo refrigerator/freezer, upright freezer, and chest freezer models in many popular sizes
- — Adjustable shelves or baskets included
- — CFC-free insulation and coolant



### PRODUCT DESCRIPTION

#### Model 3750

- 5.6 cubic ft. undercounter refrigerator/freezer.
- Adjustable thermostatic control.
- Two adjustable wire shelves, three door shelves.
- ABS plastic interior with enameled steel exterior.
- Manual defrost.
- UL listed.

#### Model 3751

- 5.6 cubic ft. undercounter refrigerator.
- Off-cycle defrost.
- Adjustable hydraulic thermostat control.
- Three adjustable wire shelves, three door shelves.
- Key lock.
- UL listed.

#### Model 3752

- 5.6 cubic ft. undercounter freezer.
- Hydraulic thermostat.
- Three fixed shelves.
- Off cycle defrost.

#### Model 3761

- 7 cubic ft. chest-style freezer.
- Includes storage basket.
- Textured enameled steel exterior.
- Key lock in lid.
- Manual defrost.

#### Model 3762

- 18 cubic ft. chest-style freezer.
- Largest chest-style model.
- Storage basket included.
- Textured enameled steel exterior.
- Key lock in lid.

#### Model 3763

- 14 cubic ft. upright two-door combination refrigerator/freezer.
- Hydraulic thermostat control with separate "cold air flow" control for freezer.
- Three adjustable shelves in refrigerator.
- Magnetic vinyl gaskets for positive door seal.
- Automatic defrost.
- UL listed.

#### Model 3764

- 10 cubic ft. upright combination refrigerator/freezer.
- Single-door model with internal freezer compartment.
- Three adjustable wire shelves.
- Manual defrost.

#### Model 3766

- 20 cubic ft. frost-free refrigerator.
- Five permanently mounted shelves plus storage basket.
- Hydraulic thermostat.
- Automatic defrost.
- Locking door comes standard.

#### Model 3767

- 20 cubic ft. upright freezer.
- Four fixed wire shelves and six door shelves.
- Door comes with built-in key lock.
- Manual defrost.
- Enameled steel interior and exterior.

#### Model 3770

- 1.9 cubic ft. bench top combination refrigerator/freezer.
- Single door with interior freezer compartment.
- One removable shelf included.
- Manual defrost.
- UL listed.

Model 3751



**48 Hour or Sooner  
Express Shipping  
Guaranteed!\***



\*see inside front cover for details

# Refrigerators/Freezers

## GENERAL PURPOSE

**Lab-Line**

### PRODUCT SPECIFICATIONS

Model #	Chamber Dimensions Inches (cm) Refrigerator Section			Chamber Dimensions Inches (cm) Freezer Section			Total Volume Cubic ft.	Exterior Dimensions* Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D		W	H	D	
3750, 3750-1	18 (46)	21 (53)	15.25 (39)	15.25 (39)	5.25 (13)	11 (28)	5.6	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3751, 3751-1	19 (48)	27.5 (70)	13 (33)		N/A		5.6	23 1/2(59)	33.3 (85)	23.75 (60)	120 (55)
3752, 3752-1		N/A		18 (46)	27.5 (70)	15 (38)	5.60	21.25 (54)	33.5 (85)	25.25 (64)	120 (55)
3761, 3761-1		N/A		26.25 (67)	28 (71)	16.5 (42)	7	31 1/2(80)	35 (89)	26.25 (67)	125 (57)
3762, 3762-1		N/A		50.25 (128)	28 (71)	22.25 (57)	18	60.25 (153)	35 (89)	29.5 (75)	200 (91)
3763, 3763-1	25 (64)	36.5 (93)	15.5 (39)	24 (61)	15.25 (39)	16 (41)	14.4	28 (71)	61.5 (156)	28 (71)	190 (86)
3764, 3764-1	19.75 (50)	32 (81)	17.5 (45)	17.25 (44)	8 (20)	12 (30)	10	23 (58)	58 (147)	23.5 (60)	150 (68)
3766, 3766-1	27 (69)	57 (145)	17.75 (45)		N/A		20	32 (81)	71 (180)	28.5 (72)	280 (129)
3767, 3767-1		N/A		27 (69)	57 (145)	18 (46)	20	32 (81)	71.5 (182)	33 (84)	240 (109)
3770	16.25 (41)	15.25 (39)	11 (28)	6.25 (16)	6 (15)	8.25 (21)	1.9	18 1/2(47)	19 (48)	19 (48)	52 (24)

\* Upon installation, we recommend clearance of 3" in rear, 2" on each side and 12" above unit.

Model 3766



# Refrigerators / Freezers

## GENERAL PURPOSE

**Lab-Line**

Model 3750



### ORDERING INFORMATION

Model #	Electrical			Temperature Range °C(°F)		US List Price
	Volts	Hertz	Watts	Refrigerator Section	Freezer Section	
3750	120	60	160	2° to 7° (36° to 45°)	-26° to -6° (-15° to 21°)	\$1209.00
3750-1	240	50	160	2° to 7° (36° to 45°)	-26° to -6° (-15° to 21°)	\$1255.00
3751	120	60	160	2° to 7° (36° to 45°)	N/A	\$1209.00
3751-1	240	50	160	2° to 7° (36° to 45°)	N/A	\$1255.00
3752	120	60	200	N/A	-26° to -15° (-15° to 5°)	\$1209.00
3752-1	240	50	200	N/A	-26° to -15° (-15° to 5°)	\$1255.00
3761	120	60	276	N/A	-21° to -4° (-5° to 24°)	\$1148.00
3761-1	240	50	276	N/A	-21° to -4° (-5° to 24°)	\$1499.00
3762	120	60	516	N/A	-28° to -18° (-18° to 0°)	\$1969.00
3762-1	240	50	516	N/A	-28° to -18° (-18° to 0°)	\$2198.00
3763	120	60	420	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)	\$1464.00
3763-1	240	50	420	0° to 10° (32° to 50°)	-20° to 0° (-4° to 32°)	\$1800.00
3764	120	60	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)	\$1321.00
3764-1	240	50	360	1° to 10° (34° to 50°)	-17° to -6° (-1° to 21°)	\$1658.00
3766	120	60	480	-5° to 15° (23° to 59°)	N/A	\$2387.00
3766-1	240	50	480	-5° to 15° (23° to 59°)	N/A	\$2632.00
3767	120	60	516	N/A	-23° to -5° (-9° to 23°)	\$1504.00
3767-1	240	50	516	N/A	-23° to -5° (-9° to 23°)	\$1805.00
3770	120	60	132	0° to 7° (32° to 45°)	-5° to 0° (23° to 32°)	\$393.00

Model 3763





# Refrigerators/Freezers

## CHROMATOGRAPHY REFRIGERATORS

**Lab-Line**

Model 3796



### READY TO USE!

- — Three models available—23, 47, and 69 cu. ft.
- — CFC-free insulation and coolant
- ★ — 60 Month Warranty on Compressor



### PRODUCT SPECIFICATIONS

Model #	Capacity Cu. Ft.	Access Ports 2" Diameter	Interior Dimensions Inches (cm)			Exterior Dimensions* Inches (cm)			Shipping Weight Lb. (kg)
			W	H	D	W	H	D	
3790, 3790-1	23	1	24 (61)	61.5 (157)	27 (69)	27.25 (70)	78.25 (199)	29.5 (75)	311 (141)
3793, 3793-1	47	2	50.75 (129)	61.25 (156)	24.5 (63)	54 (138)	78.25 (199)	29.5 (75)	495 (225)
3796, 3796-1	69	2	75 (191)	61.5 (157)	27 (69)	78 (199)	79.25 (202)	29.5 (75)	678 (308)

\*Exterior height does not include leveling feet. Add 6" for feet.

### PRODUCT DESCRIPTION

- Temperature uniformity of  $\pm 0.5^\circ$  at  $4^\circ\text{C}$ .
- Adjustable thermostat control from  $1^\circ$  to  $10^\circ\text{C}$ .
- Triple tempered glass doors.
- 2" foamed-in-place polyurethane insulation.
- Two 2" diameter PVC access ports for gassing applications.
- White vinyl exterior coating is durable and chemically resistant.
- Interior construction is .040 aluminum with baked-on enamel finish.
- Vapor-proof outlets.
- Off-cycle defrost.
- Optional chromatography support rods and casters available.

### ORDERING INFORMATION

Model #	Electrical		Duplex Outlets	Chromatography			US List Compressor Price
	Volts	Hertz		Rods	Shelves	Doors	
3790	120	60	1	N/A	3	1 (swinging)	1/3 hp \$3583.00
3790-1	240	50	1	N/A	3	1 (swinging)	1/3 hp \$3843.00
3793	120	60	2	1	6	2 (sliding)	1/2 hp \$5585.00
3793-1	240	50	2	1	6	2 (sliding)	1/2 hp \$5760.00
3796	120	60	2	2	9	3 (sliding)	1/2 hp \$7020.00
3796-1	240	50	2	2	9	3 (sliding)	1/2 hp \$7342.00

### ACCESSORIES

Model #	Description
C	Casters
CR	Chromatography Rods
DT	Digital Temperature Readout
P	Access Ports
V	Vapor-proof Duplex Outlet

When ordering, add accessory letter to end of model number, i.e.: 3790C

# Refrigerators / Freezers

## LARGE CAPACITY FREEZERS

**Lab-Line**

Model 3785



### REACH IN FOR EASY ACCESS!

- — Storage of biological materials and other nonflammable substances
- — Three models available—from 23, 49 and 72 cu. ft.
- — CFC-free insulation and coolant
- ★ — 60 Month Warranty on Compressor



### PRODUCT DESCRIPTION

- Automatic defrost.
- Key locks on all doors.
- Adjustable thermostat control from -17°C to -23°C.
- Temperature uniformity of  $\pm 0.5^\circ\text{C}$  at  $-23^\circ\text{C}$ .
- .040 polished aluminum exterior.
- 22 gauge stainless steel floor.
- Covered interior corners for ease of cleaning.
- Foamed-in-place insulation.
- PVC coated adjustable shelves.
- Compressor bottom-mounted for easy access.

Model 3788



### PRODUCT SPECIFICATIONS

Model #	W	Interior Dimensions Inches (cm)			Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
		H	D	W	H	D		
3782 3782-1	23.75 (60)	61.5 (156)	26 (66)	27 (69)	78.25 (199)	29.5 (75)	286 (130)	
3785 3785-1	50.5 (128)	61.5 (156)	26 (66)	54 (137)	78.25 (199)	29.5 (75)	425 (193)	
3788 3788-1	74.75 (190)	61.5 (156)	26 (66)	78 (198)	79.25 (201)	29.5 (75)	598 (271)	

### ORDERING INFORMATION

Model #	Electrical		Capacity (cu. ft.)	Number of Shelves (sq. ft.)	Solid Doors	Compressor	US List Price
	Volts	Hertz					
3782	120	60	23	3 (13)	1	1/2 hp	\$4621.00
3782-1	240	50	23	3 (13)	1	1/2 hp	\$4621.00
3785	120	60	49	6 (29)	2	3/4 hp	\$6217.00
3785-1	240	50	49	6 (29)	2	3/4 hp	\$6217.00
3788	240	50	72	9 (43)	3	3/4 hp	\$8619.00
3788-1	240	50	72	9 (43)	3	3/4 hp	\$8619.00

# Refrigerators / Freezers

## LARGE CAPACITY REFRIGERATORS

**Lab-Line**

MODEL 3784



### OUR LARGEST COLD STORAGE UNITS!

- — Storage of biological materials and other nonflammable substances
- — Six models available—from 23, 49, and 72 cu. ft.
- — CFC-free insulation and coolant
- ★ — 60 Month Warranty on Compressor

### PRODUCT DESCRIPTION

- Available in glass or solid door styles.
- Key locks on all doors.
- Adjustable thermostatic control from 1° to 10°C.
- Temperature uniformity of  $\pm 0.5^{\circ}\text{C}$  at 4°C.
- .040 polished aluminum exterior.
- 22 gauge stainless steel floor.
- Foamed-in-place insulation.
- PVC coated adjustable shelves.
- Compressor bottom-mounted for easy access.
- UL listed.

### ACCESSORIES

Model #	Description
B	Digital Alarm with Remote Contacts
L	6" Adjustable Stainless Steel Legs
DT	Digital Temperature Readout
S	Stainless Steel Interior and Exterior

When ordering, add accessory letter to end of model number, i.e.:3782DT

### PRODUCT SPECIFICATIONS



Model #	Interior Dimensions Inches (cm)			Exterior Dimensions Inches (cm)			Shipping Weight Lb. (kg)
	W	H	D	W	H	D	
3780 3780-1	23.75 (60)	61.5 (157)	26 (66)	27 (69)	78.25 (198)	31.75 (81)	245 (111)
3781 3781-1	23.75 (60)	61.5 (157)	26 (66)	27 (69)	78.25 (198)	31.75 (81)	277 (126)
3783 3783-1	50.5 (129)	61.5 (157)	26 (66)	54 (138)	78.25 (198)	31.75 (81)	415 (188)
3784 3784-1	50.5 (129)	61.5 (157)	26 (66)	54 (138)	78.25 (198)	31.75 (81)	480 (218)
3786 3785-1	74.75 (190)	61.5 (157)	26 (66)	78 (198)	78.25 (198)	31.75 (81)	585 (265)
3787 3787-1	74.75 (190)	61.5 (157)	26 (66)	78 (198)	78.25 (198)	31.75 (81)	655 (297)

\*Exterior height does not include castors. Add 6" for castor height.

### ORDERING INFORMATION

Model #	Electrical		Capacity (cu.ft.)	Number of Shelves (sq.ft.)	Solid Doors	Compressor	US List Price
	Volts	Hertz					
3780	120	60	23	3 (13)	Solid	1/3 hp	\$3667.00
3780-1	240	50	23	3 (13)	Solid	1/3 hp	\$3667.00
3781	120	60	23	3 (13)	Glass	1/3 hp	\$4167.00
3781-1	240	50	23	3 (13)	Glass	1/3 hp	\$4167.00
3783	120	60	49	6 (29)	Solid	1/2 hp	\$4891.00
3783-1	240	50	49	6 (29)	Solid	1/2 hp	\$4891.00
3784	120	60	49	6 (29)	Glass	1/2 hp	\$5748.00
3784-1	240	50	49	6 (29)	Glass	1/2 hp	\$5748.00
3786	120	60	72	9 (43)	Solid	1/2 hp	\$6579.00
3786-1	240	50	72	9 (43)	Solid	1/2 hp	\$6579.00
3787	120	60	72	9 (43)	Glass	1/2 hp	\$7864.00
3787-1	240	50	72	9 (43)	Glass	1/2 hp	\$7864.00