

INSTRUMENT TRAYS



A. 1002 and 1002-2



A. 63 and 63-2



A. E1650



A. E1082



A. E2304



A. E904

A. Instrument Trays And Flat/Solid Covers

Wide assortment of heavy duty, seamless trays have covered corners for easy sterilization. Covers have recessed strap handles for ease in handling and storing.

Part #	Dimensions (in.) L x W x H	Case Lot - Wt.
43	5 ³ / ₈ x 4 x 2 ⁵ / ₈	12-5.0 lbs.
43-2	Cover - Fits 43 Tray	12-3.5 lbs.
63	6 ¹ / ₈ x 6 ¹ / ₈ x 2 ⁵ / ₈	12-7.3 lbs.
63-2	Cover - Fits 63 Tray	12-5.1 lbs.
952	8 ⁷ / ₈ x 5 x 2	12-7.3 lbs.
953	8 ⁷ / ₈ x 5 x 3	12-8.6 lbs.
952-2	Cover - Fits above two trays	12-5.8 lbs.
1002	10 x 6 ¹ / ₂ x 2	12-10.3 lbs.
1002-2	Cover - Fits 1002 Tray	12-8.5 lbs.
1202	12 ¹ / ₄ x 7 ³ / ₄ x 2 ¹ / ₄	12-16.0 lbs.
1202-2	Cover - Fits 1202 Tray	12-11.3 lbs.
E1082	10 x 8 ¹ / ₄ x 2 ¹ / ₂	6-7.2 lbs.
E1084	10 x 8 ¹ / ₄ x 4	6-8.5 lbs.
1082-2	Cover - Fits above two trays	6-6.7 lbs.
E1650	16 ¹ / ₂ x 10 x 2 ¹ / ₂	6-12.5 lbs.
E1654	16 ¹ / ₂ x 10 x 4	4-9.7 lbs.
1650-2	Cover - Fits above two trays	6-9.3 lbs.
E902	6 ⁷ / ₈ x 4 ¹ / ₄ x 2	12-6.0 lbs.
E904	6 ⁷ / ₈ x 4 ¹ / ₄ x 4	12-8.1 lbs.
902-2	Cover - Fits above two trays	12-2.9 lbs.
E2302	13 ³ / ₈ x 12 ⁷ / ₈ x 2 ¹ / ₂	6-14.3 lbs.
E2304	13 ³ / ₈ x 12 ⁷ / ₈ x 4	6-17.3 lbs.
E2306	13 ³ / ₈ x 12 ⁷ / ₈ x 6	6-20.4 lbs.
2302-2	Cover - Fits above three trays	6-10.9 lbs.
E06062	6 ⁷ / ₈ x 6 ³ / ₈ x 2 ¹ / ₂	6-3.2 lbs.
E06064	6 ⁷ / ₈ x 6 ³ / ₈ x 4	6-5.0 lbs.
E06066	6 ⁷ / ₈ x 6 ³ / ₈ x 6	6-6.6 lbs.
602-2	Cover - Fits above three trays	6-2.4 lbs.
E10062	10 ³ / ₈ x 6 ³ / ₈ x 2 ¹ / ₂	6-5.3 lbs.
E10064	10 ³ / ₈ x 6 ³ / ₈ x 4	6-7.5 lbs.
E10066	10 ³ / ₈ x 6 ³ / ₈ x 6	6-9.5 lbs.
402-2	Cover - Fits above three trays	6-3.4 lbs.

ORDER PANS AND COVERS SEPARATELY.

B. Instrument Trays And Flat/Solid Covers

Part #	Dimensions (in.) L x W x H	Case Lot - Wt.
Standard Weight (22 Ga.)		
E12062	12 ³ / ₄ x 6 ⁷ / ₈ x 2 ¹ / ₂	6-6.5 lbs.
E12064	12 ³ / ₄ x 6 ⁷ / ₈ x 4	6-9.3 lbs.
E12066	12 ³ / ₄ x 6 ⁷ / ₈ x 6	6-11.7 lbs.

Heavy Weight (20 Ga.)

S12062	12 ³ / ₄ x 6 ⁷ / ₈ x 2 ¹ / ₂	6-7.5 lbs.
S12064	12 ³ / ₄ x 6 ⁷ / ₈ x 4	6-10.8 lbs.
S12066	12 ³ / ₄ x 6 ⁷ / ₈ x 6	6-13.5 lbs.
302-2	Cover - Fits above six trays	6-4.6 lbs.

ORDER PANS AND COVERS SEPARATELY.

INSTRUMENT TRAYS

B. Instrument Trays And Flat/Solid Covers, Cont...

Part #	Dimensions (in.) L x W x H	Case Lot - Wt.
--------	-------------------------------	----------------

Standard Weight (22 Ga.)

E12102	12 $\frac{3}{4}$ x 10 $\frac{3}{8}$ x 2 $\frac{1}{2}$	6-8.8 lbs.
E12104	12 $\frac{3}{4}$ x 10 $\frac{3}{8}$ x 4	6-12.4 lbs.
E12106	12 $\frac{3}{4}$ x 10 $\frac{3}{8}$ x 6	6-14.8 lbs.

Heavy Weight (20 Ga.)

S12102	12 $\frac{3}{4}$ x 10 $\frac{3}{8}$ x 2 $\frac{1}{2}$	6-10.4 lbs.
S12104	12 $\frac{3}{4}$ x 10 $\frac{3}{8}$ x 4	6-14.5 lbs.
S12106	12 $\frac{3}{4}$ x 10 $\frac{3}{8}$ x 6	6-17.3 lbs.
1311-2	Cover - Fits above six trays	6-8.0 lbs.

Standard Weight (22 Ga.)

E20122	20 $\frac{3}{4}$ x 12 $\frac{3}{4}$ x 2 $\frac{1}{2}$	6-17.0 lbs.
E20124	20 $\frac{3}{4}$ x 12 $\frac{3}{4}$ x 4	6-23.3 lbs.
E20126	20 $\frac{3}{4}$ x 12 $\frac{3}{4}$ x 6	6-27.7 lbs.

Heavy Weight (20 Ga.)

S20122	20 $\frac{3}{4}$ x 12 $\frac{3}{4}$ x 2 $\frac{1}{2}$	6-20.3 lbs.
S20124	20 $\frac{3}{4}$ x 12 $\frac{3}{4}$ x 4	6-26.4 lbs.
S20126	20 $\frac{3}{4}$ x 12 $\frac{3}{4}$ x 6	6-33.5 lbs.
2000-2	Cover - Fits above six trays	6-15.4 lbs.

ORDER PANS AND COVERS SEPARATELY.

C. Perforated Instrument Trays – P M P Series

Constructed of heavy duty 304 18-8 stainless steel, P M P perforated instrument trays are finished to provide a tough, clean and corrosion-resistant surface for long term use in a variety of applications.

Part #	Dimensions (in.) L x W x H	Case Lot - Wt.
--------	-------------------------------	----------------

Perforated Mayo Trays

For excellent steam sterilization without moisture build-up.

10F-PMT	10 x 6 $\frac{1}{2}$ x $\frac{3}{4}$	6-3.5 lbs.
13F-PMT	13 $\frac{11}{16}$ x 9 $\frac{13}{16}$ x $\frac{3}{4}$	6-9.9 lbs.
15F-PMT	15 $\frac{1}{8}$ x 10 $\frac{5}{8}$ x $\frac{3}{4}$	6-10.0 lbs.
17F-PMT	17 x 11 $\frac{5}{8}$ x $\frac{3}{4}$	6-12.5 lbs.
19F-PMT	19 $\frac{1}{8}$ x 12 $\frac{5}{8}$ x $\frac{3}{4}$	6-14.5 lbs.
1520-PMT	21 $\frac{1}{4}$ x 16 $\frac{1}{4}$ x $\frac{7}{8}$	6-22.8 lbs.

High-Sided Perforated Trays

High sides protect fragile instruments during set-up and procedures.

952-HPT	8 $\frac{1}{8}$ x 5 x 2	12-7.3 lbs.
1002-HPT	10 x 6 $\frac{1}{2}$ x 2	12-10.3 lbs.
1202-HPT	12 $\frac{1}{4}$ x 7 $\frac{3}{4}$ x 2 $\frac{1}{4}$	12-16.0 lbs.
S12102-HPT	12 $\frac{3}{4}$ x 10 $\frac{3}{8}$ x 2 $\frac{1}{2}$	6-10.4 lbs.
E1650-HPT	16 $\frac{1}{2}$ x 10 x 2 $\frac{1}{2}$	6-12.5 lbs.
S20122-HPT	20 $\frac{3}{4}$ x 12 $\frac{3}{4}$ x 2 $\frac{1}{2}$	6-20.3 lbs.

Flash Pans

End handles allow for safe, easy and secure removal from autoclave.

952-FPP	8 $\frac{1}{8}$ x 5 x 2	12-7.3 lbs.
1002-FPP	10 x 6 $\frac{1}{2}$ x 2	12-10.3 lbs.



B. E12064 Tray



B. E12104 Tray



B. E20124 Tray



C. Instrument Trays — P M P Series

INSTRUMENT TRAYS

D. Catheter Trays and Dome Covers

Heavy duty and seamless, these trays feature dome covers with strap handles for easy handling.

Part #	Dimensions (in.) L x W x H	Case Lot - Wt.
83	8¾ x 3¼ x 1½	12-4.5 lbs.
83-2	Cover for 83 Tray	12-4.0 lbs.
170-0	17¼ x 4½ x 2	12-15.0 lbs.
170-2	Cover for 170 Tray	12-10.3 lbs.
952	8¾ x 5 x 2	12-7.3 lbs.
953	8¾ x 5 x 3	12-8.6 lbs.
952D-2	Cover for 952 and 953 Trays	12-7.0 lbs.

ORDER PANS AND COVERS SEPARATELY.

E. Sterilization Trays/Bath Sets (Hot/Cold)

Tray units can be used for either hot or cold sterilization. Covers feature recessed strap handle for ease and safety in transporting, handling and storing.

Part #	Description	Dimensions (in.) L x W x H	Case Lot - Wt.
E20124	Solid Tray	20¾ x 12¾ x 4	6-23.3 lbs.
P20122	Perforated Tray	20¾ x 12¾ x 2½	6-15.4 lbs.
2000-2	Flat, Solid Cover		6-15.4 lbs.
E20126	Solid Tray	20¾ x 12¾ x 6	6-27.7 lbs.
P20124	Perforated Tray	20¾ x 12¾ x 4	6-14.8 lbs.
2000-2	Flat, Solid Cover		6-15.4 lbs.
E12104	Solid Tray	12¾ x 10¾ x 4	6-12.4 lbs.
P12102	Perforated Tray	12¾ x 10¾ x 2½	6-8.5 lbs.
1311-2	Flat, Solid Cover		6-8.0 lbs.
E12106	Solid Tray	12¾ x 10¾ x 6	6-14.8 lbs.
P12104	Perforated Tray	12¾ x 10¾ x 4	6-10.7 lbs.
1311-2	Flat, Solid Cover		6-8.0 lbs.

ORDER PANS AND COVERS SEPARATELY.

F. Sterilization Trays (Hot/Cold)

Tray units can be used for either hot or cold sterilization. Covers feature recessed strap handle for ease and safety in transporting, handling and storing.

Part #	Capacity	Dimensions (in.) L x W x H	Case Lot - Wt.
E20062	3½ Qt.	20¾ x 6¾ x 2½	6-9.9 lbs.
E20064	5¼ Qt.	20¾ x 6¾ x 4	6-12.4 lbs.
E20066	8¾ Qt.	20¾ x 6¾ x 6	6-16.5 lbs.
2006-2	Flat, Solid Cover - Fits all above pans		6-8.2 lbs.

ORDER PANS AND COVERS SEPARATELY.

G. Oblong Instrument/Drying Trays

The user of these trays benefits from coved corners and large, open bead which allow for easy handling, stacking and cleaning.

Part #	Description	Dimensions (in.) L x W x H	Case Lot - Wt.
10F	w/Rolled Bead	10 x 6½ x ¾	6-3.5 lbs.
13F	w/Rolled Bead	13 ¹¹ / ₁₆ x 9 ³ / ₁₆ x ¾	6-9.9 lbs.
15F	w/Rolled Bead	15 ¹ / ₈ x 10 ⁵ / ₈ x ¾	6-10.0 lbs.
17F	w/Rolled Bead	17 x 11 ⁵ / ₈ x ¾	6-12.5 lbs.
19F	w/Rolled Bead	19 ¹ / ₈ x 12 ⁵ / ₈ x ¾	6-14.5 lbs.



D. Catheter Trays and Dome Covers



E. Sterilization Trays/Bath Sets



F. Sterilization Tray



G. Oblong Instrument/Drying Trays

INSTRUMENT TRAYS

G. Oblong Instrument/Drying Trays, Cont...

Part #	Description	Dimensions (in.) L x W x H	Case Lot - Wt.
1520	w/Rolled Bead	21¼ x 16¼ x ¾	6-22.8 lbs.
1525	w/Flat Flange	25 x 15 x 5/8	6-24.5 lbs.
E12101		12¾ x 10¾ x 1¼	6-6.9 lbs.
E20121		20¾ x 12¾ x 1¼	6-14.9 lbs.

H. Perforated Instrument Trays

¼" diameter holes in tray bottom and sides enhance sterilizing liquid flow. Tray edges are rounded to aid in handling.

Part #	Dimensions (in.) L x W x H	Case Lot - Wt.
P20121	20¾ x 12¾ x 1¼	6-13.8 lbs.
P20122	20¾ x 12¾ x 2½	6-15.4 lbs.
P20124	20¾ x 12¾ x 4	6-16.8 lbs.
P20126	20¾ x 12¾ x 6	6-17.4 lbs.
P12101	12¾ x 10¾ x 1¼	6-6.8 lbs.
P12102	12¾ x 10¾ x 2½	6-8.5 lbs.
P12104	12¾ x 10¾ x 4	6-10.7 lbs.
P12106	12¾ x 10¾ x 6	6-14.2 lbs.

I-J. Instrument Trays without Covers

All trays feature flat bottoms, open beads and covered corners.

Part #	Dimensions (in.) L x W x H	Case Lot - Wt.
I. with Rolled Bead		
1302	13⅞ x 9⅞ x 2¼	6-11.9 lbs.
1402	14⅞ x 10⅞ x 2¼	6-13.3 lbs.
J. with Very Slight Flange		
4V	8¾ x 4⅝ x 3⅝	12-9.0 lbs.
6V	10⅝ x 5½ x 4	12-12.8 lbs.

K. Mesh Bottom Sterilization Trays

These sterilization trays feature a strong, 18-8 stainless steel mesh bottom with ⅛" perforations to protect delicate instrument tips while providing excellent steam permeation without condensation build up. The mesh bottom is offset from the tray by four strong, impact resistant, smooth perforated round ball feet to enhance drainage and drying. Extra sanitary, retractable handles set on top of tray bead promote easy grasping and lifting as well as conserve on space.

Part #	Dimensions (in.) L x W x H	Case Lot - Wt.
E20124-IST	20¾ x 12¾ x 4	6-23.3 lbs.
E1654-IST	16½ x 10 x 4	4-9.7 lbs.
E12104-IST	12¾ x 10¾ x 4	6-12.4 lbs.

L. Perforated Cleaning Baskets

This line of cleaning baskets boasts ¼" holes on bottom and sides to provide excellent penetration of cleaning solutions. Easy-to-clean retractable handles are set on top of basket to promote safe, easy handling. The bottom is offset from the surface by four "feet" to enhance draining and drying.

Part #	Dimensions (in.) L x W x H	Case Lot - Wt.
P20122-CB	20¾ x 12¾ x 2½	6-16.4 lbs.
P20124-CB	20¾ x 12¾ x 4	6-17.8 lbs.
P20126-CB	20¾ x 12¾ x 6	6-18.4 lbs.
P12102-CB	12¾ x 10¾ x 2½	6-9.5 lbs.
P12104-CB	12¾ x 10¾ x 4	6-11.7 lbs.
P12106-CB	12¾ x 10¾ x 6	6-15.2 lbs.



H. Perforated Instrument Tray



I. Instrument Tray without Cover/with Rolled Bead



J. Instrument Tray without Cover/with Very Slight Flange



K. Mesh Bottom Sterilization Trays



L. Perforated Cleaning Baskets