

WESCO

enjoy the VU

SERIES

1000

Edu VU™ 1000 Series

The Edu VU™ 1000 series educational microscopes combine state of the art technology with affordability to meet the standards and budgets of schools ranging from primary through college.

There are currently five models in the Edu VU™ 1000 Series with many additional accessories available to further customize your viewing needs.



Edu VU 1000 Series

SPECIFICATIONS	MODELS					
	1100	1110	1120	1200	1300	
Monocular observation body	•	•	•			
Dual monocular observation body				•		
Binocular observation body					•	
10X widefield eyepiece 18mm field of view	•	•	•	• (2 each)	• (2 each)	
Reverse angle revolving quadruple nosepiece	•	•	•	•	•	
DIN 4X, 10X and 40XR achromatic objectives	•	•				
DIN 4X, 10X, 40XR & 100XR oil achromatic objectives			•	•	•	
Acid and reagent resistant 100mm X 115mm stage plate with stage clips	•					
Acid and reagent resistant 100mm X 115mm stage plate with attachable X-Y mechanical stage		•				
Acid and reagent resistant 100mm X 115mm flat top graduated mechanical stage with low position 76mm X 50mm travel X-Y controls			•	•	•	
Built-in 0.65 N.A. fixed condenser with iris diaphragm and filter holder	•	•				
Focusable ABBE 1.25 N.A. condenser with iris diaphragm and swing out filter holder			•	•	•	
7 watt fluorescent illumination	•	•	•	•	•	
Electrical cord hanger	•	•	•	•	•	
Options						
Eyepiece pointer	•	•	•	•	•	
Five hole disc condenser	•	•				
Attachable mirror	•	•	•	•	•	

WESCO



Distributed By: