

# WESCO

*enjoy the VU*

SERIES

3000

## Stereo VU™ 3000 Series

*The Stereo VU™ 3000 series stereo microscopes are designed to meet the specifications and budgets of schools as well as the demands of industry with features and accessories that serve both markets.*

*There are currently twenty models in the Stereo VU™ 3000 Series.*



# Stereo VU 3000 Series

SPECIFICATIONS	MODELS							
	3100	3200	3300	3400	3110	3210	3310	3410
Binocular .65-4.5X zoom stereo observation body	•	•	•	•				
Trinocular .65-4.5X zoom stereo observation body					•	•	•	•
10X widefield eyepiece 23mm field of view	•	•	•	•	•	•	•	•
Table pole stand with focus mount	•				•			
Universal boom stand with double focus arm		•				•		
Table pole stand with incidental variable halogen and transmitted fluorescent illumination			•				•	
Table pillar stand with incidental variable halogen and transmitted fluorescent illumination					•			•

SPECIFICATIONS	FIXED STEREO MODELS											
	3120	3130	3140	3220	3230	3240	3320	3330	3340	3420	3430	3440
Binocular fixed 1X & 2X stereo observation body	•			•			•			•		
Binocular fixed 1X & 3X stereo observation body		•			•			•			•	
Binocular fixed 2X & 4X stereo observation body			•			•			•			•
10X widefield eyepiece 23mm field of view	•	•	•	•	•	•	•	•	•	•	•	•
Table pole stand with focus mount	•	•	•									
Universal boom stand with double focus				•	•	•						
Table pole stand with incidental variable halogen and transmitted fluorescent illumination							•	•	•			
Table pillar stand with incidental variable halogen and transmitted fluorescent illumination										•	•	•

# WESCO



Distributed By: