

# WESCO

SERIES

*enjoy the VU*

4000

## Lab VU™ 4000 Series

*The New Lab VU™ 4000 Series microscopes embody highly advanced infinity corrected optics and new ergonomic mechanical design for top performance from routine laboratory work to educational and research studies.*

*There are currently six models in the Lab VU™ 4000 Series with many additional accessories available to further customize your viewing needs.*



# Lab VU 4000 Series

SPECIFICATIONS	MODELS					
	4200	4220	4230	4300	4320	4330
Seidentopf binocular observation body	•	•	•			
Seidentopf trinocular observation body				•	•	•
10X widefield eyepieces 20mm field of view	•	•	•	•	•	•
Reverse angle revolving quadruple nosepiece	•	•	•	•	•	•
Infinity corrected 4X, 10X, 40XR & 100XR oil achromatic objective	•			•		
Infinity corrected 4X, 10X, 40XR & 100XR oil plan-achromatic objective		•			•	
Infinity corrected plan phase contrast PH10X PH20X, PH40XR and PH100XR oil objectives with turret phase condenser, telescope eyepiece and G.I filter			•			•
Coaxial coarse and fine focus with tension control	•	•	•	•	•	•
150mm X 150mm flat top graduated mechanical stage with low position 76mm X 50mm travel X-Y controls.	•	•	•	•	•	•
N.A. 1.25 Abbe condenser	•	•	•	•	•	•
6 volt 20 watt halogen illumination with variable light control	•	•	•	•	•	•
Field iris diaphragm		•	•		•	•

# WESCO



Distributed By: