

WESCO

enjoy the VU

SERIES

5000

Indu VU™ 5000 Series

The Indu VU™ 5000 series industrial microscopes are designed optically and mechanically to meet the demands of industry, making them the perfect tools for quality control, assembly line and O.E.M. applications.

There are currently seven models in the Indu VU™ 5000 Series with many additional accessories available to further customize your viewing needs.



Indu VU 5000 Series

SPECIFICATIONS	MODELS							
	5501B	5601B	5601BI	6501T	5601TI	5701B	5701T	
Binocular 5:1 zoom ratio stereo observation body Magnification: .7X - 3.5X	•							
Binocular 6:1 zoom ratio stereo observation body Magnification: 1X-6.3X		•						
Binocular 6:1 zoom with coaxial illumination			•					
Trinocular 6:1 zoom ratio stereo observation body				•				
Trinocular 6:1 zoom with coaxial illumination					•			
Binocular 7:1 zoom ratio stereo observation body Magnification: .65X - 4.5X						•		
Trinocular 7:1 zoom ratio stereo observation body							•	
10X widefield eyepiece 23mm field of view	•	•	•	•	•	•	•	
Table stand with focus mount and tension control	•	•	•	•	•	•	•	
6 volt 20 watt incidental halogen illumination system	•					•	•	
Options								
Universal boom stand	•	•	•	•	•	•	•	
Universal tilting drop post	•	•	•	•	•	•	•	
Universal bonder arm	•	•	•	•	•	•	•	
Universal focusing mount	•	•	•	•	•	•	•	
Table stand with variable incidental and transmitted halogen illumination	•	•	•	•	•	•	•	
Widefield 5X, 15X, 20X, 25X and 30X eyepieces	•	•	•	•	•	•	•	
0.3X, 0.5X, 0.75X, 1.5X and 2.0X objectives	•	•	•	•	•	•	•	
Fluorescent ring light and fiber optic illumination systems	•	•	•	•	•	•	•	

WESCO



Distributed By: