

Lab-Line Histology



Slide-Warmers



Model 26025

Product Description

- Uniform heat transfer over entire heating surface.
- Anodized black heavy gauge aluminum work top resists corrosion and provides high-contrast background for viewing specimen.
- Powder coated heavy-duty steel body construction is durable and easy to clean.
- Thermostat-controller regulates temperature to 65°C.
- Model 26025 features LED display and cover.
- Model 26007CE and Model 26022CE are CE certified.
- Models 26005 and 26020 are ETL marked.
- Models 26005 & 26007CE have a slide capacity of sixty-five 1" x 3" (2.54 cm x 7.6 cm) or thirty-two 2" x 3" (5.1 cm X 7.6 cm).
- Models 26020, 26022CE, 26025 & 26025-1 have a slide capacity of sixty-six 1" x 3" (2.54 x 7.6 cm) or thirty-three 2" x 3" (5.1 x 7.6 cm).

Our Large Size Economical Slide-Warmer!

- Temperature range from slightly above ambient to 65°C
- Two sizes—accommodates up to 66 slides
- Large slide capacity with space-saving design
- ☆ 24 Month Warranty



Model 26005



Application

- Designed for cytology, histology, and other biological applications in industrial and clinical lab environments.

ORDERING INFORMATION

Model #	Electrical (50/60 Hz)		Overall Dimensions Inches (cm)			US List Price
	Volts	Watts	W	H*	D	
26005 ¹	120	160	14 (36)	14 (36)	3.5 (9)	\$489.00
26007CE ¹	240	160	14 (36)	14 (36)	3.5 (9)	\$489.00
26020 ²	120	160	25 (64)	8 (20)	3.5 (9)	\$541.00
26022CE ¹	240	160	25 (64)	8 (20)	3.5 (9)	\$541.00
26025	120	160	25 (64)	8 (20)	4.25 (11)	\$725.00
26025-1	240	160	25 (64)	8 (20)	4.25 (11)	\$725.00

¹ CE

² ETL

*Height with lid open

ACCESSORIES

Model #	Description	US List Price
26023	Metal dust cover for models 26022 and 26020	\$65.00