

NEW

50ml Centrifuge Tube Refrigerator Racks

For Temperature-Stabilized Bench-Top Storage

This heavy duty all-polypropylene rack provides efficient and extended cooling of twenty-eight 15 ml or 50 ml centrifuge tubes in a chilled bath. The deep tray, holding approximately 2 liters (67 ounces) of ice and water, allows the tubes to rest on the bottom for maximum stability and

exposure to the bath. A flat perimeter rim provides space for labeling tape. The rack is steam autoclavable at 121°C (250°F) and stackable with identical racks when full or empty. To further protect tubes during transport or storage while also maintaining stackability, an optional polypropylene cover is available. Dimensions: 32 x 20 x 15cm H (12³/₄ x 7⁷/₈ x 5⁷/₈" H).



CATALOG NO.		PRICE/EACH
F18905-0050	Rack	\$39.00
F18905-0051	Rack cover	11.00



Rimless Centrifuge Tubes

For use in centrifuges up to 50,000g without refrigeration, these polypropylene tubes are rimless, round-bottomed, and autoclavable at 121°C (250°F). **12 per bag, 12 bags per case.**



CATALOG NO.	CAPACITY	O.D.	LENGTH	PRICE/BAG
F19807-0000	45ml	29mm	104mm	\$19.35
F19808-0000	90ml	32mm	162mm	31.70



Centrifuge & Test Tube

Excellent Chemical Resistance

Graduated, conical polyethylene tubes with rim are perfect for refrigerated use and at temperatures up to 52°C (125°F) continuously. **12 per bag, 6 bags per case.**



CATALOG NO.	CAPACITY	O.D.	LENGTH	PRICE/BAG
F19841-0000	15ml	17mm	117mm	\$18.55

Poxygrid® Centrifuge Bottle Rack

Supports Conical and Round Bottom Centrifuge Bottles

This rack holds large conical and round-bottom centrifuge bottles that cannot stand up by themselves (such as those made by Falcon®, Corning®, Nunc®, and Kimble®). The 102mm (4") high rack has two levels of square grids to hold 6 bottles up to 250ml (60mm diameter). The bottom rails hold the tip of the conical

bottle. Dimensions are 208 x 135mm (8¹/₆ x 5¹/₆"). Epoxy-coated steel wire. **A**



CATALOG NO.	PRICE/EACH
F19856-0250	\$18.50

Poxygrid® 50ml Centrifuge Tube Racks

50ml Multiple Tube Holders

Made of steam autoclavable, epoxy-coated steel wire, these racks hold 50ml plastic centrifuge tubes in holes 30mm square. **Packed per each, 12 per case.**



CATALOG NO.	PLACES	DIMENSIONS	PRICE/EACH
F18794-0000	8	149 x 76 x 89mm H (5 ⁷ / ₈ x 3 x 3 ¹ / ₂ " H)	\$25.25
F18794-0001	16	152 x 152 x 89mm H (6 x 6 x 3 ¹ / ₂ " H)	26.75
F18794-0002	36	229 x 229 x 89mm H (9 x 9 x 3 ¹ / ₂ " H)	38.75

For Racks to Hold Smaller Centrifuge Tubes, See Page xxx.

Poxygrid® Single 50ml Centrifuge Tube Holder

Autoclavable Tube Holder

This holder solves the problem of where to put a single 50ml conical tube so that it will not roll or spill. The holder is formed of steel wire with epoxy coating and supports the conical bottom of the tube as well as at a half height location. It is easy to use and convenient to grip and move. Dimensions are 76mm (3") wide x 64mm (2¹/₂") high. **A**



CATALOG NO.	PRICE/EACH
F18794-0003	\$9.25

NEW

Single 50ml Tube Holders

Leave Hands Free for Pipetting and Other Critical Tasks



Non-Grip Style

This handy base securely holds a 50ml conical centrifuge tube in place, and also allows the tube to be lifted freely. **5 per bag.**



Grip Style

50ml Tube Holder, Grip style securely holds the tube in the base and moves with the tube when it is lifted from the bench top. Tubes can easily be removed by holding the base and pulling the tube out. **5 per bag.**



CATALOG NO.	DESCRIPTION	COLOR	PRICE/BAG
F18795-0000	Non-Grip style	White	\$19.25
F18795-0001	Grip style	Blue	19.25