



Economy Support Stand, Polycarbonate

A sturdy and stable molded polycarbonate base supports a 46 cm (18") tall stainless steel support rod. Two hex nuts secure the 7.9 mm ($\frac{5}{16}$ ") diameter rod to the base. Offset rod mounting to the 20.8 cm ($8\frac{3}{16}$ ") diameter base leaves area for larger vessels. Useful temperature range +135°C (275°F) to -135°C (-211°F). Autoclave at 121°C (250°F). The unit is shipped unassembled. Shipping weight 1.66 lbs.



CATALOG NO.	PRICE/EACH
F18302-0000	\$17.95



Burette Clamp - Spring Action

Spring-action opening clamp with a flat wing-head screw for tightening. The stainless steel jaws are covered with rubber tubing and have an opening range of 0 to 25mm (1"). Overall clamp length is 22cm (9").



CATALOG NO.	PRICE/EACH
H18354-0000	\$37.20

Boshead


Rod Holder

This open sided stainless steel boshead accepts rods up to 13mm ($\frac{1}{2}$ ") diameter. Rods are held at right angles to one another with flat wing head clamping screws.



CATALOG NO.	PRICE/EACH
H18350-0000	\$24.75

Use the Hot Hand® Protector to Grasp Hot or Cold Objects
See Page 122



Support Stands


Corrosion Resistant
Polypropylene

Excellent for use in severe corrosive conditions, these stands have no exposed metal parts - all surfaces are polypropylene. Measuring 20 x 30cm (8 x 12") x 5cm (2") thick, the bases have smooth surfaces for easy cleaning and weigh 2.5kg (5½ lbs) for stability. The 12.5mm ($\frac{1}{2}$ ") x 76cm (30") long rod screws into the base with a deep thread and sturdy tightening knob.



CATALOG NO.	DESCRIPTION	PRICE/EACH
H18303-0000	A. With rod at center	\$57.60
H18305-0000	B. With rod at end	59.40

For Automatic Burettes, See Page 97



Polypropylene Burette Clamps

Graduations Are Always Visible

These non-corroding, polypropylene burette clamps do not obscure graduations. Burette is held firmly in place by rubber pads and a concealed steel spring. Center clamp tightens to a standard support rod; 12.5mm ($\frac{1}{2}$ ") and 9.5mm ($\frac{3}{8}$ ") bushings included. **12 per case.**



CATALOG NO.	DESCRIPTION	PRICE/EACH
F18204-0000	A. Single	\$8.50
F18206-0000	B. Double	12.90

Funnel Holders

Supports Several Funnel Sizes

Non-corrosive high density polyethylene funnel holder and clamp includes two bushings for 12.5mm ($\frac{1}{2}$ ") or 9.5mm ($\frac{3}{8}$ ") diameter support rod. Outer ring supports 76 to 152mm (3 to 6") top diameter funnels; removable inner ring holds smaller funnels from 25 to 76mm (1 to 3") top diameter. **12 per case.**



CATALOG NO.	DESCRIPTION	PRICE/EACH
F18240-0000	Single - for one funnel	\$9.80
F18242-0000	Double - for two funnels	13.40

The Filters and Funnels Section Features a Complete Line of Scienceware® Funnels, See Pages 51-54



Multiple Chromatography Column Holders

Economical and Versatile

Flexible loops on a rigid plastic back securely hold as many as eight chromatography columns up to 30mm ($1\frac{1}{8}$ ") in diameter. Columns are held in or between the flexible loops. The back is mounted on two polyethylene blocks drilled to accommodate support stand rods. Simply slide the blocks over the support stand rods at the desired height and tighten the screws. Two support rods and stands are required to hold the rack. The rack is 311mm long x 64mm wide ($12\frac{1}{4}$ x $2\frac{1}{2}$ "). (Support stands not supplied.)



CATALOG NO.	PRICE/EACH
F18202-0000	\$63.90