

A. Flexible Arm with Weighted Base

Ideal for use with pH meters. The 23 x 16cm (9 x 6") base weighs 1.5kg (3.3lbs) and has rod holding wells with thumbscrews at the center of each edge. (Electrode Holder sold separately.)

CATALOG NO.	PRICE/EACH
F18315-2322	\$176.00

B. Flexible Gooseneck Arm

This Flexible Arm has molded acetal copolymer links with ball-and-socket joints. One end is designed to accept the F18317-1000 Electrode Holder while the other has a 12.7mm (½") diameter stainless steel pin for attachment to a clamp or base. The arm is 46cm (18") long x 25mm (1") in diameter. (Base and Electrode Holder sold separately.)

CATALOG NO.	PRICE/EACH
F18315-0000	\$65.00

C. Electrode Holder


The Electrode Holder snaps onto the end of the Flexible Arm affording 3 stepped wells for holding an electrode, a reference electrode and an ATC probe. The center well has an 18mm (¾") I.D. Side wells each have a 16.5mm (⅝") I.D. and the through holes are all 12.7mm (½") I.D. The Electrode Holder comes supplied with three reducing bushings for smaller diameters and five nylon cable ties for attaching wires to the Flexible Arm. (Base and Flexible Arm sold separately.)

CATALOG NO.	PRICE/EACH
F18317-1000	\$22.50

NEW



A. Triple



B. Clamp

28mm
(1")



C. Extender



D. Utility Tray



PiRack™ Pipettor Holder System

The sturdy PiRack™ Clamp-On Pipettor Holder supports many standard brands of pipettors and keeps them safe and within easy reach. Simply attach to any shelf 28mm (1") thick or benchtop edge using the large locking thumbscrew clamp. Add Extender to increase the number of holders at a work location. Made from polypropylene, the PiRack™ Pipettor Holder is strong enough to support the weight of even the heaviest pipettor. The Utility Tray, mounted between two clamps, is ideal for storing pipetting supplies or holding pens and pencils or other small items. It is perfectly sized for Scienceware® No Wire™ Half Racks available to hold tubes, vials, mini tubes and cuvettes. All components are chemical resistant and autoclavable at 121°C (250°F).



CATALOG NO.	DESCRIPTION	PRICE/EACH
F18954-0000	A. Triple (1 clamp and 2 extender)	\$20.00
F18954-0001	B. Clamps (2 clamps only)	20.00
F18954-0002	C. Extender (4 extenders only)	19.00
F18954-0500	D. PiRack™ Utility Tray (only)	8.95



Universal Implement Stand

Keep pipettes, thermometers and other tubular items from rolling off the bench.

This item falls into the handy gadget category. Just put it on the lab bench or desk, and it is ready to hold pipettes, stirring rods, burettes, thermometers, pencils, etc. These items will not roll away while held in the semi-circular notches. The stand is especially useful as a sterile hood tool. Molded of polypropylene, the support is ¾" (83mm) diameter x 2⅝" (59mm) high and is steam autoclavable at 121°C (250°F). **2 per bag. 12 bags per case.**



CATALOG NO.	PRICE/BAG
F18973-0000	\$13.65